

Learning Activity Lesson Plan

Author Name: Michael Loftice

Activity Title:	DNA Extraction from Kiwi
Goals/ Objectives:	Students will understand and perform laboratory protocol for extracting DNA from kiwi fruit.
Brief Overview:	In this lab, students will extract DNA from kiwi fruit using common household products. DNA will be extracted, spooled onto a glass rod. The DNA can be dyed with methylene blue and observed under a microscope, if desired.
Type of course this activity would be best suited for:	I use this lab in my Agricultural Biotechnology class.
Prerequisite knowledge:	<ul style="list-style-type: none"> • Laboratory Safety. • Cell components; cell organelles • Basic principles of genetic transfer
Time Required:	One class period; approximately 50 minutes.
Materials:	<ul style="list-style-type: none"> • 1 – teaspoon measure • 1 – ¼ teaspoon measure • 1 – glass measuring cup • 1 – plastic spoon & knife • 1 – small plastic plate • 2 – small plastic cups • 1 – #2 • coffee filter • 2 – plastic transfer pipettes • 1 – graduated cylinder • 1 – test tube • table salt • Distilled water • Clear shampoo without conditioner or clear dish washing detergent • Alcohol; 95% Ethanol preferred, but can use 95% Isopropyl (COLD) • ½ kiwi fruit, peeled • Hot water bath 50-60°C • Ice water bath
Methods:	<ol style="list-style-type: none"> 1. Place alcohol in ice bath at least one hour before lab 2. Prepare hot water bath / Cold water bath 3. Prepare solution of shampoo (1 tablespoon) and table salt (¼ teaspoon.) Add enough distilled water to make 100 ml. Stir slowly to dissolve the salt – avoid making foam. (This is enough solution for about 3 student stations.) 4. Peel and cut ½ kiwi into small pieces. 5. Place the kiwi pieces into a plastic cup. Cover with 30 ml of the solution made in step #3. Mash the kiwi against the side of the cup with the spoon. (The detergent dissolves the cell membranes and releases the DNA into the solution) 6. Place the cup into the hot water bath for about 10 minutes. (This further helps to

	<p>precipitate the fats and proteins out of the solution, leaving the DNA.)</p> <ol style="list-style-type: none"> 7. Cool the mixture in the ice water bath for about 5 minutes. 8. Using the #2 coffee filter, strain the mixture over the 2nd plastic cup. Filter until you have collected about 5 ml of liquid in the bottom of the 2nd cup. 9. Using a transfer pipette, add about 1 cm of alcohol into the test tube. (Alcohol should be as cold as possible.) 10. Using the 2nd transfer pipette, add 1 ml of the prepared kiwi liquid to the test tube. (DNA is not soluble in alcohol. The DNA will precipitate out into the alcohol layer.) 11. Let the solution in the test tube sit for about 5 minutes. DO NOT shake or stir the test tube! The white, stringy DNA will precipitate out into the alcohol layer. It has the appearance of white mucus.
<p>References : (Copyright Free)</p>	<p>Introduction to DNA Extractions Lana Hays Access Excellence Fellow (Activities Exchange) Simon Kenton H.S. 11132 Madison Pike Independence, KY 41051 AELHays@aol.com</p>
<p>Extension Activities:</p>	<ol style="list-style-type: none"> 1. DNA can be spooled onto a glass rod, placed on a slide, stained with methylene blue, and observed under a microscope. 2. <u>Students may also:</u> <ul style="list-style-type: none"> • Explain why the kiwi needed to be mashed. • Explain the role the detergent played in this experiment. • Explain the role salt played in this experiment. (Salt enables the DNA strands to come together.) • Describe the appearance of the DNA extracted in this experiment. • Identify the role of biotechnology in agriculture. 3. Students can try to improve on the protocol for this experiment OR they may try to develop a new protocol for extracting DNA from another source.
<p>Standards:</p>	<p>Texas Essential Knowledge and Skills (TEKS)</p> <p>§119.68 (c)(2)(A) Identify the role of biotechnology in agriculture. (c)(2)(B) Identify fundamental principles of cell biology and molecular genetics</p> <p>(c)(3)(B) Examine laboratory techniques for manipulating Deoxyribonucleic Acid (DNA) in genetic engineering. (c)(3)(D) Apply scientific measurements and calculations.</p>