



BIOL1414 - Introduction to Biotechnology I

Exercise workbook & Lab Guide



BIOL 1414 – 2018v2

Jack O'Grady, M.S.

ACC Biotechnology Program

BIOL1414 2018v2

Table of Contents

LAB UNIT 1: INTRODUCTION TO BIOL1414.....	3
LAB UNIT 2: ESSENTIAL TOOLS IN THE BIOTECHNOLOGY LABORATORY	25
Lab Unit 2-A: Using a Micropipette	28
Lab Unit 2-B: Calibrating & Using Basic Lab Equipment	34
LAB UNIT 3: PREPARING SOLUTIONS	43
LAB UNIT 4: DNA BARCODING – ISOLATING GENOMIC DNA.....	57
Lab Unit 4-A: Collect, Document, and identify Specimens.....	61
Lab Unit 4-B: Isolate DNA from Plant or Animal tissue.....	63
Lab Unit 4-C: Agarose Gel Electrophoresis of Genomic DNA.....	71
LAB UNIT 5: ENZYME-LINKED IMMUNOSORBENT ASSAY (ELISA).....	79
LAB UNIT 6: RECOMBINANT DNA TECHNOLOGY	93
Lab Unit 6-A: Bioinformatics of Green Fluorescent Protein	93
Lab Unit 6-B: Transformation of <i>E. coli</i> with a Recombinant GFP.....	103
LAB UNIT 7: PLASMID DNA ISOLATION	115
Lab Unit 7-A: Plasmid Mini-preparation.....	120
Lab Unit 7-B: Spectrophotometric Analysis of Plasmid DNA.....	123
LAB UNIT 8: RESTRICTION ENZYME DIGEST.....	127
Lab Unit 8-A: Bioinformatics: Restriction Enzyme Digest of Plasmid DNA	130
Lab Unit 8-B: Restriction Enzyme Digest of Plasmid DNA	131
Lab Unit 8-C: Agarose Gel Electrophoresis of Restriction Digest DNA Fragments.....	133
LAB UNIT 9: PROTEIN PURIFICATION USING COLUMN CHROMATOGRAPHY.....	141
LAB UNIT 10: DENATURING POLYACRYLAMIDE GEL ELECTROPHORESIS.....	151
LAB UNIT 11: DNA BARCODING - AMPLIFY DNA BARCODE USING PCR.....	163
Lab Unit 11-A: PCR Amplification of DNA Barcode	168
Lab Unit 11-B: Analyze PCR by Gel Electrophoresis.....	170
Lab Unit 11-C: Preparing Plates for Sequencing	174
LAB UNIT 12: PCR-BASED VNTR ANALYSIS OF HUMAN DNA	177
Lab Unit 12-A: Isolating DNA & PCR of VNTR Loci	179
Lab Unit 12-B: Analysis of VNTR Loci PCR Amplicons	181
LAB UNIT 13: DNA BARCODING SEQUENCING & BIOINFORMATICS	185
LAB UNIT 14: PREPARING FOR A CAREER IN BIOTECHNOLOGY	201

APPENDIX.....	211
EMPLOYABILITY SKILLS EVALUATION FORM	213
SOLUTION PREPARATION FORM.....	215
ELECTROPHORESIS DOCUMENTATION FORM	217
DEVIATION REPORT FORM.....	219
BIOL1414: LABORATORY REPORT RUBRIC.....	221
SOP TEMPLATE & INSTRUCTIONS	223
INVENTORY CONTROL FORM.....	229

ACKNOWLEDGEMENTS

This manual could not have been possible without the dedication of the faculty, staff and students of Austin Community College’s Biotechnology Program in addition to our valuable educational partners, and numerous grants. Some lab material provided by the following educational partners: Bio-Link, Bio-Rad, Dolan Learning Center, and Edvotek. The following ACC faculty and staff developed some lab exercises: Jack O’Grady, M.S., Linnea Fletcher, Ph.D., Angela Wheeler, M.S., Patricia Phelps, Ph.D., and Evelyn Goss.

Cover art: *Aequorea Victoria*, with permission: http://opencage.info/pics/large_19389.asp



This work is licensed under a [Creative Commons Attribution-Noncommercial-Share Alike 4.0 International License](https://creativecommons.org/licenses/by-nc-sa/4.0/). This publication can be openly shared and copied with attribution. Please attribute Jack O’Grady, Austin Community College. Derivatives welcome, email jack.ogradey@austincc.edu.

FUNDING:

This material developed with the assistance of the following grants:

1. Bio-link award number 9850325 (2014-2016), 1400721 (2013), 0903317 (2010)
2. Texas Education Agency (2004-2015)
3. National Science Foundation, Award Number 0053296 (2004-2006)



LAB UNIT 1: INTRODUCTION TO BIOL1414

Jack O'Grady, M.S., Linnea Fletcher, Ph.D., Angela Wheeler, M.S., & Patricia Phelps, Ph.D.

Welcome to your first course in biotechnology! What an exciting time to learn about biotechnology! Whether you are taking this course for general interest or you want to obtain a job working in a biotechnology lab, this course is designed to provide you with a foundation of biotechnology skills.

The objectives of the lab portion of the course are to:

- ✓ Discover the field of biotechnology, and potential careers in the biosciences
- ✓ Learn basic laboratory techniques of a biotechnology or bioscience lab
- ✓ Develop skills in using basic biosciences laboratory equipment
- ✓ Cultivate critical thinking skills
- ✓ Encourage teamwork, accountability, and taking charge of your learning
- ✓ Practice accuracy in calculations, performing experiments, and in writing scientifically
- ✓ Learn how to work in a regulated biosciences work environment

This workbook is designed to provide you with an active learning discovery experience. Please bring it to every class. It will guide through the course and tie in concepts between the lecture and laboratory. Before every class, you may be asked to read a chapter from your textbook, research information online watch a video or animation or complete homework exercises. It is important to follow your syllabus and schedule, so you know when coursework will be covered and due for grading.

In this first lab unit, you will begin your discovery of the field of biotechnology, explore the laboratory you will be working in this semester, and complete your safety training.

Before you can begin working in an ACC laboratory, you must first

1. View the ACC Science Safety video.
2. Review safety policies, and tour the laboratory with your instructor and locate emergency equipment and procedures.
3. Sign a safety contract, by which you agree to comply with safety regulations.

We hope that you enjoy your experience in this introductory course. The following is a discussion of biotechnology and a description of some of the activities that you will be doing in this course.

LET'S GET STARTED!

Reading: Biotechnology, Science for the new Millennium, Daugherty
Chapter 1: What is Biotechnology?
Chapter 2: Raw Materials of Biotechnology



Part I: What is biotechnology?

Strictly speaking, *biotechnology is the use of a living organism to create a useful product*. By this definition, biotechnology would date back to the very beginnings of civilization, when humankind first learned to cultivate crops and domesticate animals in a system of agriculture. When one thinks of modern biotechnology, however, gene splicing and recombinant organisms take center stage. The Biotechnology revolution came when scientists first learned how to isolate and clone genes, allowing for genetic engineering.

Today, the biotechnology industry has grown and expanded to affect us on a day-to-day basis. Below is an excerpt from Bio.org that discusses the ever-expanding applications of Biotechnology in healthcare, agriculture, and energy (<http://www.bio.org>).

Biotechnology: Healing, Fueling, and Feeding the World (bio.org, 2014)

At its simplest, biotechnology is technology based on biology - biotechnology harnesses cellular and biomolecular processes to develop technologies and products that help improve our lives and the health of our planet. We have used the biological processes of microorganisms for more than 6,000 years to make useful food products, such as bread and cheese, and to preserve dairy products.

Modern biotechnology provides breakthrough products and technologies to combat debilitating and rare diseases, reduce our environmental footprint, feed the hungry, and use cleaner (and less!) energy, and have safer, cleaner and more efficient industrial manufacturing processes. Currently, there are more than 250 biotechnology healthcare products and vaccines available to patients, many for previously untreatable diseases. More than 13.3 million farmers around the world use agricultural biotechnology to increase yields, prevent damage from insects and pests and reduce farming's impact on the environment. And more than 50 biorefineries are being built across North America to test and refine technologies to produce biofuels and chemicals from renewable biomass, which can help reduce greenhouse gas emissions. Recent advances in biotechnology are helping us prepare for and meet society's most pressing challenges.



Let's Explore!

**Watch this! New approaches to
Biotechnology: 3D Bio-Printing video**

<https://youtu.be/peeWHtYsmdM>

In your own words, what is 3D Bio-printing?



Figure 1-1: 3-D Bio-printing a human ear



TEST YOUR
KNOWLEDGE!

WHAT IS BIOTECHNOLOGY?

1. In your own words define Biotechnology:
2. Can you think of a biotechnology product that has improved your life? What do you think your life would be like without this product? What makes it a biotechnology product?

Modern Biotechnology

The modern definition of Biotechnology usually means the manipulation of DNA to create useful products. Deoxyribonucleic acid (DNA) is the carrier of genetic information in our cells. During this course, you will learn to isolate and manipulate DNA from several different organisms including your own!

Let's try it out! Extract your own DNA!

You can easily isolate DNA from your cheek cell using common materials found at home. Follow the instructions here at this website (Nova). Your instructor may have these materials for you in the lab today or may ask you to try it at home!

<http://www.planet-science.com/categories/experiments/biology/2012/03/extract-your-own-dna.aspx>

What you need:

500ml bottled water	2 x <u>clear</u> plastic cups or glasses
Dish soap	1 tbsp. table salt (3 salt packets)
100 ml isopropyl alcohol or ice-cold Ethanol	Food coloring

Steps:

1. Add salt to a bottle of drinking water, recap and shake to dissolve.
2. Transfer 3 tbsp. of the salt water into a clear cup or glass.
3. Gargle the salt water for 1 minute. Don't swallow it!
4. Spit the water back into the cup.
5. Add one drop of detergent to the salt water. Stir gently with a spoon or swirl. Try not to create any bubbles.
6. In a separate cup, mix the alcohol and 1-3 drop(s) food coloring.
7. Gently pour the alcohol and food coloring mixture on top of the salt-water gargle. *Tilt the salt-water cup as you pour, so the alcohol mixture forms a layer on top of the salt water.
8. Wait for 2.5 minutes. You should see white clumps and strings forming. If not, gently swirl, you should see the white strings forming. The white clumps and thread are your DNA!

Observations: Write down your observations of your DNA extraction experiment:

What's happening? When you gargle the salt water and spit it back out, some of your cheek cells become suspended in the salt water. The more vigorously you gargle, the more cheek cells will collect in the salt water. The dishwashing liquid breaks down your cheek cell membranes and causes the release of DNA into the salt water. DNA is not soluble in alcohol, so it forms an aggregate where the alcohol and salt-water layers meet. Most other substances from your cheek cells stay dissolved in the salt-water layer. The white strings and clumps you see are thousands of DNA molecules clumped together. Single DNA molecules are far too small to see with the naked eye. When you gargle the salt water, you are also collecting some bacterial cells from the inside of your mouth, so the DNA you see is a mixture of your DNA and bacterial DNA!

CAREERS IN BIOTECHNOLOGY

Biotechnology Industry in Austin, Texas

The biotechnology industry has also been steadily growing in the Austin area. Today, Austin's bioscience community encompasses over 100 companies that employ more than 7000 people in the areas of research, diagnostics, pharmaceuticals and medical devices (10). Some of these companies include XBiotech, Insys, Agilent, Asuragen, Bioo Scientific, Luminex, and MyriadRBM to name a few! Austin is also a major contributor to academic research in the biological sciences, at the University of Texas, and Texas State University.



Let's Explore!

1. Watch this video to learn more about the exciting and emerging biotech industry in Austin (BioAustin, 2015): <https://youtu.be/vVLfgaDpos>
2. Go to this website and explore careers in biotechnology: <http://biotech-careers.org/>
3. What jobs seem interesting to you?

LABORATORY SKILLS AND COMMON COURSE OBJECTIVES

The following list describes the areas of expertise that we will explore in this course. For complete Common Course Objectives refer to Biotech Program Website:

<http://sites.austincc.edu/biotech/common-course-objectives/>

Basic operations in the laboratory

Students will learn and practice procedures for safe handling and storage of hazardous materials, create and follow detailed protocols, and develop skills in recording data and research notes. Additionally, students will develop solution preparation skills, using basic lab equipment.

Instruments and Equipment

An important part of working in any laboratory is the proper utilization and calibration of instruments and equipment. Students will become learn how to use basic lab equipment, follow SOPs in their operation, and explore equipment validation industry standards; IQ/OQ/PQ. Students will master the use of micropipettes, top loading and analytical balances, electronic pH meter, gel electrophoresis techniques, and a thermocycler.

Working with DNA and proteins

An important part of a lab technician job is to demonstrate basic techniques for purifying and analyzing biomolecules. Students will isolate and analyze DNA from a variety of cells, which include plant, human cheek cells, as well as transform *E. coli* with a recombinant plasmid. Additionally, students will also learn basic methods to purify and analyze recombinant proteins.

Immunochemistry

Students will be introduced to basic techniques used to detect biomolecules using antibodies and will perform the detection and quantitation of a target protein using ELISA.

Regulatory Affairs

Students will write and follow Standard Operating Procedures (SOPs), competently complete associated forms that are commonly used in a biomanufacturing facility. Additionally, students will learn regulations that govern biological laboratories, demonstrate safety procedures and protocols for disposal of hazardous chemicals and biologicals, and will perform equipment validation (IQ/OQ/PQ) on a variety of basic laboratory equipment.

Bioinformatics

Students will use computers to document and compile information to analyze research data, generate graphs, and utilize protein and genomic databases, such as NCBI, BLAST, NEBcutter, and DNAsubway.

Part II: Documentation: The Lab Notebook & Lab Exercise Workbook

Introduction: Documentation in a lab notebook is an essential skill for any biotechnician. The Food and Drug Administration's (FDA) philosophy is, "if it isn't written down, it wasn't done." Documentation details vary from lab to lab, but it is essential for one or more of the following:

- ✓ to record what was done and observed
- ✓ to establish ownership of a patent and other legal uses
- ✓ to establish criteria used to evaluate a finished product or the process to make it
- ✓ to trace the manufacture of a product
- ✓ to create a contract between a company and consumers or between a company and regulatory agencies
- ✓ to prove that a procedure was performed correctly
- ✓ to adhere to, evaluate, and develop standard operating procedures (SOP)

Lab Notebook:

Each student will maintain a 3-ring binder lab notebook to keep your workbook and laboratory exercises, notes and data throughout the semester. *Your lab notebook is graded.*

Laboratory Notebook Contents: In most workplaces, a hard-bound laboratory notebook is required. However, in this class, you are required to take detailed lab notes in your workbook and keep them in a binder as your lab notebook. Bring this binder to every class, and you will submit your lab notebook periodically during the semester for grading. Your instructor will provide you with a grading rubric for your laboratory notebook.

For this class, your lab notebook will be a 3-ring binder which must include:

1. Title Page: Your name, course name
2. Table of Contents: The notebook should be in chronological order (lab 1, 2, 3...) and tabs to easily find each lab.
3. For each lab (in this order):
 - Completed workbook exercises for that lab
 - Any notes, and raw data recorded during the exercise
 - Any forms generated during lab
 - Any handouts from lab
 - Any data analysis (graphs, tables, charts, etc.)
4. You are encouraged to keep any lecture notes in your binder with corresponding workbook exercise.

CLASS FORMAT

BEFORE CLASS

1. Look at the schedule to see which exercise in the workbook you will be covering in class and on which date.
2. Read the assigned chapters from your textbook, and review lecture material.
3. Complete any assigned exercises before coming to class – this is critical to making the most of the lab time. If you do not prepare before the lab, you may not have enough time to complete the exercise or may make mistakes, and the lab exercises may not succeed.

Most importantly, it won't be very much fun if you don't know what you'll be doing! 😊

DURING CLASS

During lab experiments, you will take notes in pen. WRITE EVERYTHING DOWN. Yes, we mean everything. How much did you weigh out? What are the supplier and the lot number of the reagent? What balance number did you use? What color was your solution? And so on. Be sure to include any changes you made to the procedure in the lab handout, even if they were at the instructor's direction. Always show calculations.

Writing down everything improves your observational skills, helps you understand the importance of each step, and provides a record of how an experiment might have gone wrong. Each student should record notes, even when working in teams.

General guidelines for writing good lab notebooks:

- Write all parts of your lab notes and observations in ink. *Writing with a pencil is forbidden in the lab*. If you make an error, draw a single line through it and enter your correction in clear and legible writing. If you discard data for any reason, you must justify your decision to do so immediately and in writing.
- *Write legibly*. Remember, coworkers, supervisors, and perhaps, the FDA may be reading your notebook. If they cannot read your writing, your work is essentially nonexistent, or they may misinterpret an important detail about your work.
- If you tape materials such as a graph, a manufacturer's specification sheet, or instrument readout into your notebook or onto a form, tape all four sides. Then write "NWUI" ("No writing under insert") on the tape, your initials, and the date.
- Keep your records factual, concise, clear and complete in all aspects.

POST-LAB

One week after all the data has been collected you must submit your workbook exercise for grading. Make sure your workbook is easy to read and well organized. Include:

1. A title page
2. Complete workbook exercises and handouts
3. All lab experiment notes, and raw data
4. Include any forms filled out during the exercise
5. Attach graphs, tables and other ways in which your data was manipulated

Late Workbooks: Part of the skills you are learning in this class is to generate quality work while meeting expected deadlines. Late assignments will be subject to the policy outlined in the syllabus. See your instructor for further information.

Missed Classes: You are expected to attend every class. Since we will be performing graded work every class, missed classes will have grade consequences. The penalty for missed labs is outlined in your syllabus.

Lab Competency

Exploring and expanding your understanding of Biotechnology is the most important outcome of this course. However, if you're interested in having a career in this field competency in laboratory skills is even more important. Note that skill competency is not limited to lab skills, but also includes attendance, punctuality, teamwork, and tidiness. Lab skill is monitored using lab exams, Employability Skills Evaluations, and lab notebook grading.

Part III: Biotechnology Laboratory Security & Safety

OBJECTIVES

Your performance will be satisfactory when you can:

- ✓ Discuss security and safety rules for the laboratory
- ✓ Recognize the correct procedure for storing and handling hazardous materials
- ✓ Locate the lab safety equipment
- ✓ Find online Material Safety Data Sheet (MSDS) databases and use them to search for MSDS sheets for chemicals used in the lab.
- ✓ Identify potential hazards with chemicals in the lab, and how to use them safely.

INTRODUCTION TO LAB SECURITY

Biotechnology laboratories are equipped with supplies and equipment that may pose a hazard if used carelessly and it is important that you learn how to handle them properly. It is often the responsibility of a biotechnician to make sure that safety rules are followed, and anyone working in a laboratory must pay attention to what they are doing and use common sense to avoid hazardous situations.

While the ACC science safety rules are designed to protect you while working in ACC laboratories, you must become self-sufficient in protecting yourself in your future jobs in the biotechnology industry. Also, lab technicians are entrusted with ensuring compliance with safety precautions in the biotechnology workplace. For this purpose, this lab exercise will introduce you to the main components of lab safety measures and procedures that apply in a biotechnology setting.

LABORATORY SECURITY

Educational institutions and biotechnology companies use a wide assortment of highly hazardous materials. When working with these materials every day, it is easy to forget about the harm these substances can cause if they are stolen. Following the terrorist attacks of September 2001 and the "anthrax letters" sent the same month, much attention has been directed to practical measures that will keep hazardous materials (biological and chemical) out of the hands of criminals. Many new federal laws were enacted in direct response to these terrorist attacks. It is important that laboratory personnel take specific actions to prevent unauthorized entry to labs, secure highly hazardous materials against theft, and ensure compliance with new security regulations. It is essential to implement procedures necessary to provide security of all hazardous materials in their areas of responsibility. One objective is to minimize the risk of theft, especially during that five-minute window when the lab is left unattended. One easy way to increase security is to make sure that your laboratory door is locked whenever the lab is left unattended, even for a few minutes. Having multiple locked door layers, such as in our laboratory where the chemicals are locked away in a preparation room is very practical in avoiding theft of hazardous material.

Different laboratories implement various security measures, which include locking up controlled substances, balances, computers, equipment and syringes and needles. Laboratory personnel should review and assess the security of their highly hazardous materials, such as infectious agents, toxins, radioactive materials, acutely toxic chemicals, carcinogens, explosive or reactive chemicals, and compressed gases. The following guidelines were adapted from Appendix F of the CDC/NIH publication, Biosafety in Microbiological and Biomedical Laboratories. The guidelines are intended to reduce the risk of unauthorized removal of hazardous materials from your laboratory:

1. **Recognize that lab security is related to but different from laboratory safety and develop a site-specific security policy.** Security, as used in this discussion, refers to measures used to control access to the laboratory to prevent theft of materials or equipment from the lab.
 - o Assess your laboratory for hazardous materials and particular security risks.
 - o Develop and implement lab security procedures for your lab group.
 - o Train your lab group on these security procedures and assign responsibilities.
2. **Control access to areas where hazardous materials are used and stored.**
 - o Close and lock laboratory doors when no one is present. Consider the use of card-keys or similar devices when the risk warrants.
 - o Do not leave hazardous materials unattended or unsecured at any time.
 - o Lock freezers, refrigerators, storage cabinets, and other equipment where biological agents, hazardous chemicals, or radioactive materials are stored when they are not in use.
3. **Know who is in your laboratory area.**
 - o Consider using a logbook for staff to sign in and out of the lab each day or using carded access devices for this purpose.
 - o Limit laboratory access to those individuals who need to be in the lab.

- All lab workers (including students, visiting scientists and other short-term workers) should wear identification badges.
 - Restrict off-hours access to individuals authorized by the principal investigator.
 - Guests should be issued badges and escorted to and from the lab. Approach people you don't recognize who appear to be wandering in laboratory areas and ask if you can help direct them.
4. **Know what materials are brought into your lab.**
 - Know what hazardous materials are being ordered and shipped to your lab.
 - Get rid of unneeded hazardous materials.
 - Use a log to sign highly hazardous materials in and out of secure storage.
 - Take periodic inventory of all highly hazardous chemicals, biological agents/toxins, radioactive materials, and controlled substances.
 5. **Know what materials are removed from your lab.**
 - Track the use and disposal of hazardous materials.
 - Require written permission before removal of highly hazardous materials from the lab.
 - Report any missing inventory.
 6. **Have an emergency plan.**
 - Recognize that controlling access can make emergency response more difficult.
 - Evaluate emergency plans with administrators, safety and security officials and, if necessary, outside experts.
 - Review emergency plans with lab personnel.
 - Provide emergency responders with information on serious hazards.
 7. **Have a protocol for reporting security incidents.**
 - Principal investigators, in cooperation with facility safety and security officials, should have policies and procedures in place for the reporting and investigation of incidents or possible incidents, such as undocumented visitors, missing hazardous materials, or unusual or threatening phone calls.
 - Train laboratory staff on procedures.

PROPER HANDLING & STORAGE OF CHEMICALS AND REAGENTS

There is no single simple formula for working safely in the laboratory since each lab facility and each experiment presents unique challenges. We will be addressing safety issues with each experiment that we do in this course and give you some specific guidelines for safety throughout the semester.

A. MSDS (Material Safety Data Sheets)

While each chemical that you use will have unique properties, some practices will aid in treating all chemicals with the level of respect that they are due. For example, labeling each chemical is required under the law and should be thorough enough so that even a person who does not work in the lab can identify any chemical. Also, every chemical in the laboratory should have a **Material Safety Data Sheet** (MSDS) or Safety Data Sheet (SDS) on file and readily available. The MSDS is a legally required technical document, provided by chemical suppliers, that describes

the specific properties of a chemical. Besides the MSDS on file in the lab, several websites offer MSDS databases. They are all broken down to the same eight sections:

1. **Chemical identity.** The manufacturer's contact information is here, along with contacts for emergency situations.
2. **Hazard ingredients/identity.** Some reagents have multiple components, and many single-component chemicals have alternative names. These are all listed here. Concentration limits for airborne exposure to a chemical are listed here. Although these indices of toxicity are mainly of concern for production workers in factories, they are also useful for evaluation of short-term exposures.
3. **Physical-chemical characteristics.** This list of physical properties tells you whether the chemical is solid or liquid and how volatile it is.
4. **Fire and explosion hazard data.** This is of particular interest in cases where fire-fighting methods must be selected.
5. **Reactivity data.** This information is essential in determining the proper handling and storage of chemicals. By knowing the reactivity patterns of a chemical, you know what substances or conditions from which you must isolate the chemical.
6. **Health hazards.** The best source of specific toxicology data is given here, such as symptoms of acute damage from exposure and some recommended emergency procedures. If a chemical has been tested for **carcinogenicity** (cancer-causing potential), that information is listed here.
7. **Precautions for safe handling and use.** This describes how to deal with spills.
8. **Control measures.** Specific recommendations for personal protective equipment (PPE) are given here.

B. NFPA Ratings (National Fire Protection Association)

Another quick assessment of a chemical's health hazards that is usually available on its container is a rating by the National Fire Protection Association (NFPA). To learn more, see their website: <http://www.nfpa.org/codes-and-standards>

A color-coded diamond shape lists numbers are rating a hazard as:

Blue for health hazard	Red for flammability	Yellow for reactivity
0 – normal material	0 – will not burn	0 – stable
1 – slightly hazardous	1 – flash point > 200° F	1 – unstable if heated
2 – hazardous	2 – flash point > 100° F	2 – violent chemical change
3 – extreme danger	3 – flash point < 100° F	3 – shock and heat may detonate
4 – deadly	4 – flash point < 73° F	4 – may detonate

The **uncolored station** of the NFPA diamond is for specific hazards:

OX – oxidizer compound, **ACID** – acidic compound, **ALK** – basic compound, **CORR** – corrosive compound, and **W** – use NO WATER

C. General Safety Precautions in Handling Hazardous Chemicals in the Lab

There are four routes to exposure to hazardous chemicals that you should keep in mind while handling them:

1. **Inhalation:** avoid by the use of fume hoods and masks
2. **Skin & eye contact:** avoid by the use of lab coats, gloves, and goggles
3. **Ingestion:** avoid eating or drinking in the lab, or leaving the lab without removing gloves and washing hands
4. **Injection:** dispose of broken glass and needles properly

General Safe Handling Procedures for Chemicals:

- ✓ Treat all chemicals as if they were hazardous until you learn otherwise
- ✓ Label all containers with contents, including concentrations and date
- ✓ If a hazardous material is contained, label it with a warning
- ✓ Think through your experiment ensure you will not be combining incompatible chemicals
- ✓ Clean your bench top before and after use
- ✓ Wash hands, often and ALWAYS before leaving the lab
- ✓ Take off lab coats and gloves before leaving the lab
- ✓ Always remove gloves before touching phones, doorknobs, light switches, etc.
- ✓ Ensure proper waste disposal and labeling.

Here are some specific tips for handling the different types of hazardous chemicals:

- ◆ **Flammables:** Do NOT heat unnecessarily, and never in the presence of a flame or source of a spark. In general, only open containers in fume hoods. When storing more than 10 gallons of flammable liquids, a special explosion-proof storage cabinet is required.
- ◆ **Corrosives:** Wear **personal protective equipment (PPE)** such as lab coats, goggles and gloves, and always add strong acids or bases to water when making solutions. Neutralize slowly to avoid rapid generation of heat and gases. Strong acids and bases should never be stored together.
- ◆ **Reactive chemicals:** Wear PPE such as lab coats, goggles, and gloves, and know the reactive properties of the chemical. Always store oxidizing chemicals away from flammable materials.
- ◆ **Toxic chemicals:** Wear PPE such as lab coats, goggles, and gloves, and know the toxic properties of the chemical. When working with a dry powder, wear a mask to avoid breathing the dust. Be aware of the waste disposal procedures.

Here are some of the most common hazardous chemicals that you will encounter in the biotechnology lab:

Carcinogens	– formaldehyde	Mutagens	– ethidium bromide
Neurotoxins	– acrylamide	Teratogens	– formamide
Nephrotoxins	– acetonitrile	Hepatotoxins	– chloroform
Corrosives	– phenol, strong acids & bases		



Let's Explore

There will be many potentially hazardous chemicals you will come into contact with in this course and on your job if you move into this field. MSDS are readily available on the internet. You can reference many chemical manufacturer websites such as Fischer Scientific and Sigma. Or try <http://www.msds.com/>. Let's take a look at some of the MSDS sheets for chemicals you may come in contact.

1. Learn more about MSDS here: <http://www.ilpi.com/msds/> You may want to try these: learn what an (M)SDS is, what one looks like, how to read one and where to find one.
2. Look up Sodium Chloride (NaCl) and Hydrochloric Acid (HCl). Summarize the primary hazard (if any) and what to do if accidentally exposed to these chemicals

NaCl:

HCl:

D. BIOLOGICAL SAFETY: CONTAINMENT

Biohazards. Students will be working with live organisms in many biotechnology labs, so it is important to be able to assess any biological hazards that they may pose and to treat them accordingly. The Biology department organism use policy can be found here:

<http://sites.austincc.edu/biology/lab-animal-policy/>

The routes of exposure to infectious agents are the same as those of hazardous chemicals: inhalation, contact with eyes and skin, ingestion, and injection. The same general precautions should be taken in handling biological hazards as the guidelines above for managing chemical hazards, especially toxic ones. Here are some general practices to maximize biological safety:

- ◆ Limit access to the lab and adequately train all lab personnel.
- ◆ Use personal protective equipment (PPE) at all times, and keep all PPE inside the lab.
- ◆ Wash hands after removing gloves and before leaving the lab.
- ◆ Always remove gloves before touching phones, doorknobs, light switches, etc.
- ◆ Avoid touching your face with your hands or gloves.
- ◆ Keep personal items such as coats and bags out of the lab.
- ◆ Minimize splashes and aerosol production.
- ◆ Disinfect work surfaces to decontaminate after a spill and after working in the lab.
- ◆ Decontaminate all regulated waste before disposal (usually by autoclaving).
- ◆ Have an insect and rodent control program in effect.
- ◆ Use a laminar flow biological safety cabinet when available.

A majority of recorded laboratory-acquired infections are due to inhalation of infectious particles, so special precautions should be taken to avoid producing aerosols when working with pathogens. While performing activities that mechanically disturb a liquid or powder, the biotechnologist should make the following adjustments.

<u>Activity</u>	<u>Adjustment</u>
• Shaking or mixing liquids	mix only in closed containers
◆ Pouring liquids	pour liquids slowly
◆ Pipetting liquids	use only cotton plugged pipettes
◆ Removing a cap from a tube	point tubes away when opening
◆ Breaking cells by sonication in the open	sonicate in closed containers
◆ Removing a stopper or cotton plug	remove slowly
◆ Centrifuging samples	use tubes with screw cap lids
◆ Probing a culture with a hot loop	cool loop first

Disinfectants such as bleach and ethanol are used extensively to decontaminate glassware and work areas, and it is important to realize that the effectiveness of disinfectants depends on the type of living microorganisms you are encountering:

<u>Resistance Level</u>	<u>Type of Organism</u>	<u>Examples</u>
Least resistant	enveloped viruses	HIV, Herpes simplex, Hepatitis B
Slightly resistant	bacteria	<i>E. coli</i> , <i>S. aureus</i>
Medium resistance	fungi	<i>Candida</i> species, <i>Cryptococcus</i>
Highly resistant	non-enveloped viruses	Poliovirus, Mycobacteria, <i>M. tuberculosis</i>
Most resistant	spore	<i>B. subtilis</i> , <i>Clostridium</i> species

E. DISPOSAL OF HAZARDOUS CHEMICALS & BIOLOGICAL MATERIALS

The disposal of hazardous chemicals is subject to state and federal regulations and is ultimately overseen by the Environmental Protection Agency. Highly toxic chemicals are regulated at low levels, and less toxic chemicals can be disposed of through city sewer systems at higher levels. Biological hazards should be contained in autoclave bags made of a high melting point plastic that is sealed and autoclaved at high temperatures and pressures to kill any live organisms completely.

In our laboratory, specific hazardous chemical and biological waste disposal will be discussed at the start of every lab. Chemical waste disposal containers can be found in the fume hood. Always keep the fume hood on and the lids to the chemical waste disposal containers on the containers.

Safety Assignment

1. Watch the safety video as directed by your instructor. It is also available online here: <https://youtu.be/uDhSSWtY3kg>
2. Your instructor will provide you with a safety training sheet and will go through this with you as a class to provide lab-specific safety training. ***Write down all the information in the blanks and keep it available in your lab notebook while working in the ACC laboratory.***
3. You must sign an ACC Safety Contract to acknowledge that you understand the safety rules and agree to abide by them.
3. Each lab group will perform a safety check today by filling in the Inspection sheet on the next page.
4. Each week a student will be an assigned safety officer. The safety officer will perform the following safety checklist. Your instructor will assign groups weekly and provide the groups log sheets. Teams will give a 1 min summary to the class of any safety violations they found. http://www.austincc.edu/sci_safe/docs/weekly_safety_inspection_log.pdf
At the end of the last class, they will give a summary of any safety issues that came up in the safety evaluation and provide the class with the checklists to include with their lab report.
5. Next week's safety officer(s) is: _____

Semester Safety Inspection Log

Campus: _____ Room number: _____

Safety Inspector(s): _____

Initial each item if acceptable or comment on necessary improvements. Immediately report problems with any of these items

- Aprons/lab coats clean and usable?
- Fire blankets intact and usable?
- First aid kits stocked?
- Spill kits accessible?
- Goggles & safety glasses used in the lab?
- Safety posters & signs readable?
- Eye-wash station and shower unimpeded and functional?
- "Notice to Employees" present?

Comments:

Date: _____

Signature: _____ printed name: _____

Signature: _____ printed name: _____



Part IV: Lab Organization & Equipment Identification

During this course, you will learn to use, calibrate and troubleshoot many pieces of equipment used in biotechnology labs, and you will be making a variety of reagents. Before you get started, it is important to learn basic functioning of the lab. This helps with efficiency, cleanliness, and keeping a safe work environment. The ACC Biotechnology program regards lab etiquette as an important part of the curriculum. Showing courtesy to students, staff, and instructors who share the work area by caring for equipment, leaving a clean workspace, and removing biological and chemical hazards are considered practicing good lab etiquette.

Each student or group should perform the following before leaving lab after every class:

1. Ensure that any solutions you have made are labeled appropriately according to SOP SOL-001 and that you have created a solution prep form for each. Store them properly in the provided storage location for your class.
2. Replace any equipment, supplies, or reagents that you have used to their proper storage place, provided that other students are not still using them.
3. Clean your work area, which includes removing all items from your lab bench and wiping it with a paper towel wet with bench cleaner.
4. Rinse out dirty glassware before putting where directed. Do not use a sharpie on equipment.
5. Return equipment and reagents to cart, and move rolling cart(s) into the prep room.
6. Turn off all equipment used during the lab by the class. Make sure the scales and scale area are CLEAN.
7. You are required to assist other students who are still working before you leave the lab.

Broken Glass: Please use the blue broken glass boxes for broken glass disposal (used slides and cover slips, Pasteur pipettes, broken glassware). There are a dustpan and broom available to use for sweeping up broken glass. ONLY GLASS GOES IN THE GLASS WASTE – NO PIPETS, PAPER TOWELS, OR USED GLOVES.

Equipment:

- ✓ An equipment locator can be found in the designated file cabinet in the lab room. Use this document to locate supplies and equipment.
- ✓ SOPs written specifically for the equipment in our department are found in the SOP packet that we have provided each semester. You will receive one hard copy at the beginning of the semester for all your classes. Please KEEP this SOP booklet with you.
- ✓ Turn off all equipment before you leave
- ✓ If any equipment is not functioning properly or appears damaged, fill out a deviation report and notify your instructor. Do not return broken equipment this is a safety hazard.



SAFETY & LAB EQUIPMENT ORIENTATION

1. With your lab partner, explore the Biotechnology Department Laboratories and Prep Room.
2. Using the Equipment Locator Key try to find the following safety-related materials as well as the equipment and materials that you will be using throughout the semester.
3. If you don't know what a piece of equipment looks like, it may be hard to find. Go on the internet and look it up! Try Google Images!
4. If you still can't find something, ask another group. And if you still can't find it, ask your instructor.
5. Turn these completed in sheets with your lab report.

Safety Related Materials/equipment	Room stored in	Location
Eye Wash Stations (all of them!)		
Fire Extinguishers (all of them!)		
Fire Blankets		
Emergency Gas shut off Valve		
Large & Small Glass Waste		
Biohazard Waste		
Liquid Chemical Waste		
General Chemicals		
Spill Kit		
Broom/dustpan		
Material Safety Data Sheets		
First Aid Kits		

Laboratory Materials/equipment	Room stored in	Location
Micropipettes		
Micropipette tips		
1.5mL microcentrifuge tubes		
Picofuge		
Microcentrifuge		
Electrophoresis power supply		
Horizontal Electrophoresis chambers		
Parafilm		
Weigh boats		
Standard pH buffers		
Unopened glove boxes		
Test tube racks		
Graduated cylinders		
Erlenmeyer Flasks		
Freezers (-20°C)		
Refrigerator (4°C)		
Shaker Incubator (37°C)		
NanoDrop Spectrophotometer		
Hotplate/stir plate		
Stir bars		

LABORATORY UNIT 1 ASSIGNMENT

1. **Obtain 3-ring Binder:** This is your lab notebook for this course and will contain the lab manual, pre-lab, and post-lab reports. We recommend a 2-3 inch binder.

For your workbook 1 submit your binder for grading which must include:

- a lab notebook title page,
 - table of contents,
 - the lab manual,
 - and staple lab unit 1 together and submit with your binder.
 - Lab Unit 1 must include:
 - ✓ title page
 - ✓ workbook
 - ✓ completed safety worksheet
 - ✓ completed safety inspection log
 - ✓ completed equipment orientation worksheet
2. **Obtain a sharpie:** You will need a permanent marker to use in the lab.
 3. **Obtain safety equipment:** You will need safety glasses, or goggles rated Z87 (or Z87.1) and closed-toed shoes.
 4. **Complete all exercises in unit 1.**
 5. **Watch Safety Video:** The ACC Safety Committee has produced a video explaining safety rules and regulations. You must watch this video. We will do this in class together. It is also available online here: <https://youtu.be/uDhSSWtY3kg>
 - a. **Sign Safety Contract:** You must sign an ACC Safety Contract to acknowledge that you understand the safety rules and agree to abide by them.
 - b. **Fill-in Safety Worksheet:** Your instructor will provide you with a safety training sheet and will go through this with you as a class. This will provide room-specific safety training. ***Write down all the information in the blanks and keep it available in your lab notebook while working in the ACC laboratory.***
 - c. **Complete safety Inspection Log:** Include this in your workbook.
 - d. **Fill-in Orientation Sheets:** Using the Equipment Locator Key provided, explore the laboratory with your lab partner. Fill in the attached Orientation Sheets. ***if there is equipment you don't know, look up pictures on the internet!*** Include these sheets with your report.



LAB UNIT 2: ESSENTIAL TOOLS IN THE BIOTECHNOLOGY LABORATORY

Jack O'Grady, M.S., Linnea Fletcher, Ph.D., Angela Wheeler, M.S., Patricia Phelps, Ph.D.

OBJECTIVES

Your performance will be satisfactory when you can:

- ✓ Identify essential lab equipment and describe their function
- ✓ Choose appropriate glassware and equipment for accurate measurement
- ✓ Understand the role of the reagents you use in the laboratory
- ✓ Correctly operate a micropipette
- ✓ Demonstrate how to calibrate and use a pH meter and top loading scale

During your training in the ACC Biotechnology program, you will learn to use, calibrate and troubleshoot many pieces of equipment used in biotechnology labs, and you will be making a variety of reagents. You are required to keep a list of the equipment that you learn to use and a brief description of the purpose of the machine. For example, a PCR machine is used to amplify a particular section of DNA.

Concerning the equipment, to use it you need to know its location in the laboratory. Please locate the following items in the lab. Your lab instructor or lab tech may provide a guided tour of the lab pointing out all the equipment listed. Alternatively, your instructor may provide you with an equipment locator key.

LET'S GET STARTED!

Reading:

Chapter 3.1: The Basic Skills of the Biotechnology Workplace!

Measurement of Volume. There are many different types of glassware used in a biotechnology lab. The glassware used will depend on the purpose of its usage as different glassware have various levels of accuracy of measurement! In the next few labs, students will learn how to choose the correct glassware for the appropriate situation.

Erlenmeyer flasks are primarily for the preparation or storage of solutions, not an accurate volume adjustment. Although there are volume markings on these flasks, they are not calibrated and should not be relied upon for exact volume measurements.

Beakers are used for preparing solutions, such as dissolving powdered reagents in water. Beakers are useful when a pH adjustment requires access to the solution by a pH probe. The volumetric markings on beakers are not reliable – and is not an accurate measuring tool.



Graduated cylinders are calibrated with sufficient accuracy for most volume measurements when preparing solutions. For example, the calibration of most 100 mL graduated cylinder can be relied upon to measure to within ± 0.6 mL accurately. Graduated cylinders are most frequently used to bring solutions to a final volume.

Volumetric flasks are used to measure one particular volume with the highest degree of accuracy and are used to make standard solutions for analytical assays. For example, the calibration of a 100 mL volumetric flask can have an accuracy of ± 0.1 mL. Note, there are no graduated markings in a volumetric flask – you can only bring solutions to one final volume in that flask.

Pipets are glass or plastic devices that are routinely used to measure and transfer liquids by drawing the liquid into the tube with a bulb or mechanical pump.

- **Pasteur pipets** are small glass tubes used with a **bulb** to transfer volumes as little as a single drop and as large as a few milliliters. They are not graduated and are not used to measure volumes.
- **Beral pipets (transfer pipets)** are plastic pipettes with a bulb at one end used for transfer of liquids. Sometimes they have calibration marks, which have a low level of accuracy. They are often disposable, sterile and individually wrapped.
- **Serological pipets** are graduated tubes used to measure anywhere from 1 to 50 mL. When the liquid has drained from this pipet, the final drop in the tip is transferred with a puff of air. These are known as TC or to contain pipets. Serological pipets are most frequently used in bioscience labs with pipet controller (or pipet-aid), mechanical devices that aspirate liquid into and out of the pipet.
- **Mohr, or “to deliver,” pipets** are similar to blowout pipets, but do not require a puff of air to deliver the desired volume accurately. They can be identified by the label “TD” on the top. These are sometimes called TD serological pipets.
- **Volumetric pipets** are not graduated but are calibrated to deliver a single, highly accurate, volume.
- **Micropipettes** are mechanical devices with disposable plastic tips, which deliver with a high degree of accuracy adjustable microliter volumes of liquid. There are several micropipette devices available of varying sizes, such as a 0.5-10ul, 2-20ul, 20-200ul, and 100-1000ul.
- **Repeater Pipettes** are mechanical devices that can be set to deliver, repeatedly, a precise microliter volume. The liquid is aspirated and dispensed from a disposable tip.
- **Multichannel micropipettes** can deliver the same volume from as many as 12 tips simultaneously.

Measurement of Weight. Instruments for weighing materials are called balances, and most laboratories have more than one type of balance, depending on the amount of material being measured and the degree of accuracy required.

- **Electronic balances** usually have a digital readout, and weighing dishes can be **tared** to read zero mass before using. Most balances used for the preparation of solutions have a sensitivity of ± 0.01 g, (“**top loading balance**”) but **analytical balances** can be sensitive to ± 0.01 mg or less. Electronic balances require routine maintenance and recalibration.

Measurement of pH. Most solutions prepared in the biological laboratory must have a carefully controlled pH. Buffers are prepared by adjustment to a specific pH with strong acid and base solutions, using a meter to monitor the pH. A **pH meter** is a voltmeter that measures the electrical potential between two electrodes. One electrode is in contact with your solution, and the other is in contact with a reference solution. Usually, both of these electrodes are combined in a single pH probe that you place in your solution. These meters can read to the nearest 0.1 pH unit, but require frequent calibration with reference buffers of known pH.

Solution Preparation. Solution preparation involves mixing liquids and dissolving solids in liquids. There are many specialized devices in addition to balances, volume measuring devices, and pH meters engaged in these processes.

- **Magnetic stirrers** come in the form of a box with a magnet inside attached to a motor that spins the magnet. When a vessel containing a magnetic stir bar is on top of the magnetic stirrer, the stir bar rotates and stirs the contents of the container.
- A **vortex mixer** rotates the bottom of a tube rapidly; setting up a vortex in the liquid that quickly mixes the contents.

Centrifugation. Many pieces of equipment are used to centrifuge biological samples for analysis.

- A **preparative centrifuge** has a balanced rotor that holds vessels and spins at high speed, up to 20,000 rpm. High-speed centrifugation will cause insoluble particles such as cells, and sometimes subcellular components, to form a pellet at the bottom of the vessel. Rotors are available that hold vessels as small as a few milliliters to as large as a liter. These centrifuges are often refrigerated so that heat-sensitive material, such as cells and proteins, are not damaged due to the high heat generated during centrifugation.
- A **tabletop or clinical, the centrifuge** is not refrigerated and spins at a much slower speed than a preparative centrifuge. Rotors for clinical centrifuges hold tubes with a capacity of 15 mL or less and a 'swing-bucket' rotor. These are frequently found in medical laboratory testing labs.
- A **microcentrifuge** holds 1.5mL microcentrifuge tubes with liquid and can centrifuge at high speeds to separate liquids and particulates in solutions. There are a variety of rotor sizes available, and there are refrigerated centrifuges available as well.
- A **picofuge** is a fixed low-speed microcentrifuge, which spins much slower than a microcentrifuge and functions mostly to move liquid from the sides and top of a microcentrifuge tube to the bottom.

Lab Unit 2-A: Using a Micropipette

OBJECTIVES

Your performance will be satisfactory when you can:

- ✓ Properly handle and operate a micropipette
- ✓ Correctly select the appropriate micropipette to measure a particular volume
- ✓ Successfully and confidently operate a microcentrifuge
- ✓ Learn and practice Good Laboratory Practices (GLP)

Purpose: The purpose of this exercise is to become familiar with some of the essential tools of the Biotechnician. In this lab, you will learn how to handle and operate a micropipette properly, learn to select the appropriate micropipette to measure a particular volume correctly.

Safety: In this laboratory exercise we will not be using any hazardous materials, operations or dangerous equipment. However, it's part of Good Laboratory Practice to always wear close-toed shoes, lab coat, safety glasses and gloves when working in a laboratory.

Follow these **Good Laboratory Practices (GLP)** and make them a habit for every lab:

1. Keep your work area clear of unnecessary items
2. Keep everything you need within reach
3. Gather all materials and set up disposal before you begin working
4. Label each container BEFORE you fill it
5. Change gloves often to avoid contamination and never wear your gloves out of the lab
6. Never perform protocols from memory; always read every step every time you perform a procedure, and then check it off as it is completed
7. Always cap bottles of stock solutions and chemicals when finished

Introduction to the micropipette: The micropipette is one of the biotechnician's most frequently used tools. There are different brands of micropipette such as Brinkmann, Labsystems or Rainin. Each of the brands has their devotees, but all of them work in the same way.

All micropipettes are essentially long tubes with a handle and an adjustable piston inside. A disposable tip is placed on the bottom of the micropipette or pipette or micropipet (you will see it written many different ways). This tip is the only piece that is inserted into the liquid. In or near the handle is a screw/knob/button that adjusts the volume of the micropipette by moving the piston up and down. On top, there is a plunger button for filling and for dispensing the liquid. There is often a second button on top for ejecting the tip. Micropipettes come in a variety of sizes; for example, 1-20 uL, 20-200 uL, and 100-1000 uL. When choosing which micropipette to use, the rule of thumb is to *select the smallest size pipette that can deliver the desired volume.*

A biotechnician is usually issued a personal set of the micropipette and is responsible for cleaning and verifying them on a regular basis. During this semester you and your lab partner will be assigned a set of pipettes. It is your responsibility to care for and maintain them.





Let's Explore!

Watch this video! https://youtu.be/uEy_NGDfo_8

Describe below how to use a micropipette to measure 200ul of liquid accurately. Be specific, and include tips you learned in the video.

Your turn! Let's use micropipettes!

Materials

Per Pair:	Per person:
<ul style="list-style-type: none"><input type="checkbox"/> 50 ml conical containing sugar solution<input type="checkbox"/> 50 ml conical containing deionized water<input type="checkbox"/> 50 ml conical containing colored water<input type="checkbox"/> Set of 3 micropipettes<input type="checkbox"/> Box of 0.5-10 μL tips<input type="checkbox"/> Box of 20-200 μL tips<input type="checkbox"/> Box of 100 – 1000 μL tips<input type="checkbox"/> Picofuge<input type="checkbox"/> Plastic tip waste beaker	<ul style="list-style-type: none"><input type="checkbox"/> 1.5 mL microcentrifuge tubes (13)<input type="checkbox"/> Microcentrifuge tube rack

PROCEDURE

Part I: Organizing Your Work Space.

1. Before you begin working clean the bench top with the cleaner provided by spraying a paper towel and wiping the bench with the wetted towel.
2. Collect everything you will need for the lab. Check off each material as you retrieve it and return it to your bench.
3. Note, each person in your group will perform each of the measurements. Please obtain your own microcentrifuge tube rack and microcentrifuge tubes.
4. Locate the set of micropipettes you and your lab partner will be using. Ensure there is a complete set, with the same set numbers on them. If you do not have a complete numbered set, let your instructor know.

Part II: Familiarizing yourself with the micropipette

1. Components of a Micropipette

Examine the set of micropipettes for your lab group. Notice the volume ranges for the set you have (top of micropipette). Familiarize yourself with the components of the micropipette.

 <p>PLUNGER BUTTON</p> <p>TIP EJECTOR BUTTON</p> <p>VOLUME ADJUSTMENT KNOB</p> <p>DIGITAL VOLUME INDICATOR</p> <p>PLASTIC SHAFT</p> <p>TIP EJECTOR ARM</p> <p>DISPOSABLE TIPS</p>	  <p>This is a p20 micropipette. What volume is this micropipette set to dispense? _____ μl</p>
---	--

2. **Tips for Optimal Micropipetting Technique.** Of all the factors contributing to the performance of a micropipette, *the skill of the operator is the most critical!* Here are some tips to help improve your micropipette technique. These tips are demonstrated in an excellent video here: <http://www.artel-usa.com/tip1.aspx>

- Pre-wet the tip.** Aspirate and entirely expel liquids three times before dispensing. This increases the humidity within the tip and reduces sample loss due to evaporation.
- Dispense liquids at ambient temperature.** Allow liquids to ambient temperature before dispensing. Humidity and pressure are temperature-dependent and therefore will affect the volume dispensed from the tip.
- Examine the tip before and after dispensing.** Remove droplets from the side of the tip, and ensure there are no air bubbles in the tip before dispensing. After dispensing touch the tip to the side of the container if there is a droplet attached to the side of the tip.

- d. **Use standard mode.** Depress the plunger to the first stop, immerse into liquid, aspirate by releasing the plunger slowly. Remove tip from the liquid and depress the plunger to the second stop to dispense the entire contents.
- e. **Pause consistently after aspiration.** Pause for one second before removing from the liquid. If aspirating a viscous liquid, pause for at least 3 seconds or until you can see no further liquid moving up the tip.
- f. **Pull the pipette straight out.** This is particularly important for volumes less than 50 μ l. If you hold your tip at an angle, it can alter the volume aspirated.
- g. **Minimize handling of the pipette and the tip.** Touching the tip and tube with a gloved hand will warm them up resulting in inaccurate aspiration. Touching the tip with a bare hand can lead to contamination.
- h. **Immerse the tip to the proper depth.** Immerse a p1000 5-6mm, and smaller micropipettes 2-3mm only. Too little immersion will lead to aspiration of air, too much can result in liquids clinging to the side of the tip, or touching the bottom causing incorrect aspiration volumes.
- i. **Use the correct pipette tip.** Use the manufacturer's recommended tip for the micropipette. Remember to use barrier tips when using biohazards or to avoid cross-contamination of your experiment is necessary (such as working with PCR).
- j. **Use consistent plunger pressure and speed.** Depress and release the plunger smoothly and slowly and consistently.

Part III: Practicing operating a micropipette. Follow the steps below, checking off each step as you go. Each student will perform these tasks independently.

GLP Tip: Never lay a micropipette down with a filled tip or hold it upside down or sideways. The liquid will not leak out if you hold it upright but it may enter the instrument if you hold it upside down, and contamination will result.

1. Practice setting the volume on the micropipette. Look at the top of the micropipette to identify its measuring range. Remember that the highest value listed on the top is the largest volume you can measure on that pipette. On a 100 to 1000 μ L micropipette, the most significant measurable amount is 1000 μ L; on a 20-200 micropipette, it is 200 μ L. Likewise, the smaller value in the range is the lowest measurable volume; on a 2-20 μ L micropipette, the lowest measurable volume is two μ L.
2. Set a 100-1000 μ L micropipette to 500ul, a 20-200 μ L micropipette to 150ul, and a 2-20 μ L to 20ul. What are these values in mL?

500 μ L _____ mL

150 μ L _____ mL

20 μ L _____ mL

3. Obtain 3, graduated 1.5 mL microcentrifuge tubes and place them in a microcentrifuge tube rack. Close the lids, and label the tops: 500ul, 150ul, and 20ul. Each student should have their own rack and perform these steps. Obtain a tube of colored sugar solution – you may share with your lab partner.
4. Micropipette 500 μ L of *colored sugar solution* into a 1.5 mL centrifuge tube as follows:
 1. What micropipette will you use for this? _____
 2. Set micropipette to 500ul.
 3. Place appropriate tip on the end of the micropipette and close tip box lid.
 4. Using one hand, hold the micropipette and press down on the plunger with your thumb to the first stop.
 5. Open sugar solution with your other hand and bring up to eye level.
 6. Submerge the end of the tip just under the surface of the colored sugar solution and slowly release thumb aspirating up into the tip the colored solution.
 7. Notice that you have no bubbles in the tip. If you do, eject tip into waste container and repeat.
 8. Open 1.5mL centrifuge tube, bring to eye level, and dispense liquid into the tube. Close lid. Eject tip into a waste container.
 9. Verify the volume on the side of the tube. If it is incorrect, repeat. Keep tube.
5. Repeat this measuring 150ul. Keep tube.
6. Repeat this measuring 20ul. Keep tube.
7. Have your lab partner verify your volumes.
8. Show instructor all three tubes before you continue!

GLP TIPS!

- a. Always keep your micropipette tip box lids closed, and your microcentrifuge tubes closed at all times. Open only when using at that moment. This avoids external contamination into your tubes (such as splashing, or dust, or your hair or skin cells...).
- b. Always label the top so that it can be read without removing the tube from the rack, and orient the tubes in the same direction so that you won't confuse letters like "H" and "I" and numbers like "6" and "9". Also, only use a permanent marker, such as a Sharpie, that will not erase or bleed if it gets wet. If your tubes are to be stored or mixed in a microcentrifuge, label with your initials.

Part IV Testing Your micropipetting Skill

Each student will perform this task independently.

1. Retrieve a microcentrifuge rack, 10 graduated 1.5mL tubes and colored solution.
2. Label the top of tubes with a Sharpie 1-10.
3. Using the table below, measure the following amounts into the indicated tubes. Keep lids closed at all times!
4. Mix the contents by 'vortexing' briefly.
5. Use the pecofuge to 'pop-spin' the liquid back down into the tube.
6. Check the graduations on the side of the tube to ensure you have aliquoted the liquid correctly.
7. You can check the accuracy of your measurements by setting a micropipette to the total volume that is supposed to be in the tube and slowly withdrawing all of the solution from several tubes. Your pipetting was accurate if you leave no solution behind and have no air bubble in your tip.
8. Have your lab partner check your tubes.
9. Have your instructor check your work before discarding tubes in general trash.
10. Return all equipment and 50mL conical tubes of colored liquid to the cart.

Tube #	Contents	Tube #	Contents
1	5 μ L blue	7	100 μ L clear
2	10 μ L blue		20 μ L blue
3	100 μ L blue	8	500 μ L clear
4	1000 μ L blue		20 μ L blue
5	5 μ L clear 20 μ L blue	9	1000 μ L clear 20 μ L blue
6	20 μ L clear 20 μ L blue	10	500 μ L clear 500 μ L blue

Table 2-1: Testing micropipetting skill

References:

1. Ellyn Daugherty, *Biotechnology: Science for the New Millennium*. 2012. EMCParadigm Publications. ISBN: 978-0-76384-284-0
2. Seidman & Moore, *Basic Laboratory Methods for Biotechnology: Textbook & Laboratory Reference*, 2nd edition. 2009. Prentice Hall. ISBN: 0321570146
3. "10 Tips to improve your pipetting technique". 2009. Artel. www.artel-usa.com

Lab Unit 2-B: Calibrating & Using Basic Lab Equipment

OBJECTIVES

Your performance will be satisfactory when you can

- ✓ Calibrate a pH meter using standard pH buffers and an SOP
- ✓ Use a pH meter to determine the pH of an unknown solution correctly
- ✓ Calibrate an electronic balance using standard weights and an SOP
- ✓ Use an electronic balance to obtain the desired mass of a substance
- ✓ Learn the measuring accuracy of glassware – determine the appropriate glassware to use to measure a specified volume accurately

INTRODUCTION

Different pieces of lab equipment are designed to measure properties such as temperature, pH, mass, and volume to varying degrees of accuracy. If the temperature markings on the side of a thermometer are not set accurately, the instrument's measurements will not be accurate. The accuracy of these markings is due to the **calibration** of the thermometer.

Calibration, Verification & Validation

Micropipettes are calibrated by the manufacturer before they are sent to you, but they do become less accurate the more they are used. Therefore, the performance of a micropipette should be verified periodically. GMP and ISO laboratories have written policies for performance evaluation for micropipettes.

Calibration is a process that compares a known (the “standard” device) against an unknown (the target device in question). During the calibration process, the offset between these two devices is quantified, and the target device is adjusted back to tolerance (if possible). A calibration report usually contains both “as found” and “as left” data. When a micropipette is considered to be out of calibration, it is typically sent to the manufacturer for recalibration.

Verification is simply the process of “verifying” that a device is in tolerance (within an acceptable range). Verification usually results in “as found” data. If the device is not in tolerance, it is sent for recalibration.

Validation is a detailed process of confirming that the instrument is installed correctly, that it is operating effectively, and that it is performing without error. Validation is broken into three different tests: the installation qualification (IQ), the operational qualification (OQ), and the performance qualification (PQ).

Some equipment must be periodically calibrated because the settings are not as immovable as lines on a graduated cylinder or thermometer. The calibration of instruments such as pH meters, electronic balances, and micropipettes can be rendered inaccurate by factors such as movement, humidity, dirt, electrical field changes, and many others.

PROPER USE BASIC LAB EQUIPMENT

Use of a serological pipet and electronic pipet-aid. Serological, or “blowout,” pipets are graduated tubes used to measure anywhere from 0.1 to 50 mL. They are typically made of plastic and single-use (disposable) and have the top end plugged with cotton to prevent contamination and overflow. When the liquid has drained from this pipet, the final drop in the tip is transferred (pushed out) with a puff of air. These are known as TC or to contain pipets.

When filling a pipet, bring the container with the liquid and the pipet to eye level. The tapered end is held beneath the surface of the liquid at all times. The liquid is drawn into the pipet by suction until the level is just above the volume of liquid to be delivered, then bring the level down to the meniscus of the desired volume. When reading the volume, ALWAYS view the pipet dead-on at eye level with the pipet held vertically, perpendicular to the ground. Pipets are designed to be used with a hand pump or bulb, of which there are many varieties. Most biotech labs now use electronic devices called pipet-aids. This helps save time, carpal tunnel in the wrist, and improve accuracy.



Let's Explore

Watch the following video on how to use a serological pipet correctly:

<https://youtu.be/aei-tU1ZIkE>

With a lab partner, go to this website and test your serological pipet and micropipetting skill!

<http://www.lsteam.org/iet/pipetting/>

Briefly describe below how to use a serological pipet to measure 5mL of liquid accurately. Be specific.

Calibrating and Operation of an Electronic Balance. The standards used to calibrate electronic balances are objects of known mass. For balances that measure to $\pm 0.01\text{g}$, the standard is usually a 200-gram weight. These balances are used to measure amounts over 0.05 g. When you place the 200-gram weight on balance in calibration mode, the balance recognizes the weight as 200 grams, and will then use that information to measure other masses. Refer to SOP in SOP booklet on the calibration and operation of our specific balance model (there are several!).

General guidelines to follow when using balances are as follows:

- The amount being weighed determines what type of balance is used in this laboratory (top-loading or analytical). Typically < 1g is best weighed on an analytical balance.
- Use clean spatulas to weigh out material.
- Never put excess chemicals back into their original containers. Always ask your instructor how to discard of excess chemicals.
- Leave a CLEAN balance. Chemicals left on balance will corrode it.
- Turn off the balance when not in use.

Calibration and Operation of a pH Meter: Read the SOP provided in the SOP booklet of the module of pH meter used in your lab. The definition of pH is the hydrogen ion (H⁺) concentration of a solution. By definition, any solution with a pH < 7 is acidic and any solution with a pH > 7 is basic. In this laboratory, pH is measured by using a pH meter. The pH meter measuring system consists of a voltmeter that measures voltage, two electrodes and the sample that is being measured. When the two electrodes are immersed in a sample, they develop an electrical potential (voltage) that is measured by the voltmeter. **Review the SOP for “Operation and Maintenance for a pH meter” before class.**



Let's Explore!

Watch this video on calibrating the Acument pH Meter we use in our labs. Follow along with your SOP booklet and make additional notes in your SOP booklet <https://youtu.be/UfdKhja6u2I>
Draw a flowchart below of how to calibrate a pH meter.

Calibration of Common Laboratory Equipment

Purpose: The purpose of this lab is to calibrate two pieces of equipment; pH meter and a top loading scale. The calibrated pH meter will be used to determine the pH of an unknown solution. The calibrated scale will be used to determine the most appropriate, and an accurate measuring tool to measure 10 mL.

Safety: In this laboratory exercise we will not be using any hazardous materials, operations or dangerous equipment. However, it's part of Good Laboratory Practice to always wear close-toed shoes, lab coat, safety glasses and gloves when working in a laboratory.

Materials

Large waste beaker	pH meter with calibration SOP
Solution of unknown pH	top loading balance with calibrated SOP
1, 10 mL graduated cylinder	200 g standard mass for balance calibration
1, 50-mL beaker	3 pH standard buffers for calibration
Large weigh boat	3, 30mL beakers
Color solutions	Squirt bottle
Kimwipes	10 & 5 mL serological pipet
Pipet-aid (automatic pipet)	Transfer disposable pipet
1, 50mL Erlenmyer flask	

PROCEDURE

A. Calibrate a pH meter. Check off each step as you move along.

1. Collect all materials and set up your workbench with your lab partner.
2. Select a pH meter to calibrate and record its model & unit number. _____
3. Turn on and calibrate the pH meter using the SOP booklet provided.
4. Always rinse the pH electrode with a wash bottle of distilled water, catching the rinse liquid in a labeled waste beaker.
5. Wipe the tip of the electrode with a Kimwipe before inserting into another liquid.
6. After calibration is completed – Did you receive the 'Good Electrode' message? If not, clear attempt and repeat the calibration one more time. If you did not get a 'good electrode' message, ask your instructor for assistance.
7. After calibration is complete, wash probe, tap dry and use to measure the pH of an unknown solution. Swirl the solution slightly to ensure that it makes good contact with the electrode before recording a reading.
8. Record the pH _____
9. Check with instructor if pH is correct. If not, verify the pH probe by measuring the pH of buffer 7. Is it close to 7.0? If not, repeat experiment!

B. Calibrating and Using an Electronic Balance. Check off each step as you move along.

1. Retrieve all materials needed and set up your workbench with your lab partner.
2. Select a balance and record the model number. _____
3. Follow the instructions in the SOP provided for the balance model to calibrate the balance.



4. Verify calibration by re-measuring 200g weight.
5. Place a weigh boat on the balance pan and press the tare button. What does TARE do?

6. Place a 5 mL pipet on an automatic pipet-aid.
7. Open the lid to the colored solution and bring up to eye level.
8. Draw up 4 mL of colored water using a 5 mL pipet using an automatic pipet-aid.
9. Dispense carefully into tared weigh boat.
10. Record the mass of the orange water. Each lab partner will repeat this measurement 3 times, tare the weigh boat in between each measurement:
Attempt 1: _____ g.
Attempt 2: _____ g.
Attempt 3: _____ g.
Average: _____ g.
11. At sea level, 1.00 mL of water weighs 1.00 g. Is your measurement accurate? _____

C. Measuring Accuracy of Glassware. This section may be completed by both lab partners together.

GLP Tip: Always inspect glassware before use. This is especially important if you are handling hazardous materials or heating the glassware. Show your instructor any chipped, cracked or broken glassware. Broken glassware is disposed of in Broken Glassware Waste box using a dustpan and broom.

1. Place weigh boat on scale and tare.
2. Measure 10 mL of water in a 50 mL beaker using the lines on the beaker for your measurement.
3. Pour the water into the weigh boat and record the mass below. Don't forget your units!
4. Repeat three times.
5. Calculate the average mass of 10 mL of water measured with a beaker.
6. Repeat step 1-5 using a 10 mL graduated cylinder
7. Repeat step 1-5 using a 50 mL Graduated Cylinder.
8. Repeat step 1-5 using a 10 mL pipet.
9. Repeat step 1-5 using a 50 mL flask.
10. Clean up your lab station. Pour the contents of the waste beaker down the sink with plenty of water, and return dirty glassware to the cart.
11. Spray and wipe off your lab bench with a paper towel.
12. Return all leftover reagents and equipment where you got them.
13. Ensure the balances, pH meters, and the areas around them are clean and dry and turn off the equipment.

RESULTS

	50 mL Beaker	10 mL Grad Cylinder	50 mL Grad Cylinder	10 mL Pipet	50 mL Flask
Measurement 1					
Measurement 2					
Measurement 3					
AVERAGE:					

Table 2-1: Measurement of 10mL of liquid using various pieces of lab equipment.

Lab Unit 2: Basic Tools in the Biotechnology Laboratory

1. Why should you avoid touching the micropipette tips with your bare fingers?
2. What happens if you push the plunger to the second stop before drawing up the liquid?
3. Why should you always keep the lids on the microcentrifuge tubes closed?
4. Why is it important to verify or calibrate lab equipment before use?

5. The definition of a gram is the mass of 1 mL of pure water at 20°C (about room temperature) and 1 atmosphere of pressure. Referring to the 4mL pipetting exercise:

What should be the average mass of the water your group measured? _____ g

What was the average mass of the water your group measured? _____ g

6. Summarize the average mass of the 10mL of water you measured (include units) below:

50 mL Beaker _____ 10 mL graduated cylinder _____

50 mL graduated cylinder _____ 10 mL pipet _____

50mL Flask _____

Circle the one you predict to measure most accurately. Put a star beside the one that did measure most accurately.

7. Was your predicted most accurate glassware the same as your measured most accurate glassware? Explain any discrepancies.

CONCLUSION: Write a short conclusion statement for this laboratory exercise.

How successful were you at using the micropipette? In what way(s) can you become more skillful? Which glassware was more accurate in your experiment? Is this an expected result?



LAB UNIT 3: PREPARING SOLUTIONS

Jack O'Grady, M.S., Linnea Fletcher, Ph.D., Angela Wheeler, M.S., Patricia Phelps, Ph.D.

OBJECTIVES

Your performance will be satisfactory when you can:

- ✓ Correctly prepare a solution of a given molarity
- ✓ Calculate and perform parallel and serial dilutions and distinguish between the two
- ✓ Use a microcentrifuge to pellet a precipitate
- ✓ Prepare a graph using MS Excel, and determine linear regression, R^2 value

INTRODUCTION

A common task for any Biotechnician is solution preparation. This is an essential skill that must be mastered to be a successful technician in the workplace. You will be expected to prepare many different kinds of solutions correctly 100% of the time. An incorrectly prepared solution can cost a biotech company a lot of money and time. In a pharmaceutical company, an improperly prepared solution can harm someone.

LET'S GET STARTED!

Reading:

Chapter 3.2-3.6: The Basic Skills of the Biotechnology Workplace!

Additional Resources:

1. The optional Seidman textbook has several chapters on solution making, with many step-by-step instructions and practice questions.
2. The ChemCollective has many online tutorials: <http://www.chemcollective.org/tutorials.php>

Let's start with the basics - What is a solution? It is defined as a solute (smaller amount) dissolved in a solvent (larger amount). The concentration of a solution frequently must be known to a high degree of accuracy. An incorrectly prepared solution can destroy months of hard work or cost companies thousands of dollars. Therefore, companies usually have an SOP (Standard Operating Procedure) for the preparation of each solution to minimize mistakes. All calculations are recorded in the lab notebook, even if a calculator is used. Critical calculations are double-checked by another person (and sometimes triple-checked). The exact mass and volume of reagents used are recorded in the notebook. This information, along with the date and the preparer's name or initials, is recorded on a solution preparation form and on a label on the bottle itself; these forms are provided in the Appendix.

There are several critical aspects to making solutions that should be followed at all times.

- Check and recheck each calculation. It is best if two people make a calculation independently and then cross check their answers.
- Read each reagent bottle twice, once before using and once afterward, to ensure that the right reagent is used.
- Complete a Solution Prep form for every solution you prepare. Solution prep forms should include detailed and relevant information. See Appendix instructions and example.
- Label each bottle before filling. Use tape and a permanent marker and labeling SOP.
- Record any changes observed, no matter how trivial. Your notes can be used to trace back a problem to its source quickly.

Cleaning Glassware. Properly cleaning glassware is one of the most critical aspects of the job of a biotechnician. Improper cleaning can have disastrous and costly consequences for a company! Cleaning glassware for the lab is a lot more involved than just sticking it in a dishwasher. You will need to consider the type of glassware itself, what the dirty glassware was used for, and most importantly, what is the purpose of the glassware. Note, the water source is a vital part of cleaning glassware!

The five common steps to washing glassware (or plastic ware) are as follows.

1. **Pre-rinse:** Soak or pre-rinse all glassware after use, to prevent the contaminants from drying onto the glassware.
2. **Contaminant Removal:** Wash using approved detergents and solvents along with scrubbing will help with contaminant removal. Typically, a lab glass detergent such as Alcon is used. As with all washing, hot water is preferable to cold water. Often lab brushes are used to help wash debris off glassware.
3. **Rinse:** The rinse step is essential in removing the detergent and cleaning solvents. Many SOP's specify that glassware is rinsed 3 to 5 times in tap water.
4. **Final Rinse:** Always use purified water for the last rinse, and is usually performed 1-3 times.
5. **Drying:** This is either done in the air upside down on a rack or by heat. ***Never hand dry!*** Clean, dry glassware is often stored covered in a closed cupboard to avoid contamination.

A glassware washing SOP is included in your SOP booklet. If you are required to wash glassware in lab class, you MUST follow this SOP!

Why is glassware washing a critical part of a biotech company operation?

**TEST YOUR
KNOWLEDGE!**

PART I: PREPARING A MOLAR SOLUTION

Molarity. Molarity is the most common unit of concentration in the biotechnology lab. The molarity of a solution is defined as the number of moles of solute per liter of solution. The symbol for molarity is M, but it can also be written as moles/Liter, or mol/L. A mole of any element always contains 6.02×10^{23} (*Avogadro's number*) atoms. Because some atoms are heavier than others, a mole of one element weighs a different amount than a mole of another element. *The weight of a mole of a given element is equal to its atomic weight in grams.* Consult a periodic table of elements to find the atomic weight of an element. For example, one mole of the element carbon weighs 12.0 g.

Example: Using a periodic table, calculate the molar mass of chromium oxide (CrO_2). The atomic weight of chromium is 52.00, and that of oxygen is 16.00. You must count the oxygen twice because there are two per formula unit of chromium oxide.

ANSWER: $52.00 + 2(16.00) = 84.00 \text{ g/mol}$

Practice:

Using a periodic table, calculate the molar mass of potassium sulfate (K_2SO_4).

Molecular Mass K_2SO_4 : _____ Units? _____

We can't directly measure moles, but we can measure mass. To calculate the mass of a chemical needed to prepare a given volume of a solution of desired molarity, you must convert the number of moles to mass, using the molar mass as a conversion factor.

$$\begin{array}{rclcl} \text{Mass} & = & \text{molarity} & \times & \text{volume} & \times & \text{molar mass} \\ \underline{\text{?}} \text{ g} & = & \underline{\quad} \text{ moles/liter} & \times & \underline{\quad} \text{ L} & \times & \underline{\quad} \text{ g/mole} \end{array}$$

Don't forget to convert mL to L, if necessary.

Example: How will you prepare 100 mL of 1 M NaOH (MW 40.0 g/mol)?

Calculations: $\text{g} = 40 \text{ g/mol} \times 1.0 \text{ mol/L} \times 0.1 \text{ L} = 4 \text{ g}$

Protocol:

1. Place weigh boat on an electronic balance and press tare.
2. Weight 4.0 g of NaOH and pour pellets into 100 mL beaker
3. Add 70 mL of deionized water and a stir bar.
4. Place on a stir plate and stir until dissolved.
5. Pour into 100 mL graduated cylinder and bring to 100mL with diH₂O. Cover and mix.
6. Transfer to a labeled bottle.
7. Fill in a solution preparation form.

Practice: What mass of NaCl is needed to prepare 100mL of a 0.25M NaCl solution (58.44 g/mol)?

PART I: PREPARATION OF STOCK SOLUTIONS

PURPOSE: The purpose of this lab is to prepare 2.0M CaCl₂ and 2.0 M MgSO₄ stock solutions.

SAFETY: Reagents used in this experiment may cause eye irritation if splashed in the eye. Use PPE, including gloves, lab coat, and safety glasses. For disposal, dump down the drain with lots of water separately. When combined these reagents will form a precipitate and may clog the drain.

MATERIALS

2, 30 or 50-mL beakers
2, 25-mL graduated cylinders
labeling tape
permanent marker
50 mL conical (2)

Calcium chloride – CaCl₂ • 2H₂O
Magnesium sulfate – MgSO₄ • 7H₂O
Top loading balance with SOPs
2, spatulas
2, weigh boats
stir plate & stir bar

PROCEDURE

In this lab procedure, you and your lab partner will split the preparation of both solutions; one will prepare a 2.0M solution of calcium chloride, and the other will prepare a 2.0M solution of magnesium sulfate. Before you begin, perform calculations, and check calculation with your lab partner!

Before class calculate the molecular weight of the two reagents you will be using today:

a. Calcium chloride – CaCl₂ • 2H₂O _____ g/ mol

b. Magnesium sulfate – MgSO₄ • 7H₂O _____ g/ mol

When you get to class, and you are collecting your reagents – VERIFY the molecular weight you calculated by verifying the molecular weight of the reagent written on the side of the bottle. Are they the same? Different?

TEST YOUR
KNOWLEDGE!

Protocol for the preparation of 25mL of a 2M CaCl ₂ solution	Experimental Notes
<ol style="list-style-type: none"> 1. Gather all materials needed & label beakers and tubes using tape. 2. Place weigh boat on the scale and press tare. 3. Weigh _____ g of CaCl₂ 4. Pour into 50mL beaker 5. Add 15mL of diH₂O to the beaker along with a stir bar. 6. Place on a stir plate and dissolve. 7. Transfer solution to a 25mL graduated cylinder and BTV (Bring to volume). 8. Cover with Parafilm and mix by inversion 9. Transfer to a labeled 50mL conical tube. 10. Fill in solution preparation form. 	<p>Scale Model & #: _____</p> <p>Actual grams CaCl₂ weighed: _____ g</p> <p><u>Show Calculations:</u></p>

Your Turn! Write a protocol below for the preparation of 25mL of 2.0 M MgSO₄:

Protocol for the preparation of 25mL of a 2M MgSO ₄ solution	Experimental Notes
<ol style="list-style-type: none"> 1. Gather all materials needed & label beakers and tubes using tape. 2. Place weigh boat on the scale and press tare. 3. Weigh _____ g of MgSO₄ 4. 	<p>Scale Model & #: _____</p> <p>Actual grams MgSO₄ weighed: _____ g</p> <p><u>Show Calculations:</u></p>



PART II: PREPARING PARALLEL DILUTIONS

Parallel dilutions consist of adding additional solvent (usually water) to a solution to reduce its concentration. There are many ways of expressing dilution factors.

1. Combining one part food coloring with nine parts water dilutes the food coloring to 1/10. This means that there is one part food coloring in 10-parts total volume. The denominator an expression such as 1/10 is the total volume of the solution.
2. The food coloring dilution above can also be referred to as 1:9 food coloring to the water. The colon (:) means "to." A 1:10 food coloring to water dilution would be 1/11, not 1/10 because the total number of parts is 11.
3. Frequently, stock solutions in biotechnology labs are concentrated and must be diluted before using. A buffer that is ten-fold more concentrated than the usable concentration is referred to as a 10X solution. One must dilute a 10X solution by a factor of 10 (by adding 1 part of the 10X stock to 9 parts of solvent) before using.
4. In dilutions, parts can be of any unit. If you combine 1 mL food coloring with 1mL water, you are using the same dilution factor (1:1 or ½) as the person who combines 1 ounce of food coloring with 1 ounce of water. If you combine one ounce of food coloring with one liter of water, the dilution factor is not 1:1, because the units are not the same.
5. To dilute a more concentrated stock solution to a less concentrated solution the following formula is used:

$$C_1V_1 = C_2V_2$$

Where,

C_1 = original concentration (of stock solution)

C_2 = final concentration (of diluted solution)

V_1 = original volume (to be taken from stock solution)

V_2 = final volume (of diluted solution)

Example: Calculate how many mL of a 1.0 M stock solution of NaCl are needed to prepare 100 mL of a 0.050 M solution (also referred to as 50 mM).

$$C_1V_1 = C_2V_2$$

$$(1.0 \text{ M})(? \text{ mL}) = (0.050 \text{ M})(100 \text{ mL})$$

$$\text{mL} = \frac{(5.0 \text{ M})(\text{mL})}{1.0 \text{ M}}$$

$$\text{mL} = 5.0$$

PART II: PREPARATION OF PARALLEL DILUTIONS

PURPOSE: The purpose of this lab is to prepare parallel dilutions using CaCl_2 and MgSO_4 stock solutions. These solutions will be used to determine if there is a linear relationship to the ratio of CaCl_2 and MgSO_4 with a mass of precipitate produced.

SAFETY: Reagents used in this experiment may cause eye irritation if splashed in the eye. Use PPE, including gloves, lab coat, and safety glasses. For disposal, dump down the drain with lots of water separately. When combined these reagents will form a precipitate and may clog the drain.

MATERIALS

Each group

One, 10 mL graduated cylinder
1, 5 mL serological pipet
1, 10 mL serological pipet
100 – 1000 μl micropipette
microcentrifuge tubes (9)
microcentrifuge tube rack

Class Shares

Top loading balance
Microcentrifuge
2, 15 mL conical tubes
transfer pipets
 CaCl_2 & MgSO_4 stocks prepared in Part A
permanent marker

PROCEDURE: Parallel Dilutions

*Calculate the volume of the 2M CaCl_2 stock solution made in Part A needed to make 10.0 mL of a 1.50 M solution. Prepare the calculations in your exercise book BEFORE you come to the lab. When you get to the lab, verify your calculations with your lab partner before you begin.

Protocol : Preparing 1.5M & 0.5M CaCl_2	Experimental Notes
<p>Part I: Preparing solutions: 1.5M CaCl_2</p> <ol style="list-style-type: none">1. Label a 10mL graduated cylinder 1.5M CaCl_22. Using a _____ mL serological pipet, add _____ mL of 2.0M CaCl_2.3. BTV to 10mL with dH_2O.4. Cover and mix.5. Transfer to a labeled disposable 15mL conical. <p>Part I: Preparing solutions: 0.5M CaCl_2</p> <ol style="list-style-type: none">1. Label a 10mL graduated cylinder 0.5M CaCl_22. Using a _____ mL serological pipet, add _____ mL of 2.0M CaCl_2.3. BTV to 10mL with dH_2O.4. Cover and mix.5. Transfer to a labeled disposable 15mL conical.	<p>Calculations:</p> <p>A: 1.5 M CaCl_2 (using 2.0 M stock)</p> <p>B: 0.5 M CaCl_2 (using 2.0 M stock)</p>



Part II: Mixing solutions & form precipitate

1. Label 9 microcentrifuge tubes 1-9.
2. Weigh each tube and record the empty tube weight in a data table below.
3. Pipette 500 μl of 2 M magnesium sulfate stock solution into each tube.
4. Add the following to the indicated tubes:
 Tube 1, 2, 3 —500 μl of 2 M CaCl_2
 Tube 4, 5, 6 —500 μl of 1.5 M CaCl_2
 Tube 7, 8, 9 —500 μl of 0.5 M CaCl_2
5. Place tubes in a microcentrifuge and ensure they are balanced!
6. Spin the tubes in a microcentrifuge for 5 minutes at 10,000rpm.
7. Discard the water by using a micropipette or transfer pipet. Do not disturb the pellet.
8. Weigh the tubes with pellets and record in the table below.
9. Perform calculations as outlined in the table.
10. Discard the stock solutions down the sink with plenty of water unless your instructor asks you to save them.

RESULTS:

Tube Number	CaCl ₂ Concentration (M)	Pre-weight (g)	Post weight (g)	Weight of PPT (g)
1	2			
2	2			
3	2			
Average				
4	1.5			
5	1.5			
6	1.5			
Average				
7	0.5			
8	0.5			
9	0.5			
Average				

Table 3-1: Weight of precipitate formed when various concentrations of CaCl_2 are added to 2.0M MgSO_4 .

PART III: PREPARING SERIAL DILUTIONS

Serial dilutions are dilutions made from other dilutions. They are made for one of the following reasons:

1. A series of dilutions with the same dilution factor is desired.
2. The final concentration desired is so small that the original volume (C_1) cannot be accurately measured.

For serial dilutions:

$$\text{Dilution factor} = (V_1 + V_2) / V_1$$

Where V_1 is the volume of the solution being diluted
 V_2 is the volume of solvent used to dilute the solution

(Note: V_2 is also the ending volume of the diluted solution)

Example, how will you prepare 10, 5-fold dilutions of a reagent with an end volume of 40mL in each of the ten tubes?

$$\text{Dilution Factor} = (V_1 + V_2) / V_1$$

$$DF=5$$

$$V_2 = 40 \text{ mL}$$

Find V_1

$$5 = (V_1 + 40\text{mL}) / V_1$$

$$5V_1 = (V_1 + 40\text{mL})$$

$$4V_1 = 40\text{mL}$$

$$V_1 = 10 \text{ mL}$$

PROTOCOL:

1. Using a 50mL graduated cylinder dispense 40mL of diH₂O into each of 10 labeled beakers
2. Using a 10mL pipet, dispense 10mL of stock solution into beaker 1 and mix.
3. Using a clean pipet, transfer 10mL of the diluted mixture in beaker 1 into beaker 2 and mix.
4. Continue until beaker 10 – discard 10mL out of beaker 10 to leave 40mL remaining.

IMPORTANT TO NOTE: No, really, this is important, pay attention to this!

- a. The V_1 & V_2 in this formula is NOT the same as the V_1 & V_2 in $C_1V_1 = C_2V_2$.
- b. In this formula, V_2 is both the volume of solvent used in each dilution and the final volume of that dilution. Why? Once you have made a dilution by adding V_1 mL of solution to V_2 mL of water, you remove V_1 mL of that dilution to make the next one. Thus, you always end up with V_2 mL in each dilution but the very last one.

PART III: PREPARATION OF SERIAL DILUTIONS

PURPOSE: The purpose of this lab is to prepare serial dilutions of HCl and determine the relationship between HCl concentration and pH. These solutions will be used to determine if there is a linear relationship between pH and HCl concentration ranging from 1M to 0.001M HCl.

SAFETY: Concentrated solutions of hydrochloric acid are extremely corrosive. The risk of damage through contact lessens as the acid is diluted; very dilute solutions are only mildly corrosive. The concentrated acid releases dangerous quantities of hydrogen chloride vapor; inhaling this can be extremely harmful. Contact with the eyes or skin can cause severe permanent damage. Always wear safety glasses and nitrile gloves when handling concentrated HCl. Do not allow the acid or a solution of it to come into contact with your skin. The concentrated acid must always be used in an area equipped with proper ventilation. Dispose of waste in WASTE BEAKER IN CHEMICAL FUMEHOOD ONLY.

MATERIALS

Each group

Four, 100-mL beakers
One, 50-mL graduated cylinder
3, 5- or 10-mL serological pipets
labeling tape
permanent marker

Class Shares

1.00 M HCl
pH meters with calibration SOP
pH standard buffers and three 30 or 50-mL beakers
250-mL waste beaker washing pH probe
HCL waste beaker (Chemical hood)

PROCEDURE: Serial Dilutions

*Calculate the volumes of hydrochloric acid (HCl) and water (solvent) needed to prepare a serial dilution of a 50 mL final volume, of each of 0.1 M, a 0.01 M, and 0.001 M HCL solution. Prepare the calculations below **BEFORE** you come to class. Verify your calculations with your lab partner.

CALCULATIONS:

What formula will you use?

Protocol	Experimental Notes
<ol style="list-style-type: none"> 1. Collect all materials before you begin. Both you and your lab partner can perform this experiment together. 2. Using tape, label 3, 100 mL beakers with the concentrations of HCl. 3. Using a 50mL graduated cylinder add 50 mL of diH₂O to each of 100mL beaker. 4. Using a _____ mL pipet, dispense _____ mL of 1M HCl to the beaker labeled 0.1M HCl. 5. Add a stir bar, place on stir plate, and mix. 6. Using a clean _____ mL pipet, transfer _____ mL of 0.1M HCl into the 0.01M HCl beaker and mix. 7. Using a clean _____ mL pipet, transfer _____ mL of 0.01M HCl into the 0.001M HCl beaker and mix. 8. Calibrate a pH meter using calibration SOP 9. Measure pH of each solution. 10. Create a data table below and record data. Remember to label data table with a title. 11. <u>Pour your acid solutions into a labeled waste container.</u> Clean up your work area and the area around the pH meters. 	<p>pH meter: _____</p> <p>Temperature: _____</p>

RESULTS:

Using a ruler, create a table below for your data. Include a descriptive title at the bottom.

Table 3-2: (Provide a title):

PART IV: GRAPHING DATA

As a biotechnician, you will often make scatter diagrams and line graphs to illustrate the data that you collect. Scatter diagrams are commonly used to show the relationship between two variables. For example, in an absorbance spectrum, the variables would be the wavelength of light and the amount of light absorbed. Although this data is recorded in a table, a scatter diagram can illustrate more visually the relationship between the two datasets. Most companies use MS Excel to plot graphs; therefore during this semester, you will prepare your graphs using this software. Your instructor will demonstrate in class.

How do you know which variable is to be on the x-axis, and which is to be on the y-axis? The x-axis should be the **independent variable** or the parameter that you selected to vary. The y-axis should be the **dependent variable** or the data that you obtained from your measurement.

A “best fit” line for a scatter diagram. If you look at your scatter plot and the middle points on the graph are close to forming a straight line; it is reasonable to conclude that the relationship between the independent and dependent variables is linear. The straight line defines this linear relationship. The most valid best-fit straight lines that illustrate a linear relationship are determined using a type of statistical analysis, called linear regression analysis.

Best-Fit Linear Regression. Linear regression is generated using a computer program or scientific calculator. A correlation coefficient (R^2 value) is a mathematical relationship of linearity of the data. The closer to the value of 1, means the data is linear. The closer to the value of 0.85, indicates the data may be random. Depending on the purpose of the data analysis, you may require an R^2 value to be above 0.95.

Always follow these guidelines when preparing a graph showing experimental data:

1. A graph should be given a descriptive title to explain the experimental data.
2. Adjust both axes so that the completed graph will nearly fill the page. It is not necessary to have a 0,0 axis unless this is a data point.
3. Both axes should be labeled and marked with appropriate units of measurement.
4. The x-axis should show the independent variable; variable that the experimenter chooses and can change. The y-axis should show the dependent variable; the one that the experimenter observes.
5. Draw a line on a scatter plot to illustrate a potential relationship. If there is a linear relationship between the independent and dependent variable ($R^2 > 0.95$), draw a best-fit straight line through the points that are consistent with the linear relationship. If there is no linear relationship, leave the scatter plot as is, or draw curved lines between the data.

Graphing a semi-log plot. Your variables may not have a linear relationship, in which case a straight line cannot represent your data. In the biological sciences, the relationship may be exponential rather than linear. This means one value doubles for each single-unit increase in the other value. For example, each time a cell divides, the number of cells is doubled. This means that if you repeatedly count the number of cells in culture over a given interval of time, the cell count will not rise linearly with time, but rather exponentially with time. On a semi-log graph,

the X-axis is typically linear (each increment is spaced equally and represents an equal unit of measurement), but the Y-axis is exponential (each increment is NOT spaced proportionately and does NOT equal the same unit of measure).

Graphing with Microsoft Excel. In the workplace, you will never graph by hand. Using a computer-based graphing program is an essential skill of any biotechnician. During this semester you will have many opportunities to hone this skill by practicing graphing different types of data using various types of graphs. The most common graphing program used in this field is MS Excel. You may find tutorials for the use of Microsoft Excel spreadsheets and graph at various places online. A thorough tutorial on graphing data will be provided in class. Bring your data to class! You may also attend Biotechnology Program Open Labs, and the Learning Lab to get individual tutoring assistance with graphing. YouTube is an excellent source of information for a specific version of MS Excel. Check it out! For Macs: https://youtu.be/ANz1TY_qg9w and for PCs: <https://youtu.be/MM7dCOuhBs8>

Lab Unit 3 – Preparing Solutions Assignment

1. Using MS Excel, and the data generated from your experiment, create a graph to compare the molar concentration of CaCl_2 (x-axis) to the amount of precipitate produced as shown by the weight of your pellet (y-axis). Include linear regression and R^2 value on the graph. Print out and attach to your report.
2. What is the equation of the line? _____ R^2 value: _____
3. What does the R^2 value tell you about your data?
4. Using MS Excel and the data generated in your experiment, create a graph to compare the molar concentration of HCl and pH. Graph the $-\log_{10}$ HCl concentration on the X-axis and pH on the Y-axis. Include linear regression and R^2 value on the graph. Print out and attach to your report.
5. What is the equation of the line? _____ R^2 value: _____
6. How accurate were your dilutions? How do you know?

CONCLUSION: Write a conclusion statement for your experiment. Did you see a linear relationship in your MgCl_2 pellets? Was there a linear relationship between molarity of HCl and pH?

LAB UNIT 4: DNA BARCODING – ISOLATING GENOMIC DNA

Jack O’Grady, M.S. & Cold Spring Harbor DNA Learning Center

This lab has been modeled from the Cold Spring Harbor DNA Learning Center DNA Barcoding 101 Manual: “Using DNA Barcodes to Identify and Classify Living Things.” 2012
(<http://www.dnabarcoding101.org/>)

OBJECTIVES

This laboratory demonstrates several important concepts of modern biology. During this activity, you will:

- ✓ Collect and analyze sequence data from plants or animals – or products from them.
- ✓ Use DNA sequence to identify & explore relationships between species.
- ✓ Collect plants, animals, or products in your local environment or neighborhood.
- ✓ Extract and purify genomic DNA from tissue or processed material.
- ✓ Analyze genomic DNA preparation using agarose gel electrophoresis

NOTE: This is a research-based lab which will span three labs throughout the semester.

LET’S GET STARTED!

Reading:
Chapter 4: Introduction to studying DNA

INTRODUCTION

Taxonomy, the science of classifying living things according to shared features, has always been a part of human society. Carl Linnaeus formalized biological classification with his system of binomial nomenclature that assigns each organism a genus and species name.

Identifying organisms has grown in importance as we monitor the biological effects of global climate change and attempt to preserve species diversity in the face of accelerating habitat destruction. We know very little about the diversity of plants and animals – let alone microbes – living in many unique ecosystems on earth. Less than two million of the estimated 5-50 million plant and animal species have been identified. Scientists agree that the yearly rate of extinction has increased from about one species per million to 100-1,000 per million. This means that thousands of plants and animals are lost each year. Most of these have not yet been identified. Classical taxonomy falls short in this race to catalog biological diversity before it disappears. Specimens must be carefully collected and handled to preserve their distinguishing features. Differentiating subtle anatomical differences between closely related species requires the subjective judgment of a highly trained specialist – and few are being produced in colleges today.

Now, DNA barcodes allow non-experts to objectively identify species – even from small, damaged, or industrially processed material. Just as the unique pattern of bars in a universal



product code (UPC) identifies each consumer product, a “DNA barcode” is a unique pattern of DNA sequence that identifies each living thing. Short DNA barcodes, about 700 nucleotides in length, can be quickly processed from thousands of specimens and unambiguously analyzed by computer programs.

The International Barcode of Life (iBOL) organizes collaborators from more than 150 countries to participate in a variety of “campaigns” to census diversity among plant and animal groups. Including ants, bees, butterflies, fish, birds, mammals, fungi, and flowering plants and, additionally, within ecosystems – including the seas, poles, rain forests, kelp forests, and coral reefs. The 10-year Census of Marine Life, completed in 2010, provided the first comprehensive list of more than 190,000 marine species and identified 6,000 potentially new species.

There is a surprising level of biological diversity, literally in front of our eyes. For example, DNA barcodes showed that a well-known skipper butterfly (*Astraptes fulgerator*), identified in 1775, is ten distinct species. DNA barcodes have revolutionized the classification of orchids, a sophisticated and widespread plant family with an estimated 20,000 members. The urban environment is also unexpectedly diverse; DNA barcodes were used to catalog 54 species of bees and 24 species of butterflies in community gardens in New York City.

DNA barcodes are also used to detect food fraud and products taken from conserved species. Working with researchers from Rockefeller University and the American Museum of Natural History, students from Trinity High School found that 25% of 60 seafood items purchased in grocery stores and restaurants in New York City were mislabeled as more expensive species. One mislabeled fish was the endangered species, Acadian redfish. Another group identified three protected whale species as the source of sushi sold in California and Korea. However, using DNA barcodes to identify potential biological contraband among products seized by customs is still in its infancy.

Barcoding relies on short, highly variable regions of the mitochondrial and chloroplast genomes. With thousands of copies per cell, mitochondrial and chloroplast sequences are readily amplified by polymerase chain reaction (PCR), even from very small or degraded specimens. A region of the chloroplast gene *rbcl* – RuBisCo large subunit – is used for barcoding plants. The most abundant protein on earth, RuBisCo (Ribulose-1, 5-bisphosphate carboxylase oxygenase) catalyzes the first step of carbon fixation. A region of the mitochondrial gene *COI* (cytochrome c oxidase subunit I) is used for barcoding animals. Cytochrome c oxidase is involved in the electron transport phase of respiration. Thus, the genes used for barcoding are involved in the key reactions of life: storing energy in carbohydrates and releasing it to form ATP.

This laboratory exercise utilizes DNA barcoding to identify plants or animals – or products made from them. First, a sample of tissue is collected, preserving the specimen whenever possible and noting its geographical location and local environment. A small leaf disc, a whole insect, or a sample of muscle tissue is suitable sources. DNA is extracted from the tissue sample, and the barcode portion of the *rbcl* or *COI* gene is amplified by PCR. The amplified sequence (amplicon) is submitted for sequencing in one or both directions.




The sequencing results are then used to search a DNA database. A close match quickly identifies a species that is already represented in the database. However, some barcodes will be entirely new, and identification may rely on placing the unknown species in a phylogenetic tree with near relatives. Novel DNA barcodes can be submitted to the database at the Barcode of Life Data System (BOLD) (www.boldsystems.org) at the University of Guelph.

OVERVIEW OF EXPERIMENTAL METHODS


I. COLLECT, DOCUMENT, AND IDENTIFY SPECIMENS




II. ISOLATE DNA FROM PLANT OR ANIMAL TISSUE



III. AMPLIFY DNA BY PCR



IV. ANALYZE PCR PRODUCTS BY GEL ELECTROPHORESIS



SEQUENCE PCR PRODUCT AND ANALYZE RESULTS





Let's Explore

1. Read the introduction to this lab and watch “The History of DNA Barcoding” video BEFORE class: https://youtu.be/Z_CIP6-E4VY
2. Write down some ideas on a class project you’d like to propose.

3. Answer the following questions:
 - a. What is DNA barcoding?

 - b. Why is DNA barcoding important?

 - c. What are a few examples of practical applications of DNA barcoding?

Bring your ideas to class! What are some ideas your classmates shared?

Write a short paragraph to describe the DNA Barcoding project you have chosen.

Lab Unit 4-A: Collect, Document, and identify Specimens

INTRODUCTION:

The DNA isolation and amplification methods used in this laboratory work for a variety of plants and animals – and many products derived from them. Your collection of specimens may support a census of life in a specific area or habitat, an evaluation of products purchased in restaurants or supermarkets, or may contribute to a more significant “campaign” to assess biodiversity across large areas. It may make sense for you to use sampling techniques from ecology. For example, a quadrant samples the plant and animal life in one square meter (or ¼ square meter) of habitat, while a transect collects samples along a fixed path through a habitat.

Use common sense when collecting specimens. Respect private property; obtain permission to collect in non-public places. Respect the environment; protect sensitive habitats, and collect only enough of a sample for barcoding. Do not collect specimens that may be threatened or endangered. Be wary of poisonous or venomous plants and animals. Consult your instructor if you are in doubt about the safety or conservation status of a potential specimen. You will also need a small sample for classical taxonomic analysis and to act as a reference sample if you plan to submit your data to GenBank.

Do not take more sample than you need. Only a small amount of tissue is needed for DNA extraction – a piece of plant leaf about 1/4 inches in diameter or a piece of animal tissue the size of a pencil eraser. Minimize damage to living plants by collecting a single leaf or bud, or several needles. When possible, use young, fresh leaves or buds. Flexible, non-waxy leaves work best. Tougher materials, such as pine needles or holly leaves, can work if the sample is kept small and is well ground. Dormant leaf buds can often be obtained from bushes and trees that have dropped their leaves. Fresh frozen leaves work well. Dried leaves and herbarium samples are variable. Avoid twigs or bark. If woody material must be used, select flexible twigs with soft pith inside. As a last resort, scrape a small sample of the softer, growing cambium just beneath the bark. Roots and tubers are a poor choice because high concentrations of storage starches and other sugars can interfere with DNA extraction.

REAGENTS, SUPPLIES, & EQUIPMENT

To share

- ◆ Collection tubes, jars, or bags
- ◆ Tweezers, scalpel, and scissors
- ◆ (Smart) phone with a camera or digital camera with GPS
- ◆ Field guide or taxonomic key

PROTOCOL:


1. Collect specimens, according to a strategy or campaign outlined by your instructor. Fresh, healthy, actively growing plant leaves work well with this protocol. “Field Techniques Used by Missouri Botanical Garden” has many useful methods for collecting and preparing plant specimens: <http://www.mobot.org/mobot/molib/fieldtechbook/pdf/handbook.pdf>



2. Use your phone or digital camera to photograph your specimen in its natural environment, or where it was obtained or purchased.
 - a. Take wide, medium, and close-up views.
 - b. Include a person for scale in wide and medium shots. Include a ruler or coin for scale in close-ups.
3. A global positioning system (GPS)-enabled phone or camera stores latitude, longitude, and altitude coordinates along with other metadata for each photo. Visualize or extract this geotag information if appropriate. The class may share their collection location by dropping a pin on a project Google map.
4. Use a field guide or taxonomic key to identify your specimen as precisely as possible: kingdom > phylum > class > order > family > genus > species. Taxonomic keys for local plants or animals are often available online, at libraries, or from universities, natural history museums, and botanical gardens.
5. Check to see if your specimen is represented in the Barcode of Life Database, BOLD (www.boldsystems.org):
 - a. Search by entering genus and species names in the search bar at top right. If the species is represented in the database, the “Taxonomy Browser” will list the number and sources of specimen records.
 - b. Click on “Download Public Sequences” for a FASTA file of available barcode sequences.
 - c. Click on “Taxonomy Browser” at the top left to explore barcode records by the group.
6. Use clean (wash them!) tweezers, scalpel, or scissors to collect a small sample of tissue. You should collect young leaves and flowers from plants.
7. See your instructor about when to bring in your sample. A fresh sample is best. However, you may freeze your sample at -20°C until you are ready to begin DNA extraction.

Lab Unit 4-B: Isolate DNA from Plant or Animal tissue

This universal DNA extraction method uses Promega's Wizard Genomic DNA Purification Kit (Cat # A1120). The background and introduction information come from Bio-Rad's Cloning and Sequencing Manual.



DNA Isolation Overview

The fundamental steps of DNA purification are sample lysis and then purification of the DNA from contaminants. There is a myriad of protocols available for isolating DNA from organisms in the molecular lab. The more “classical” methods have remained virtually unchanged for decades, and the more modern methods involve kits that are commercially available. The best method for any particular application depends on these fundamental considerations:

- ✓ Where the DNA is isolated from will determine the cell lysis techniques used.
- ✓ The purity requirements for the intended use of the DNA being isolated will determine how many purification steps will be involved.
- ✓ The type of DNA being isolated: genomic DNA has different physical properties from those of plasmid DNA.

The isolation of nucleic acids is a common practice in the molecular lab. Necessary steps to nucleic acid isolation may include:

1. The disruption of the cell membrane, and cell wall when necessary, by mechanical, chemical and enzymatic treatment.
2. Enzyme degradation is used for selective isolation of DNA (by RNase treatment) or RNA (by DNase treatment).
3. The separation of nucleic acids from other cytoplasmic components by combinations of these steps:
 - a. Phenol extraction of proteins, followed by selective precipitation of nucleic acids under high salt and cold alcohol treatment.
 - b. Selective precipitation of nucleic acids under high salt or cold alcohol treatment.
 - c. Selective adsorption onto a chromatographic matrix in a centrifuge (a “spin column”) followed by desorption by a unique buffer system.

Cell Lysis

The successful isolation of DNA requires methods that prevent nuclease degradation of the DNA. Some buffer constituents used to promote lysis and denaturation of nucleases include detergents, enzymes to inactivate DNases, denaturants such as guanidinium, and additional organism-specific components such as enzymes for cell wall, or RNase to remove contaminating RNA molecules.

In general, animal tissues are efficiently lysed, because they have no cell wall, and a gently detergent treatment usually is sufficient to break open cells. Yeast and microbial cells, on the other hand, have rigid cell walls that must be weakened enzymatically before the cell will release its DNA. In the case of bacteria, lysozyme enzyme is added, while in the case of yeast a more complex mixture of enzymes must be used to degrade cell wall polymers. Plant cell walls are abraded mechanically by grinding frozen plant tissue, often with glass beads or sand and a mortar and pestle.

Remove Cellular Debris

This step is typically performed with centrifugation. In this protocol, protein precipitation solution is added to aid in reducing the contaminating proteins and nucleases in the supernatant.

Isolate DNA from other cell components

The method you select for your application depends on the size and source of the DNA to be isolated. Genomic DNA can frequently be rendered insoluble and quickly precipitated by addition of alcohol to the mixture.

PCR from Genomic DNA

For PCR to be successful, the DNA extracted needs to be relatively intact. The best sources for DNA extraction includes young green leaves, but fruit, roots, or germinating seeds should also suffice. *It is better to use tissue that is still growing, as the nucleus-to-cytoplasm ratio will be more favorable, cells walls will be thinner, and the amount of potentially harmful secondary products will be less.* There are two features of plants that make DNA extraction different from animals. First, plants have a tough cell wall made of cellulose that has to be penetrated. Second, a significant part of every plant cell is a vacuole that contains acids, destructive enzymes (including nucleases), and unique secondary compounds (products produced from pathways that are not part of primary metabolism) that potentially damage DNA. To minimize contaminants from the vacuolar contents, salts and other inhibitors have been added to the lysis buffer.

Determining the Quality & Quantity of Nucleic Acids

Once nucleic acid has been isolated, it is essential to analyze the quality of the preparation. The stability of the nucleic acid and its performance in subsequent enzymatic steps is affected by its purity and how much you use is determined by its concentration. Both the amount of nucleic acid isolated and its purity are affected by the type of tissue that it is isolated from, the amount of tissue used, and the isolation technique used.

Three general techniques, which are typically used in labs to evaluate the concentration and purity of an isolated nucleic acid such as genomic DNA:

Spectroscopic analysis of ultraviolet absorption at 260 nm will provide total nucleic acid concentration, and a ratio of A_{260}/A_{280} (contaminating protein) and A_{260}/A_{230} (contaminating carbohydrates) will give you an idea of quality. The NanoDrop, a low-volume spectrophotometer, is most commonly used as it requires microliters of sample volume.

The concentration of a nucleic acid sample (either RNA or DNA) can be determined using UV spectrophotometry. Both RNA and DNA absorb UV radiation at 260nm, thanks to the nitrogenous bases, making it possible to detect and quantify at concentrations as low as 2.5ng/ μ l. You can calculate the concentration of the DNA in your sample as follows:

DNA concentration (μ g/ml) = (OD 260) x (dilution factor) x (50 μ g DNA/ml)/(1 OD260 unit)

In contrast to nucleic acids, proteins have a UV absorption maximum of 280 nm, due mostly to the tryptophan residues. Purity determination of contaminants can be recognized by the calculation of a “ratio” of these two absorbance readings. The ratio A_{260}/A_{280} is used to estimate the purity of nucleic acid since proteins absorb at 280 nm. Pure DNA should have a ratio of approximately 1.8, whereas pure RNA should give a value of roughly 2.1. Absorption at 230 nm reflects biological contaminants of the sample such as carbohydrates, peptides, salts, or proteins. Also, many chemicals commonly used in nucleic acid preparations, such as phenol, EDTA, and SDS, can be detected by their absorbance at 230 nm.

Qualitative analysis using **agarose gel electrophoresis** is one of the most common methods that provide valuable information about size, quality, and relative concentration.



Let's Explore

Watch this video on DNA extraction before coming to class: <https://youtu.be/90Gsxal57-g>

Draw a brief flowchart of DNA isolation from plants below:

PREPARING FOR AN EXPERIMENT

Starting from this lab exercise, you will learn to create your own experimental design. There are several important aspects to a good laboratory protocol. The first is understanding the purpose of your experiment. The purpose statement sets up the goal of the experiment, and most of the time answers an experimental question.

What is the experimental question in this (DNA Barcoding) experiment? (Write your ideas below)

Writing a Purpose Statement. For each lab, read the entire laboratory exercise and write a short, 2-3 sentence purpose statement. It should begin “The purpose of this lab...” and should include the experiment of the lab exercise. For example, in one of the lab exercises, you will isolate your own genomic DNA to perform a VNTR analysis. “*The purpose of this lab is to perform a VNTR analysis of human genomic DNA.*” Although you learn how to use PCR to perform this analysis, the purpose is not to learn how to perform PCR.

Your Turn! Create a Purpose Statement for this DNA Barcoding lab exercise:

Share this statement with your lab partner and provide each other editing advice.

Another important aspect of an experimental protocol is to understand any safety hazards that may arise during the experiment. This may include safe handling of a chemical such as what PPE do you wear or using a chemical fume hood. It may also include the safe disposal of reagents, such as autoclaving biohazard material.

Writing a Hazard Communication Statement. After you have read the entire laboratory exercise, provide a comprehensive list of the potential hazards (chemical, biological, physical...) that you may be exposed to during this exercise. List precautions that you will take for each hazard. For example, if you use Hydrochloric Acid, you may state the following: “*Hydrochloric Acid (HCl) DANGER: Corrosive. Avoid contact with skin and eyes, Avoid inhalation of fumes and mist. Do not mix with caustics or other reactive materials. Wear PPE (gloves, protective eyewear, and close-toed shoes)*”

Your Turn! Does this laboratory exercise have any special hazards? List all below:

Writing a Procedure. The lab handouts include a lot of background material and other information in the procedural steps for your instruction in these techniques. Your procedure should not include this type of information and should be limited only to the actual steps taken in a procedure without explanation. You should read the instructions in your manual and extract just the action required of you during lab. Your procedure should be numbered steps. Thus, you will create a document that is easier to follow during the lab session. So far you've been writing small procedural protocols. Starting in lab 7, you will begin to write your own based on the instructions provided in the manual.

The composition of the protocol is a skill that you must master. It is sometimes difficult to gauge the amount of detail that a protocol needs. A protocol that is too long and detailed is too cumbersome to use routinely, while one lacking sufficient detail will not lead to uniformity when different people perform the procedures. In this course, we will guide you through these decisions by providing you with a lab protocol to follow. In general, a protocol that needs the most detailed information is used by a large number of people, is infrequently used so that the users will not remember exactly how it is done and involves especially sensitive or critical steps of a process.

REAGENTS, SUPPLIES, & EQUIPMENT

Creating a comprehensive list of equipment and reagents will not only save you time but help you organize your experiment. Your list should be comprehensive enough that it will take you 5 minutes to collect your items and not have to return searching for items. Being prepared for your experiment will ensure experimental success.

<p>Each person:</p> <ul style="list-style-type: none"> • Container with crushed ice • DNA rehydration solution 100µL) • 70% ethanol (600ul) • Isopropanol (600ul) • 2 microcentrifuge tubes (1.5 ml) • Nuclei lysis solution (600ul) on ice • Protein Precipitation solution (200ul) • Plastic pestles • Tissue specimen • Microcentrifuge tube rack 	<p>Share:</p> <ul style="list-style-type: none"> • Microcentrifuge • Picofuge • Water bath or heating block at 65°C • Vortex • Micropipettes and barrier tips
--	--

PROTOCOL

Protocol	Experimental Notes
<p>PART I: DNA EXTRACTION PROTOCOL</p> <ol style="list-style-type: none">1. Obtain plant or animal tissue ~10-20 mg or ¼ inch diameter from your sample.2. Place sample in a clean 1.5 mL tube labeled with your initials.3. Add 100 µL of nuclei lysis solution to the tube.4. Twist a clean plastic pestle against the inner surface of 1.5 mL tube to forcefully grind the tissue for 1 minute. Use a clean pestle for each tube if you are doing more than one sample.5. Add 500 µL more nuclei lysis solution to the tube and mix by vortexing.6. Incubate the tube in a water bath or heat block at 65°C for 15 minutes. Then stand tube at room temperature for 5 minutes.7. Add 200 µL of protein precipitation solution to each tube. Vortex tubes for 5 seconds and incubate on ice for 5 minutes.8. Place your tube and those of other groups in a balanced configuration in a microcentrifuge, with cap hinges pointing outward. Centrifuge for 4 minutes at maximum speed to pellet protein and cell debris.9. Label a clean 1.5 mL tube with your initials. Use a fresh tip to transfer 600 µl of supernatant to the clean tube. Be careful not to disturb the pelleted debris. Discard old tube containing the pelleted debris.10. Add 600 µL of isopropanol to the supernatant in the tube. Close cap, and mix by rapidly inverting tubes several times.11. Centrifuge for 1 minute at maximum speed to pellet the DNA.	<p>Name of Plant: _____</p> <p>Tissue Sampled: _____</p> <p>Mass Sample: _____ g</p>

12. Carefully pour off the supernatant from the tube, and add 600 μL of 70% ethanol. Close cap, and invert tube several times to *“wash” the pellet*.
13. Centrifuge the tube for 1 minute at maximum speed.
14. Quickly, and carefully pour off the supernatant and pop-spin in picofuge for 3 seconds. Use a 200ul micropipette with a fresh tip to remove any remaining ethanol, being careful not to disturb the pellet.
15. Air-dry the pellet for 5-10 minutes to evaporate remaining ethanol. Do not air dry longer or resuspension will be difficult.
16. Add 100 μL of the DNA rehydration solution to each tube, and dissolve the DNA pellet by pipetting up and down several times. Incubate the DNA at 65°C for 30 minutes, or overnight at 4°C *if* the DNA does not easily resuspend.
17. Store your sample on ice or at -20°C if continuing another day.

Time: _____ min

Temperature: _____ degC

Time: _____ min

Stored: _____

PART II: SPECTROPHOTOMETRIC ANALYSIS OF DNA

Determine the concentration of your DNA on the low-volume (NanoDrop) Spectrophotometer

1. Using the SOP provided in your SOP booklet set the NanoDrop spectrophotometer to read DNA (nucleic acid).
2. Blank the NanoDrop with 1.5 μL of fresh rehydration solution (not your heated aliquot). Why not use water? _____
3. Measure 1.5 μL of your DNA preparation. Record data.
4. Obtain a freezer box for the whole class to store DNA. Label the box with the class name and instructor name. Ensure your tubes are clearly labeled. If you do not move on to part C, store the freezer box at -20°C until next class.

DNA concentration:
_____ ng/ μL

A260/A280: _____





Lab Unit 4-C: Agarose Gel Electrophoresis of Genomic DNA

Agarose gel electrophoresis is a molecular biology technique used to separate a mixture of nucleic acid fragments according to size in an agarose matrix. DNA or RNA samples are separated by applying an electric field to move the largely negatively charged molecules toward the positive electrode. The separated nucleic acids can be visualized using a nucleic acid binding stain and exposed to UV light.

Agarose. The gel matrix is composed of agarose molecules in supercoiled bundles that are aggregated into three-dimensional structures with channels and pores through which the nucleic acid fragments move. The large pore size and gel strength make agarose suitable for nucleic acid molecules above 50 base pairs. The size of the nucleic acid fragments will dictate which pore size is most appropriate; large fragments are best separated on 0.8% gel, 1000-3000bp size fragments on 1%, and small fragments (100-500bp) on 2-3% gels. The speed of migration through the matrix is also proportional to the voltage applied, the percentage of gel and size of the fragment. For mini-gel analysis, the voltage is typically set to 5V/cm between the electrodes. A mini-gel is typical run around 100V.

Because DNA molecules are double-stranded (and the same shape), they migrate through the gel matrix separating according to size alone. Smaller molecules travel faster through the gel matrix than larger molecules at a rate that is inversely proportional to the log₁₀ of the number of base pairs. Therefore, by running a known molecular weight ladder at the same time as your nucleic acid sample, you can determine the size of the nucleic acid fragments using a linear regression analysis.

SYBR Green stain. SYBR Green is a very suitable dye for detecting nucleic acid fragments in agarose gels, as it binds nucleic acids and emits fluorescent radiation on UV illumination. Because SYBR green is relatively non-toxic, it can be disposed of down the sink or in the regular trash; however, it is always a good practice to wear nitrile gloves and safety glasses when working with any chemical.

Gel Buffers. Two buffers are commonly used when preparing agarose gels for separating DNA fragment; TAE (with sodium acetate) and TBE buffers (with boric acid). Both have benefits and limitation in use. Both buffers are prepared from 50X stocks, which can be stored at room temperature.



Let's Explore

Learn more about gel electrophoresis, here!! <http://learn.genetics.utah.edu/content/labs/gel/>

PRE-LAB EXERCISE: PREPARE AN SOP OF 50XTAE AGAROSE GEL RUNNING BUFFER:

Review from your lecture: *A Standard Operating Procedure (SOP) is used for operations that are standard, meaning, laboratory tasks that you repeat the same way and expect the same outcome no matter who is performing the procedure, when or how many times.* These are common in every workplace, but especially regulated biotechnology labs such as biomanufacturing departments.

COMPONENTS OF A STANDARD OPERATING PROCEDURE (SOP)

1. Title
2. ID number, approval date, revision number, date of revision
3. Signatures of preparers/approvers.
4. Statement of Purpose
5. Scope describes when the procedure is relevant
6. A Statement of Responsibility describes who performs task
7. Safety Statement. Any known hazards?
8. Materials. Lists the essential items necessary to perform the procedure.
9. Calculations required. Leave space with an example
10. Procedure. Details the tasks and activities required to successfully perform or conduct the operation described, written in numbered steps.
11. References to other documents/forms.
12. Instructions on how to document that the procedure was performed

Each group must prepare their own **100 mL bottle of the 50XTAE buffer** that they will use throughout the semester. Perform calculations below and verify with your lab partner. Include a calculations section in your SOP. You may hand-write them into your SOP.

The recipe for 50XTAE is as follows:

- 2 M Tris base (molecular wt = 121.14 g/mol)
- 0.05 M EDTA (from a 0.5 M stock provided)
- 5.72% glacial acetic acid (100% acetic acid provided in fume hood – caution!)

PART I: PREPARE 50X TAE STOCK BUFFER

1. Use the template provided by your instructor and create an SOP for preparing 100mL of 50X TAE using the reagents listed above. You may work with a lab partner, but generate an SOP to submit with your lab report. *show your calculations* You may write them by hand if typing them is too complicated.
2. Exchange your 50X TAE SOP with your lab partner (or another lab group) for review. Check proper format and calculations. Make comments directly on the SOP. You will submit this annotated copy with your lab report as well as a final copy to reflect any corrections needed.
3. Print out a clean final copy of your SOP. Use this SOP to prepare a 100mL 50XTAE solution.
4. Label your reagent according to the labeling SOP and ask your instructor where it will be stored for the semester.
5. Fill in a solution preparation form to include with your report.
6. Prepare an inventory control form for your 50XTAE solution. Keep all inventory control forms with the solutions (your instructor may provide a binder for this). Every time you use an aliquot from this bottle, you must document and log it on the inventory control form.

PART II: PREPARE 1% AGAROSE GEL WITH SYBRsafe

Materials required:


Agarose	Balance, weigh boats
1X gel running buffer	SYBR Safe stock solution, 10,000X
Gel-casting apparatus	imaging system
Gel electrophoresis box	Erlenmeyer flask
Power supply	Saran Wrap
Microwave oven or hot plate	Micropipettes and tips
Molecular weight DNA ladder (1kB DNA ladder, Promega G7541)	

1. Before you begin – check all calculations with your lab partner and other lab groups.
2. Set up a gel casting apparatus as directed by your instructor.
3. **Prepare 400 mL of 1X TAE buffer. Calculation:** How will you prepare 400mL of 1XTAE from the 50XTAE stock you made?

Protocol: Write a step-wise procedure to prepare 400mL of 1XTAE buffer from a 50X Stock.


4. Weigh out _____ g of agarose and add it to a 125 mL Erlenmeyer flask.

Calculation: How much agarose will you need to make 30mL of a 1% agarose gel?


5. Add 30mL of 1X gel running buffer in the Erlenmeyer flask and place the flask on a level surface and carefully mark the glass at the fluid level with a permanent marker (DO NOT mark the white marking area with a permanent marker! It cannot be removed). Insert a Kim wipe in the top to prevent evaporation.
6. Heat the mixture until all agarose has dissolved using a microwave. Typically, 1 min on med-high works well. Swirl until all particles are dissolved.
 CAUTION! The flask will be hot! Use hot hands to remove flask and swirl. Do not hold over the top of your face while swirling.
7. Observe the fluid level about the mark you made on the flask. If a significant amount of water has evaporated, carefully add water to return to the level of the mark and swirl the solution.
8. Cool the solution to 50 – 60°C while continuing to swirl until it is cool to the touch. Do not cool too long or it will solidify in the flask.
9. Add _____ ul of SYBR Safe stock solution (10,000X) to a final concentration of 1X.

Calculation: How much SYBRsafe will you add to your flask to achieve a final concentration of 1X starting with a 10,000X stock?

10. Pour the gel immediately into your casting tray. Allow the gel to form completely; typically, 20 minutes at room temperature is sufficient.



NOTE: When you are done pouring the molten agar into the gel cast tray, allow any molten agar residue to solidify in the glassware for 10 min. Remove excess solidified agarose with a spatula and dispose of in the trash can, & place glassware in the dirty dishes bin.



NOTE: If you make a mistake with this solution DO NOT dump down the sink. The powdered and molten agarose will solidify in the sink and clog it. Ask your instructor where to dispose of the waste.

11. Carefully, remove the comb and bumpers from the gel, place the gel in the electrophoresis chamber, and add a sufficient volume 1X TAE gel running buffer to completely cover the surface of the gel. This is usually the remaining volume of buffer.
12. Add 5ul of 5-6X load dye to a clean 1.5mL centrifuge tube. Add 10-20ul (0.5-1ug) of genomic DNA to this tube and flick to mix and pop spin in a picofuge.
13. The molecular weight marker should already be prepared for you. Collect an aliquot (6ul) to run on the gel with your samples.
14. Load 6ul of the molecular weight marker in lane 1, and all (25ul) of the DNA sample into lane 2 (don't skip lanes).
15. Place the cover on the electrophoresis chamber in the correct orientation and connect leads to a power supply. Set the power supply to approximately 80 volts and allow to electrophorese until the blue tracking dye is about 3/4th the way to the bottom of the gel. This may take almost one hour!
16. Fill out a gel electrophoresis documentation form to include with your report.
17. Capture an image using the gel documentation system. Affix a copy of your gel image to the gel electrophoresis documentation form and turn this in with your lab report.
18. Save the remaining genomic DNA in a labeled microcentrifuge tube, return to freezer box, and store at -20°C.



4. We will continue the DNA barcoding experiment in future lab exercises. Briefly, outline the remaining steps to DNA barcoding using a flowchart.

5. Write a conclusion statement for your experiment:

LAB UNIT 5: ENZYME-LINKED IMMUNOSORBENT ASSAY (ELISA)

Jack O'Grady, M.S.

This lab is based on Bio-Rad's ELISA ImmunoExplorer Kit (Cat#166-2400EDU). The introduction, instructions, and analysis is copied in part or whole from the instruction manual.

OBJECTIVES

Your performance will be satisfactory when you can:

- ✓ Understand the mechanism for an ELISA assay
- ✓ Perform an ELISA
- ✓ Graph data using MS Excel, perform linear regression
- ✓ Utilize an ELISA standard curve to determine the concentration of an unknown protein sample

LET'S GET STARTED!

Reading:

Chapter 5: Introduction to studying protein

Chapter 6: Identifying a potential Biotechnology product

Chapter 7: Spectrophotometers and Concentration Assays

INTRODUCTION

Immunology is the study of the immune system and how the body protects itself against disease. Over 100 years ago, biologists found that animals' internal immune systems respond to invasion by "foreign entities" or antigens. When an invader enters the body, it provokes an immune response that begins with the production of proteins called antibodies.

Like magic bullets, antibodies seek out and attach themselves to invading entities (antigens), flagging the invaders for destruction by other cells of the immune system. The antigenic invaders may be any molecules foreign to the body, including components of infectious agents like bacteria, viruses, and fungi. Today, antibodies have become vital scientific tools, used in biotechnology research and to diagnose and treat disease. The number of different antibodies circulating in the blood has been estimated to be between 10^6 and 10^{11} , so there is usually an antibody ready to deal with any antigen. In fact, antibodies make up to 15% of your total blood serum protein. Antibodies are very specific; each antibody recognizes only a single antigen.

You are about to perform an ELISA (enzyme-linked immunosorbent assay). The ELISA relies on antibodies to detect the presence of antigens in liquid samples. Because they are antibody-based, ELISAs are called immunoassays. ELISAs can detect minute amounts of disease agents in samples such as body fluids (before the body has had a chance to mount an immune response). Smallpox virus is an example of a disease agent that can now be detected using an ELISA. If exposure is detected and treated with a vaccine within 2–3 days, patients do not develop




smallpox. Other applications for ELISA include testing for West Nile virus, HIV coat protein p24, SARS virus, anthrax spores, hormones such as hCG in pregnancy tests, illegal steroids in drug tests, bacteria in food safety tests, and the presence of genetically modified organisms contaminating non-GMO food.

Where Is ELISA Used in the Real World?


With its rapid test results, the ELISA has had a significant impact on many aspects of medicine and agriculture. ELISA is used for such diverse purposes as home pregnancy tests, disease detection in people, animals, and plants, detecting illegal drug use, testing indoor air quality, and determining if food is labeled accurately. For new and emerging diseases like severe acute respiratory syndrome (SARS), one of the highest priorities of the US Centers for Disease Control and Prevention (CDC) and the World Health Organization (WHO) has been to develop an ELISA that can quickly and easily verify whether patients have been exposed to the virus. Over-the-counter kits that are based on the same principles as this ELISA activity include home pregnancy and ovulation tests, and tests for the presence of illegal drugs like marijuana and cocaine.

Some tests give positive or negative results in a matter of minutes. For example, home pregnancy dipstick tests detect levels of human chorionic gonadotropin (hCG), a hormone that appears in the blood and urine of pregnant women within days of fertilization. The wick area of the dipstick is coated with an anti-hCG antibody labeled with a pink compound (step 1). When the strip is dipped in urine, if hCG is present it will bind to the pink antibody, and the pink hCG-antibody complex will migrate up the strip via capillary action (step 2). When the pink complex reaches the first test zone, a narrow strip containing an unlabeled fixed anti-hCG antibody, the complex will bind and concentrate there, making a pink stripe (step 3). The dipsticks have a built-in control zone containing an unlabeled secondary antibody that binds unbound pink complex (present in both positive and negative results) in the second stripe (step 4). Thus, every valid test will give a second pink stripe, but only a positive pregnancy test will give two pink stripes.



How Are Antibodies Made?

When exposed to antigens, all mammals generate an immune response and produce antibodies, proteins that recognize and bind tightly to the specific antigens. Each antibody recognizes only a single antigen. Animals such as goats, rabbits, and mice can be injected with an antigen and, after a period, their serum will contain antibodies that specifically recognize that antigen. If the antigen was a disease-causing agent, the antibodies could be used to develop diagnostic tests for the disease. In an immunoassay, the antibodies used to recognize antigens like disease agents are called primary antibodies.



A) Structure of IgG bound to the HIV capsid protein p24 as determined by X-ray crystallography (Harris et al. 1998, Momany et al. 1996). These structures can be downloaded from the Protein Data Bank (www.pdb.umd.edu, (Berman et al. 2000) using the PDB identification codes 1IGY and 1AFV and manipulated using free online software such as Rasmol and Protein Explorer. **B)** A commonly used representation of an antibody bound to an antigen.

Secondary antibodies recognize and bind to primary antibodies in an immunoassay. They are prepared by injecting antibodies produced by one species of animal into another species. This works because the antibodies produced by different species are different enough from each other that they will provoke an immune response. For example, if you want a secondary antibody that will recognize a primary human antibody, inject human antibodies into an animal like a rabbit. After the rabbit immune response, the rabbit serum will contain antibodies that recognize and bind to human antibodies. Secondary antibodies are frequently labeled to make them visible.

In this experiment, the secondary antibodies you will be working with are conjugated to an enzyme named horseradish peroxidase (HRP); HRP in the presence of its substrate, TMB, produces a blue color.

Controls in Immunoassays

For an immunoassay to be valid, it must include both positive and negative controls, i.e., samples that will give known results. Controls are always run side by side with experimental samples. If you do not run a positive control and the experiment provides negative results, how can you be sure the results are truly negative? What if the assay simply did not work? If a positive sample gives a negative assay result, it is called a false negative.

Conversely, if you do not run a negative control and the experiment gives all positive results, how can you be sure the results are truly positive? What if the assay was contaminated with antigen? If a negative sample gives a positive assay result, it is called a false positive.


Controls are also needed to guard against experimental error and to ensure that the assay is working correctly. There can be problems with reagents, which can degrade due to age or poor storage conditions. Operators can make mistakes by choosing the wrong reagents, causing errors in dilutions or in pipetting, or failing to remove unbound reagents. Poor record keeping is another source of false assay results. Most of these possibilities can be checked for within the assay with the appropriate controls. Now let's put this all together.

The main steps in this antigen detection ELISA are:

1. Add your sample and control samples to the wells in a microplate strip. Your samples contain many proteins and may or may not contain the antigen. Incubate for 5 minutes to allow all the proteins in the samples to bind to the plastic wells via hydrophobic interaction. This is called an immunosorbent assay because proteins adsorb (bind) to the plastic wells.



2. Add primary antibody to the wells and incubate. The antibodies will seek out the antigen from the many proteins bound to the well. If your sample contains the antigen, the antibodies will bind it tightly and remain in the well.



3. Detect the bound antibodies with HRP-labeled secondary antibody. If the primary antibodies have bound to the antigen, the secondary antibodies will bind tightly to the primary antibodies.




4. Add enzyme substrate to the wells, wait 5 minutes, and evaluate the assay results. If the antigen was present in your sample, the wells will turn blue. This is a positive diagnosis. If the wells remain colorless, the antigen was not present in your sample and the diagnosis is negative.



Like most ELISA assays, this kit relies on a Horseradish Peroxidase (HRP) conjugated antibody and the TMB (3, 3', 5, 5'-tetramethylbenzidine) substrate. TMB is a chromogen that yields a blue color when oxidized with hydrogen peroxide (catalyzed by HRP) that has significant absorbance peaks at 370 nm and 652 nm. The color then changes to yellow with the addition of acid with

maximum absorbance at 450 nm. The relative amount of protein in the well will be directly proportional to the amount of signal that is obtained at 450 nm.

This kit uses an indirect ELISA format, in which a sample is added to a microtiter plate well. After coating the well with the sample, the unbound sample is removed, the wells are washed and a primary antibody to the target protein is added. Unbound primary antibody is washed away and a secondary antibody, which recognizes this primary antibody, binds. The secondary antibody is conjugated to HRP, is added (the “detection antibody”). The signal is generated by reaction with the TMB substrate as described above. The intensity of the signal (measured at 450 nm) is directly proportional to the amount of target protein in the sample. Dilutions of a standard are used to construct a standard curve, from which the concentration of target protein in the samples are determined by interpolation.



ELISA antibodies can recognize antigens directly or by recognizing another antibody that recognizes the antigen

QUANTITATIVE ELISA

A set of known standards is analyzed along with the samples at the same time using the same reagents. The range of known standard concentration is determined by both target antigen concentration as well as dynamic range of ELISA kit used.

A standard curve is constructed by plotting the absorbance obtained from each reference standard against its concentration in ng/mL, as shown below. The following figure is a typical antigen standard curve. Using the equation generated by the standard curve, the target antigen concentration (ng/mL) of an unknown sample can be determined.

To determine the total amount of target antigen in the sample, use the following equation and solve for the x value:

$$y = mx + b, \text{ rearrange: } x = (y-b)/m$$

Where:


b=the y-axis intercept of your standard curve

m=the slope of your standard curve

x=the antigen concentration

y=the 450 nm Absorbance of the sample

For example: If the OD reading of the sample is: 1.248, using the above equation and the Standard Curve shown above the antigen concentration would be: 8.342 ng/mL




Sample Standard curve for ELISA

ELISA DIAGNOSTIC TEST FOR CANCER DETECTION

Despite advances in detection and treatment, cancer remains a deadly disease with more than 1.63 million new cases diagnosed each year in the United States (American Cancer Society, 2012), and approximately 7.6 million people die worldwide (Center for Disease Control, 2013)

Early cancer detection is key to decreased total numbers of death from cancer. Protein biomarkers are showing promise in aiding early detection system. Carcinoembryonic Antigen (CEA) is one such antigen that has shown promise in assisting early detection, disease progression monitoring as well as determining response to pharmacologic intervention. The figure below is data reproduced from Bayer Diagnostics/Siemens white paper.



Any number of conditions may cause an increase in CEA antigen levels in serum. Upregulation of CEA may be an indication of many potential problems. Additional diagnostic testing helps pinpoint the actual problem. *CEA serum levels for healthy nonsmokers is less than 3ng/mL. The CEA level for smokers is elevated and puts them at risk of cancer at up to 15 ng/mL. With cancer patients, CEA levels are significantly elevated in serum and usually, exceed 20 ng/mL.*

In this simulation serum ELISA lab we will determine the CEA levels of two patients ("A" and "B"). You will run a set of CEA antigen standards and determine the concentration of your patients from a standard curve.

References:

1. MaxDiscovery™ GAPDH ELISA Kit Manual – 3401-01 from ©BIOO Scientific Corp
2. Elynn Daugherty, Biotechnology: Science for the New Millennium. 2012. EMC Paradigm Publications. ISBN: 978-0-76384-284-0
3. Bio-Rad's ELISA Immuno Explorer Kit (Cat#166-2400EDU)



Let's Explore

Watch the following video on using ELISA to detect and quantitate specific antigens. Draw a single well to show how Indirect ELISA works. <https://youtu.be/RRbuz3VQ100>

TITLE: *(write a fully descriptive title to your experiment – Write this here and on your title page!)*

PURPOSE: *(write a purpose statement to your experiment)*

SAFETY: *(review the procedure below. Write a safety statement your experiment)*

MATERIALS: Review the protocol below and devise a detailed list of reagents and equipment you will need to complete this ELISA assay. Verify your list with your lab partner!

PROCEDURE: Watch this video <https://youtu.be/849HN1ueUhs> and draw a simple flowchart below summarizing the ELISA procedure you will be performing:


PROTOCOL

NOTE: Check off each step as you complete them! Write down observations as you move through the experiment.

PART I: ELISA ASSAY:


You will perform the ELISA ASSAY on BOTH strip plates at the same time! Why?

1. Label the outside wall of each well on one 12-well strip with the numbers 1–12. Label the first three wells of a second 12-well strip with a "+" for the positive controls, the next three wells with a "-" for the negative controls, the next three wells with the initials of one of your sample tubes, and the last three wells with the initials of the second sample tube.



2. Use a pipet to add 50 μ l of PBS from the yellow tube labeled "PBS" to wells labeled #2 through #12.
3. Add 100 μ l from the yellow tube labeled "1,000 ng/ml AG" to the well labeled #1.
4. Perform serial dilution from well #1 through well #11 in the following manner:
 - a. Pipet 50 μ l out of well #1 and add it to well #2. Pipet up and down gently three times to mix the sample in well #2.
 - b. Using the same pipet tip, transfer 50 μ l from well #2 to well #3 and mix the sample in well #3.
 - c. Using the same pipet tip, transfer 50 μ l from well #3 to well #4 and mix the sample in well #4.
 - d. Repeat this transfer and mixing step, moving to the next well each time. STOP when you reach well #11; discard the 50 μ l of solution from well #11 into a waste container.
5. In the second microplate strip, use a fresh pipet tip to transfer 50 μ l of the positive control (+) from the violet tube into the three "+" wells.
6. Use a fresh pipet tip to transfer 50 μ l of the negative control (-) from the blue tube into the three "-" wells.
7. Use a fresh pipet tip to transfer 50 μ l of each of your team's samples into the appropriately initialed three wells.


4. Wait 5 minutes while all the proteins in the samples bind to the plastic wells.
5. Wash the unbound sample out of the wells:
 - a. Tip the microplate strip upside down onto the paper towels so that the samples drain out, then gently tap the strip a few times upside down on the paper towels. Make sure to avoid splashing sample back into wells.



- b. Discard the top paper towel.
- c. Use a transfer pipet filled with wash buffer from the beaker to fill each well with wash buffer taking care not to spill over into neighboring wells. The same transfer pipet will be used for all washing steps..




- d. Tip the microplate strip upside down onto the paper towels so that the wash buffer drains out, then gently tap the strip a few times upside down on the paper towels.
 - e. Discard the top 2–3 paper towels.
6. Repeat wash step 5.
 7. Use a fresh pipet tip to transfer 50 μ l of primary antibody (PA) from the green tube into all 12 wells of the microplate strip.



8. Wait 5 minutes for the primary antibody to bind.
9. Wash the unbound primary antibody out of the wells by repeating wash step 5 **two** times.



10. Use a fresh pipet tip to transfer 50 μ l of secondary antibody (SA) from the orange tube into all 12 wells of the microplate strip.




11. Wait 5 minutes for the secondary antibody to bind.
12. Wash the unbound secondary antibody out of the wells by repeating wash step 4 **three** times.



The secondary antibody is attached to an enzyme (HRP) that chemically changes TMB (the enzyme substrate), turning it from a colorless solution to a blue solution. Predict which wells of your experiment should turn blue and which should remain colorless and which wells you are not sure about.

13. Use a fresh pipet tip to transfer 50 μ l of enzyme substrate (SUB) from the brown tube into all 12 wells of the microplate strip.



6. After 5 min, add 50 μ l of 0.18M sulfuric acid. This will turn the blue wells to a bright yellow color.
7. Analyze in a plate reader set to 450nm.
8. Using the Excel spreadsheet of your data, determine the CEA concentration of your patient samples A & B and assess the risk of cancer for each patient by comparing your results to the literature provided in the introduction.

Part II: Quantitative ELISA Analysis

Standard Concentration (ng/mL)	Absorbance (450nm)	Absorbance – Blank Absorbance
0		
1		
2		
4		
8		
16		
31		
63		
125		
250		
500		
1000		
Negative control		
Patient A		
Patient B		

- Using MS Excel, graph your standard curve. Plot (absorbance-blank) on the y-axis and CEA concentration on the x-axis. Label both axes (with units) and give the graph and appropriate descriptive title.
- Using MS Excel, determine the equation of the line *for the linear points of the curve* as follows. Note the graph may have a sigmoidal (“S”) shaped curve. Notice and high and low points on the graph that may not be in the linear range and remove them. Keep the “no analyte” (0, 0) point.
 - Input data into MS Excel in two columns, with the CEA concentration on the left column, and the corresponding absorbance (-blank) on the right column.
 - Select data, click on insert, click on scatter plot
 - Right click on one of the points, select, add a trend line
 - Click on Linear, set intercept at 0,0, click on display equation on chart and display R² value on chart
- The R² value is a correlation coefficient that will tell you how well your data correlates linearly to each other. A R² value >0.95 is acceptable, but >0.98 is preferred.
- Click on the chart elements (Right corner of graph box, “+” and select axis, axis title, and chart title. And complete the corresponding titles. Print out a copy of your graph and include with your report.
- Using the equation of the line, determine the CEA concentration of your sample.



LAB UNIT 6: RECOMBINANT DNA TECHNOLOGY

Jack O'Grady, M.S., Linnea Fletcher, Ph.D., Angela Wheeler, M.S., Patricia Phelps, Ph.D.

INTRODUCTION TO RECOMBINANT DNA TECHNOLOGY:

Recombinant DNA technology has dramatically impacted and advanced numerous fields – from medicine to agriculture. Recombinant DNA technologies have helped create genetically altered plants, which can withstand different environmental conditions and reduce the use of pesticides to yield crops that are more abundant. They have also helped produce advances in medicine, such as treatments for cancer, recombinant therapeutic protein drugs, and vaccines. One interesting use of recombinant DNA technology is in the field of animal husbandry. Recombinant DNA technologies about animal husbandry have led to the development of transgenic animals and clones, including the first successfully cloned animal, Dolly, the sheep. To learn more, visit Khan Academy: <https://www.khanacademy.org/science/biology/biotech-dna-technology>

Lab Unit 6-A: Bioinformatics of Green Fluorescent Protein

OBJECTIVES

Your performance will be satisfactory when you can:

- ✓ Locate scientific publications and biological databases of DNA and protein sequences on the Internet (GenBank, BLAST, NEBcutter, ClustalQ)
- ✓ Retrieve and compare sequence information from databases – BLAST, Jmol, ClustalW
- ✓ Compare evolutionary relatedness and draw phylogenetic trees from sequence comparisons
- ✓ Download the GenBank submission for pGLO plasmid
- ✓ Use NEBcutter to generate a restriction map of a plasmid
- ✓ Use BLAST to identify pGLO proteins and use Jmol to explore the 3-D structure of GFP
- ✓ Use Clustal software to align different jellyfish species the GFP gene, and generate a phylogram.

INTRODUCTION TO BIOINFORMATICS

The **Human Genome Project** (HGP), completed in 2003, was a 13-year project coordinated by the U.S. Department of Energy and the National Institutes of Health. During the early years of the HGP, the Wellcome Trust (U.K.) became a primary partner; additional contributions came from Japan, France, Germany, China, and others. For more information:

http://www.ornl.gov/sci/techresources/Human_Genome/project/hgp.shtml

Human Genome Project goals were to:

1. identify all the approximately 20,000-25,000 genes in human DNA,
2. determine the sequences of the 3 billion chemical base pairs that make up human DNA,
3. store this information in databases,
4. improve tools for data analysis,
5. transfer related technologies to the private sector, and
6. address the ethical, legal, and social issues that may arise from the project



Though the HGP is now complete, analyses of the data will continue for many years. Follow this ongoing research on the HGP page:

http://www.ornl.gov/sci/techresources/Human_Genome/project/timeline.shtml

An essential feature of the HGP project was the federal government's long-standing dedication to the transfer of technology to the private sector. By licensing technologies to private companies and awarding grants for innovative research, the project catalyzed the multibillion-dollar U.S. biotechnology industry and fostered the development of new medical applications.

Another goal of the Human Genome Project is to sequence the genomes of other species of interest, such as model organisms used by biologists, pathogens of medical importance, and crop plants of agricultural importance. This goal has also been exceeded, and today the genomes of tens of thousands of species have been sequenced. Check out the list here, at the NCBI website: <https://www.ncbi.nlm.nih.gov/genome/browse/#!/overview/>

The Human Genome Project has been a catalyst for change in the way biologists approach the study of living things. Biologists today using the sophisticated laboratory technology for sequencing DNA are collecting data faster than they can interpret it. A new field called **bioinformatics** is developing for the storage and management of the data stored in these rapidly growing databases, as well as for the use of a computer as a general tool for discovering how living things work.

The power of bioinformatics approach for the discovery of genes has been proven with the completion of the yeast genomic sequencing project in 1996. Once the genome was fully sequenced, bioinformatics approach could be used to scan for and identify genes. The genes discovered this way could be compared to a large number of genes that had already been discovered through more classical molecular and genetic techniques. The results were remarkable. Before the yeast genome was sequenced in 1996, an international collaboration of scientists studying the genetics of this model organism had identified an impressive 2,000 genes by conventional genetic analysis. When the yeast genome sequencing was completed, bioinformatics searches for similarities of DNA sequences from other organisms were able to locate an additional 2,000 genes. Meaning, in less than one year, a single laboratory using DNA sequencing and computer searches of sequence data could both duplicate and double the gene discovery of a 20-year international effort.

Once a gene has been identified, many new questions can be asked: what kind of protein does it code for and what is its function? How does it interact with other molecules of the cell? Is it expressed at all times as a so-called “housekeeping gene,” or is it a developmentally regulated gene? Is the expression tissue-specific? Is it expressed in response to an environmental factor? These questions are the same questions that have been asked by molecular and cell biologists for decades, usually by studying one gene, and its protein or proteins, at a time. With the copious amounts of information coming from the genomics project, however, biologists can ask the same questions about more complex systems. Instead of asking about one protein at a time, biologists can now ask questions about hundreds of proteins at a time, looking for patterns of



structure and patterns of expression. Looking at the proteins on a genomic scale is a new field now called **proteomics**. When a new gene has a sequence that has been found to be homologous to a gene in a database that has already been characterized, sometimes many of these questions about protein structure and function can be answered quickly by the bioinformatics approach. For example, the 2,000 new genes discovered by the yeast genomic sequencing project, discussed above, matched genes of other organisms whose function had already been determined.

Bioinformatics approach is playing an increasingly important role in protein structure studies. Although the final structure of a protein is determined by the amino acid sequence of that protein, we have yet to model the correct folding of a protein by its amino acid sequence by computer. There is progress, however, in achieving this so-called “holy grail” of proteomics. As our database of protein structures grows, it is easier to predict protein structures based on similarities in amino acid sequence. Also, we have discovered by analysis of sequence databases that there are certain conserved protein families with high sequence homology in part, if not all of the amino acid sequence. Computer programs can currently predict protein structures by homology modeling when the sequence homology is as low as 25%. Therefore, if the amino acid sequences agree by more than 25%, the computer program can accurately predict the secondary and tertiary structure of the amino acid sequences.

GenBank. In this exercise, you will use a computer to access GenBank, the database repository of all DNA and protein sequences housed at the National Center for Biotechnology Information (NCBI) at the National Institutes of Health (NIH). Researchers and scientists submitted the sequences in GenBank. A unique number identifies each submission and is called an **accession number**. In the first part of this bioinformatics activity, you will use the accession number of the **pGLO** plasmid (also known as pBAD-GFPuv) to search GenBank and retrieve the plasmid sequence for pGLO.

Bioinformatics and Gene Structure. Once a sequence is determined, scientists can use computer programs to search the DNA for genes. Genes have a **start codon** (ATG) that signals the start of the protein and a **stop codon** which signals the end of the protein and for the ribosome to release the mRNA. These two particular codons must be ‘in the frame.’ This means that each codon (every three nucleotides) must code for a specific amino acid, from the start codon until the stop codon is reached. This is referred to as an **open reading frame** or ORF. Since not all ORFs encode a protein, you will use a common research tool called **BLAST**. This putative protein sequence for each ORF will be queried against the protein database in GenBank using BLAST. You will determine which sequence most closely matches the ORF.

Once the gene and ORF have been identified, an essential piece of information to clone or manipulate the gene using recombinant DNA technology is restriction mapping. In this lab, you will use the free software **NEBcutter** to determine the restriction map of the pGLO plasmid. The Protein Data Bank (**PDB**) is a valuable research database which hosts 3-D protein structures identified by a PDB-ID number. This software program also hosts **Jmol**, which you will use to view the 3-D structure of the GFP protein.

Part I: Using GenBank to search for and download the pGLO plasmid sequence

1. Open MS Word and save the open document with the lab name. Create a title at the top of the document. Remember to save your work as you go!
2. Open a web browser and go to the National Center for Biotechnology Information (NCBI) website. (www.ncbi.nlm.nih.gov/)
3. Click the **analyze** button and then nucleotide on the pull down menu.
4. The accession number for pGLO plasmid is **U62637**. Type this accession number in the search window at the top (with the nucleotide dropdown). Press Search.
5. Find the FEATURES section of the accession (left-hand side). How many coding sequences (genes) are reported for this plasmid sequence? _____
6. Find the ORIGINS section of the accession. Report the first 15 bases of the plasmid sequence.

Part II: Using NEBcutter to determine the restriction site map and ORFs of pGLO Plasmid

1. NEB has an excellent array of recombinant DNA cloning online tools. Go to <http://nebcloner.neb.com/> and explore. Click on the traditional cloning workflow, and Watch the video on molecular cloning. Draw a flowchart of cloning workflow below:

2. In a new tab, go to NEB website: www.neb.com and click on NEBcutter under Tools & Resources (<http://nc2.neb.com/NEBcutter2/>)
3. Type in the Accession number for pGLO plasmid (**U62637**) into the “GenBank Number” box. Click on the **Circular** button at **The Sequence is** section. This is important! Leave the other settings at default and click **submit**.
4. Once you press submit the window will reload and show an image of the plasmid with restriction sites and with the ORFs identified. Right-click on the restriction map, select **copy and paste it into your MS Word document** with an appropriate heading.
5. Click on one of the ORFs, a, b, or c. This ORF will open in a new window with the amino acid sequence of this gene. Copy and paste this gene sequence into your MS Word document with an appropriate heading.
6. To identify this gene, copy the sequence starting at “>” and paste it into the NCBI Protein BLAST database here: <https://blast.ncbi.nlm.nih.gov/Blast.cgi> Click on **Protein Blast**, then copy sequence into “Enter Accession Number.”
7. Scroll down and click **BLAST** at the bottom of the page. Note that this may take quite a few minutes depending on how busy the server is. Be patient. The window will refresh multiple times before it is complete.
8. When results are loaded, copy the conserved domain image “**superfamily**” at the top of the page (move cursor over top image, right click, and copy image) and paste into your MS word document with an appropriate heading.
9. Scroll down to **Sequences producing significant alignments** (under the multicolored chart). The sequences at the top of the table are the results which are the closest match to the query sequence.
10. Click on the **accession number** for the TOP RESULT only (on the right). Copy and paste this amino acid sequence (ONLY!) into your result MS word document. DO NOT copy the entire document – this is hundreds of pages long!
11. Go back to NEBcutter and repeat with the remaining ORFs.
12. Fill in the following table:

ORF in Plasmid	Top Search Result (protein name & organism)

Part III: Using JSmol to view a 3-D model of GFP

1. Go to the **Protein Database** website <https://www.rcsb.org/>
2. The PDB ID# for GFP is **3i19**. Enter this in the search field and click on **Go**.
3. Right click on 3D image of GFP protein structure, press copy, and paste into your MS Word document with appropriate heading.
4. Click on **“3D View Structure”** and explore the various 3D protein features of GFP. Note: As a reminder, you should be doing this work in an updated browser, such as Chrome. NOT MS explorer.
5. Observe the protein from all sides and answer the following questions:

- a. GFP is composed of several structures that comprise a barrel – what is the name of this secondary structure? How many are needed to make the barrel?
- b. Do you see any alpha helices? How many are there?
- c. Describe the ‘active site’ for fluorescence.

Part IV: Using Clustal Omega to compare GFP sequences from different species of jellyfish

GFP is a protein found in many species of jellyfish. It was originally isolated in *Aequorea victoria*. **ClustalQ** is an alignment tool that can be used to compare sequences of similar proteins.

1. Return to the **NCBI** website (www.ncbi.nlm.nih.gov/)
2. Copy and paste the FASTA format of each of the following GFP sequences into your MS Word document and press save! Searching the accession numbers below, click on FASTA format, and copy and paste the entire file from “>” onward.
 - a. AAC53663.1: pGLO plasmid
 - b. AAA27722: GFP from *Aequorea victoria*
 - c. AAN41637: GFP from *Aequorea coerulescens*
 - d. AAK02062: GFP from *Aequorea macrodactyla*
3. Copy and paste each FASTA file one after another. Start each on a new line, but do NOT insert a paragraph spacing or line between sequences.

4. Go to the **Clustal** website: <http://www.ebi.ac.uk/Tools/msa/clustalo/>
5. Copy and paste all four FASTA files and paste them into the search window provided. Ensure **Protein** and **Clustal with numbers** is selected. Click **submit**. It may take a few minutes to complete the request.
6. Click on **Show Colors**. Copy and paste your aligned sequences into your MS word results document. The colors represent different chemical properties of the amino acids:
 - o Red: small hydrophobic
 - o Blue: Acidic
 - o Magenta: Basic
 - o Green: Hydroxy, sulfhydryl, amine, and G
 - o Gray: Unusual amino acids

The GFP plasmid was generated by genetically modifying the native form of GFP to obtain a brighter and more stable fluorescent molecule. Looking at the consensus sequence and the different colored amino acid identifies the differences between the native GFP from *Aequorea victoria* and the GFP that has been genetically modified in the pGLO plasmid. What amino acids were changed? How do you think those changes modified the structure?

7. Click on the **Phylogenic Tree**, to observe the relationship between these species. Copy and paste it into your MS Word document (you may need to do a print screen for this tree).
8. Using your phylogenic tree results, what two GFP proteins are most alike? What two are most different? Discuss what this means – evolutionarily speaking?

9. Save and print out a copy of your results from your MS Word file and include with your report.



Lab Unit 6-B: Transformation of *E. coli* with a Recombinant GFP

Adapted from Bio-Rad's pGLO Transformation Kit (Cat#166-0003EDU)

http://www.bio-rad.com/LifeScience/pdf/Bulletin_9563.pdf

OBJECTIVES

Your performance will be satisfactory when you can:

- ✓ use sterile technique
- ✓ make *E. coli* cells competent for transformation
- ✓ transform *E. coli* with plasmid DNA
- ✓ select for recombinant clones on antibiotic selection plates
- ✓ analyze and troubleshoot the results of a transformation experiment

INTRODUCTION

DNA recombination or molecular cloning consists of the insertion of DNA fragments from one type of cell or organism into replicating DNA of another type of cell. The cell is said to be **transformed**, and many copies of the inserted DNA can be made in the cell. If the inserted fragment is a functional gene coding for a specific protein, many copies of that gene and translated protein could be produced in the host cell if there is a promoter preceding the site of insertion. This process has become necessary for the large-scale production of proteins (*Bacillus thurengiensis* toxin, insulin, human growth hormone, Factor VIII, etc.) that are of value in agriculture, medicine, and other sciences.

While transformation is a relatively rare event under natural conditions, it is possible to manipulate conditions to make transformation frequencies higher in the laboratory. For example, plasmids can be used as vectors to carry fragments of DNA into bacterial cells. **Plasmids** are closed, circular DNA molecules that are capable of autonomous replication within a host cell. There are many naturally occurring plasmids, but the plasmids used in the biotechnology laboratory are those that have a high copy number in host cells. After the host cell has been transformed with a **high copy number** plasmid, the plasmid will multiply and be maintained at levels of hundreds to thousands of copies within each cell.

Plasmids have been genetically engineered to contain a cluster of restriction enzyme sites within a short region of the plasmid called a **multiple cloning site**, or MCS. This allows for insertion of DNA fragments produced from a restriction digest to be incorporated into plasmid DNA at its multiple cloning sites after digestion with the same restriction enzyme. After allowing the sticky ends of fragments of target DNA to anneal to the complementary sticky ends of the plasmid, the DNA insert is fixed in place with covalent bonds by DNA ligase, forming a recombinant DNA (rDNA) plasmid. In this experiment, the plasmid contains a gene coding for a green fluorescent protein (GFP) isolated from a jellyfish (*Aequorea victoria*).

The plasmids used to transform bacterial cells also have a **selectable marker** gene. This gene codes for a protein that allows the scientist to distinguish cells that have been successfully transformed by plasmid DNA from those that have not. The most common selectable markers



are antibiotic resistance genes, which allow for selection of transformants by growth on media containing the antibiotic. Non-transformed cells will die and transformed cells will survive under these conditions. In this lab exercise, a plasmid with an antibiotic resistance gene for ampicillin (amp) is used; the ampicillin resistance gene codes for an enzyme that destroys the ampicillin in the surrounding growth media.

Making bacterial cells competent for uptake of DNA

The efficiency of transformation can be improved by carefully managing conditions before and during the transformation. For example, the choice of bacteria and plasmid can affect the efficiency of transformation, because many plasmids have a narrow host cell range and will only transform bacterial cells of a single species. Transformation frequencies are considerably higher when using fresh bacterial cells taken from actively growing cultures. Also, bacterial cells can be made **competent** for DNA uptake by pretreatment with chloride salts of divalent cations such as calcium, followed by a cold-shock and a heat-shock step. The metal ions and temperature changes affect the structure and permeability of the cell wall and cell membranes such that DNA molecules can pass through. Cells that are allowed to recover in non-selective growth media and at their optimal growth temperature following transformation also have higher transformation frequencies. The recovery time allows the transformed cell to amplify the plasmids and to express the antibiotic resistance gene required for survival on the antibiotic-containing selection medium.

Green Fluorescent Protein (GFP)


The green fluorescent protein (GFP) is a fluorescent protein naturally occurring in the Pacific jellyfish *Aequoria victoria* that has been successfully cloned into some organisms from bacteria to mice. Although originally chosen for its novelty of causing the transgenic organisms to glow green, GFP has been successfully used as a marker for transformation. Recent studies have created gene fusion in which the GFP gene is fused to genes of target markers on either the N- or C-terminus of the protein that they encode. The GFP becomes a marker for the intracellular location of the target gene product, tracking its migration by fluorescence microscopy into the nucleus, mitochondria, secretory pathway, plasma membrane or cytoskeleton. GFP can also be used as a reporter of gene expression levels as well as a measure of protein-protein interactions. Therefore, GFP is a very useful tool for both geneticists and cell biologists.

To date, the GFP gene has been introduced and expressed in many biological systems including yeast, bacteria, fish, plant, fruit fly, and many types of mammalian cells including human cells. This was such an important scientific discovery that Martin Chalfie, Osamu Shimomura, and Roger Y Tsien were awarded the 2008 Nobel Prize in Chemistry in 2008 for their roles in the discovery and development of GFP.

The green fluorescent protein is a medium-sized protein of 238 amino acids and a molar mass of 27,000 Daltons. In spectrophotometry, it shows a major absorption peak at 395 nm and a minor absorption peak at 475 nm. The characterizing molar extinction coefficients are 30,000 and 7,000 M⁻¹cm⁻¹ respectively. **Fluorescence at 509 nm** is not energy requiring and depends on the

amino acids serine-65, tyrosine-77, and glycine-67. This trimer forms a fluorescent chromophore after translation by cyclization and oxidation reactions.

Once isolated, the GFP is stable across a wide range of temperatures and pH. It is very resistant to denaturation, requiring treatment with 6 M guanidine hydrochloride at 90°C or pH of <4.0 and >12.0. Furthermore, it can renature completely within minutes following many denaturing protocols, including sulfhydryl reagents such as 2-mercaptoethanol.



Overall Shape of GFP Monomer (Carson, M, 1987. *J. Mol. Graphics* 5:103-106.)

GFP consists of a dimer, each made of a barrel-shaped cylinder made primarily of β pleated sheets on the outside and α -helices on the inside, a structure that is unique among proteins. This structure produces a compact domain that surrounds and protects the fluorophore located at the center of each cylinder as shown in Fig. 2.1. The N-terminal region of the protein acts as a “cap” on the end of the protein, further protecting the core fluorophore. When this cap is disrupted, the fluorescence may be easily quenched. The dimers are probably held together with the hydrophilic interactions of the pleated sheets on the outside of the cylinders.

GFP PLASMID – pGLO

pGLO contains several genes that enable replication of the plasmid DNA and expression of the fluorescent trait (phenotype) in bacteria following transformation. Some of the essential genes include:

- **GFP** — the jellyfish gene that codes for the production of Green Fluorescent Protein.
- **amp^r** — A gene that encodes the enzyme beta-lactamase, which breaks down the antibiotic ampicillin. Bacteria containing the amp^r gene can be selected by placing ampicillin in the growth medium.
- **ori** — The origin of pGLO plasmid DNA replication.
- **araC** — A gene that encodes the regulatory protein that binds to the pBAD promoter. Only when arabinose binds to the araC protein is the production of GFP switched on.
- **pBAD Promoter** — A specific DNA sequence upstream of the GFP gene, which binds araC-arabinose and promotes RNA polymerase binding and transcription of the GFP gene.
- **Multiple Cloning Sites** — Regions of known restriction (*NdeI*, *HindIII*, *EcoRI*, etc.) sites that permit insertion or deletion of the gene of interest.



Plasmid Map of pGLO plasmid

In this laboratory exercise, you will transform *E. coli* bacteria with a plasmid that has the GFP gene inserted into it. In our next lab exercise, you will induce this plasmid to produce recombinant GFP and purify it using column chromatography.

To learn more about Bio-Rad's pGLO transformation kit lab exercise, watch the following video: <https://youtu.be/c40UudFILGw> Summarize the procedure in a flow chart below:

References:

pGLO Bacterial Transformation Kit (Cat#166-0003EDU) Manual. 4006097 Rev E



pGLO Transformation Protocol:

LABORATORY SAFETY

The *Escherichia coli* strain used in this experiment is not considered a pathogen, but *E. coli* bacteria colonize the intestinal tracts of animals. Although it is rarely associated with any illness in healthy individuals, it is good practice to follow simple safety guidelines in handling and disposal:

1. Gloves and goggles should be worn at all times.
2. Wipe down the lab bench with antibacterial cleaner before starting the lab and before leaving the laboratory.
3. All materials, including plates, pipettes, loops, and tubes that come in contact with bacteria should be autoclaved before disposal in the garbage.
4. Wash hands thoroughly with soap and water after removing gloves.

MATERIALS

Each group:

Micropipettes and plugged tips
Sterile microcentrifuge tubes
Microcentrifuge tube racks
Floating microcentrifuge racks
Gloves
Sterile inoculating loops
Ice bath
2 LB Petri plates
2 LB/amp Petri plates
1 LB/amp/ara plate
Biohazard bags and stands
LB broth
Cold 50 mM CaCl₂ transformation solution

Class shares:

Benchtop cleaner
37° C incubator
42° C water bath & floats
Hot gloves
Longwave UV light source
Fresh starter culture of *E. coli*
Petri Dishes
(Lyophilized *E. coli*)
pGLO plasmid solution (0.08ug/ul)
Petri dishes for preparing agar plates
Agar
Transformation solution
Ampicillin vial
Arabinose vial

PROCEDURE

Part I: Preparing Plates (day 1)

Note: Part of this part may have already been completed. *Ask your instructor.*

1. Label Petri plates with a sharpie as demonstrated by your instructor. There will be 16 plates LB, 16 plates LB/amp and eight plates LB/amp/ara.
 - a. One stripe means LB (only)
 - b. Two stripes means LB-amp
 - c. Three stripes means LB-amp-ara
2. Clean the area that will be used to pour plates.
3. Prepare nutrient agar (autoclave-free) – (note: the instructor may have already prepared this in an autoclave for you before the start of the class)



- a. The agar plates should be prepared at least three days before the student experiment is performed. They should be left out at room temperature for two days and then refrigerated until they are to be used. The two days on the bench top allows the agar to dry out (cure) sufficiently to readily take up the liquid transformation solution.
 - b. To prepare the agar, add 500 ml of sterile water to a 1 L or larger sterile Erlenmeyer flask or glass bottle.
 - c. Add the entire contents of the LB nutrient agar packet. Swirl the flask to mix the agar, and autoclave-sterilize 20min, 121°C. Allow the LB nutrient agar to cool in a clean 50°C water bath.
4. Prepare arabinose and ampicillin
 - a. **Note:** Arabinose requires at least 10 minutes to dissolve - be patient. With a new sterile pipet, add 3 ml of transformation solution directly to the vial to rehydrate the sugar. Mix the vial; a vortex helps.
 - b. With a new sterile pipet, add 3 ml of transformation solution directly to the vial to rehydrate the antibiotic.
 5. Pour LB nutrient agar plates (LB, LB/amp, LB/amp/ara)
 - a. First, pour LB nutrient agar into the 16 plates that are labeled LB.
 - b. Stack the empty plates 4 to 8 high and starting with the bottom plate lift the lid and the upper plate straight up and to the side with one hand and pour the LB nutrient agar with the other.
 - c. Fill the plate about one-third to one-half (~12 ml) with agar, replace the lid and continue up the stack. Pour 16 plates in this fashion and label them as LB.
 6. Now, add the hydrated ampicillin to the remaining LB nutrient agar. Swirl briefly to mix. Pour into the 16 plates that are labeled as LB/amp using the technique utilized above.
 7. Last, add the hydrated arabinose to the remaining LB nutrient agar containing ampicillin. Swirl briefly to mix and pour into the eight plates labeled as LB/amp/ara using the technique utilized above.
 8. Allow plates to set at ambient temperature for next class. Keep plastic sleeves for plates! You may leave at ambient temperature for two days. Store at 4°C.
 9. Next class leave out one LB (only) plate per group and put the rest back into the sleeves. Store upside down (agar side up) in the fridge.

Part II: Starter Culture (day 2)

1. Rehydrate lyophilized *E. coli* by directly adding 250ul of transformation solution directly into the vial.
2. Recap, swirl to mix and allow to stand for 5 min at ambient temperature.
3. Each group prepares a started culture plate by streaking plate for isolated colonies.
 - a. Label the top of an LB (only plate) with your group initials, the date, LB, and *E. coli* starter culture. Keep a lid on the plate at all times. Remember this plate has no antibiotics! Anything from the air that lands on the plate will grow in this nutrient-rich environment.
 - b. Shake culture to mix just before use.
 - c. Using a 10ul loop, insert loop into bacteria vial and remove a loop full of rehydrated *E. coli*. The loop will have a shine to it (like when you are blowing bubbles).
 - d. Dispense the 10ul loop of culture on a labeled LB only plate, spread as on single line on the plate (see figure below). Dispose of the loop.
 - e. Obtain a clean loop and repeat streak, dispose of the loop. Repeat two more times. See diagram below and the demonstration by your instructor.
 - f. Return lid, turn upside down and incubate 24-48hrs at 25-37°C. The colonies should be 1mm wide.
 - g. NOTE: You must use a fresh starter plate. Prepare 1-2days before the transformation section.




Figure 1: Streaking for isolated colonies. This is your starter plate.
Watch this video for a detailed explanation: <https://youtu.be/0heifCiMbfY>

Part III: Transformation of *E. coli* (day 3)

BIOHAZARD! Note, all tips, tubes and gloves must be disposed of in the biohazard trash. Clean bench top before and after working. Discard gloves in biohazard trash and wash hands before leaving the lab.

1. Label the top of one closed microcentrifuge tube "+pGLO" and another "-pGLO." Label both tubes with your initials along the side.
2. Add 250 μ L of transformation solution (CaCl_2) into both tubes. Place the tubes on ice.
3. Use a sterile loop to pick up a single colony of bacteria from the starter plate. Pick up the +pGLO tube and immerse the loop into the transformation solution at the bottom of the tube. Spin the loop between your index finger and thumb until the entire colony is dispersed in the transformation solution. Place the tube back on ice. Using a new sterile loop, repeat for the -pGLO tube.
4. Examine the pGLO plasmid DNA solution with the UV lamp. Does it glow?
5. Add 10 μ L of the pGLO plasmid into the cell suspension of the +pGLO tube. Close the tube and return it to the ice. **Do not add plasmid DNA to the -pGLO tube.** Why, not?
6. Incubate the tubes on ice for 10 minutes.
7. While the tubes are incubating on ice, label your four agar plates around the bottom of the plate with the agar (not the lid). See instructor for a demonstration. Ensure you are labeling the correct plate!
 - o Label one LB/amp plate: +pGLO
 - o Label the LB/amp/ara plate: +pGLO
 - o Label the other LB/amp plate: -pGLO
 - o Label the LB plate: -pGLO
8. Heat shock the sample by transferring both tubes into the 42°C water bath for exactly 50 seconds (use a timer!). Immediately after 50 seconds, place both tubes back on ice. For the best transformation results, the change from the ice (0°C) to 42°C and then back to the ice must be immediate!
9. Incubate tubes on ice for 2 minutes.
10. Remove the tubes from the ice and place on the bench top. Open each tube, add 250 μ L of LB broth to the tube, pipet up and down to mix, and reclose the tube. Repeat with the second tube, using a clean pipet tip.
11. Incubate the tubes for 10 minutes at room temperature.

12. Flick the closed tubes with your finger to mix. Using a clean pipet tip for each sample, transfer 100 μL of the appropriate reagents onto the appropriate plates as follows:
 - o LB/amp plate: +pGLO tube
 - o LB/amp/ara plate: +pGLO tube
 - o LB/amp plate: -pGLO tube
 - o LB plate: -pGLO tube
13. Use a new sterile loop for each plate. Spread the liquid evenly around the surface of the agar by quickly skating the flat surface of a new sterile loop back and forth across the entire plate surface. Return the lid and allow liquid to absorb into the plate on the bench top. This may take 5-10 minutes.
14. Once the liquid had fully absorbed into the surface of the plate, invert, and stack your plates together. You can tape them together and label with group initials.
15. Place the inverted stack in the 37°C incubator overnight. Why do you incubate plates with the agar side facing down? The plates should be removed after approximately 18-24 hours and refrigerated until the next class period.
16. If there are extra LB (only) plates, use them to practice streaking for isolated colonies. Use your *E. coli* starter culture LB plate. Put these into the 37°C incubator overnight with your other plates.

Part IV: Analysis of transformants (day 4)

1. Observe the plates with the UV lamp. Do any plates have glowing green colonies? Why?
2. Count the number of colonies on each plate. A convenient method to keep track of counted colonies is to mark the colony with a marker on the outside of the plate as you count it. Record results in a table in the attached worksheet.
3. **Keep plates from this experiment for your next two labs!** Seal plates closed with Parafilm, or place them in a zip lock baggy so they do not dry out. Label and store them at 4°C as directed by your instructor.

Calculating Transformation Efficiency:

In many experiments, it is important to transform as many cells as possible genetically. For example, in some types of gene therapy, cells are collected from the patient, transformed in the laboratory, and then put back into the patient. The more cells that are transformed to produce the desired protein, the more likely that the therapy will work. The transformation efficiency is calculated to help scientists determine how well the transformation is working.

In the next page is a transformation calculation worksheet, which will walk you through step-wise how to calculate transformation efficiency: the number of transformed cells per microgram of DNA used.

Transformation Efficiency Worksheet: *Include this transformation worksheet in your lab report!*

The transformation efficiency is calculated using the following formula:

$$\text{Transformation efficiency} = \frac{\text{Total number of cells growing on the agar plate}}{\text{Amount of DNA spread on the agar plate (in } \mu\text{g)}}$$

1. Total number of transformed colonies on the LB/amp/ara plate: _____
2. The total amount of pGLO plasmid DNA in the bacterial cells on the LB/amp/ara plate:
Check stock of pGLO plasmid. Its concentration is: _____ ug/ul (0.08ug/ul)

Volume of pGLO plasmid used _____ ul

Multiply the concentration of pGLO plasmid by amount transferred to tube: _____ ug pGLO

3. Total volume in the transformation solution:

CaCl₂: _____ ul
pGLO plasmid: _____ ul
LB broth: _____ ul
Total Volume: _____ ul

4. Volume spread on the plate: _____ ul
5. Calculate the fraction of the total transformation solution used:
$$\frac{\text{Volume spread on LB/amp/ara plate}}{\text{Volume of total transformation solution}} = \text{_____}$$
6. Calculate the mass of pGLO plasmid on the plate:
Fraction of total transformation x total pGLO DNA in solution
_____ x _____ ug = _____ ug

7. Calculate the number of colony forming units (CFUs) per ug of DNA transformed
Transformation Efficiency:

$$\frac{\text{Number of colonies transformed}}{\text{ug of DNA on the plate}} = \text{_____} = \text{_____} * \text{transformants/ug}$$

*NOTE: Report transformation efficiency in a whole number

Lab Unit 6B - Transformation of *E. coli* with pGLO Plasmid Assignment

1. What is a bacteria colony? Is a single colony only one bacteria cell?
2. Describe how you would properly streak for a single, isolated colony. Why is this important?
3. Exogenous DNA (such as a plasmid) does not passively enter *E. coli* cells that are not competent. What treatment do cells require to be competent?
4. Fill out the following table with observations from your experiment: Show your calculations.

Experimental Conditions	# of colonies/plate	Appearance of colonies under UV light
+pGLO on LB/amp		
+pGLO on LB/amp/ara		
-pGLO on LB/amp		
-pGLO on LB		

LAB UNIT 7: PLASMID DNA ISOLATION

Jack O'Grady, M.S., Linnea Fletcher, Ph.D., Angela Wheeler, M.S., Patricia Phelps, Ph.D.

Adapted from Bio-Rad's Aurum Plasmid Mini-prep Kit (Cat# 732-6460)

OBJECTIVES

Your performance will be satisfactory when you can:

- ✓ understand the structure and function of plasmid DNA
- ✓ demonstrate safe and proper lab skills of working with biohazard material
- ✓ demonstrate the safe and proper use of lab equipment including micropipettes, microcentrifuge, and electrophoresis equipment
- ✓ demonstrate knowledge of DNA isolation techniques from live organisms
- ✓ understand the use of spectrophotometry in determining the quality and quantity of a nucleic acid preparation
- ✓ demonstrate the proper use of a low-volume spectrophotometer

INTRODUCTION

The isolation of nucleic acids is a common practice in the molecular lab. Scientists commonly perform a combination of the steps listed below to isolate nucleic acid:

1. The disruption of the cell membrane, and cell wall when necessary, by mechanical, chemical and enzymatic treatment.
2. Enzyme degradation is used for selective isolation of DNA (by RNase treatment) or RNA (by DNase treatment).
3. The separation of nucleic acids from other cytoplasmic components by combinations of these steps:
 - a. Phenol extraction of proteins, followed by selective precipitation of nucleic acids under high salt and cold alcohol treatment.
 - b. Selective precipitation of nucleic acids under high salt or cold alcohol treatment.
 - c. Selective adsorption onto a chromatographic matrix in a centrifuge (a "spin column") followed by desorption by a special buffer system.

The original *alkaline lysis method for purifying plasmid DNA* from bacterial cultures requires organic reagents and time-consuming steps to obtain high-quality DNA. *The Aurum plasmid mini-prep kit* (used in this lab) has been optimized for the rapid purification of high-quality, high-yield plasmid DNA. This kit uses the silicon dioxide exoskeleton of diatoms as the DNA binding matrix. This matrix carries a partially positive charge and therefore binds negatively charged nucleic acids. The advantages of this porous substrate include ease of resuspension, high affinity for DNA, simple and efficient processing, elution in an aqueous buffer, and an inherently large surface-to-volume ratio. All of these properties contribute to the highest purity and yields of DNA. Plasmid DNA purified with the Aurum plasmid mini-prep kit can be used directly for fluorescent sequencing, cell transfection, electroporation, and enzymatic restriction and modification.

DETERMINING THE QUALITY & QUANTITY OF NUCLEIC ACID

Once the nucleic acid has been isolated, its' quality and concentration is evaluated. The amount of the isolated nucleic acid used in subsequent steps depends on its concentration. The stability of the nucleic acid and its performance in subsequent enzymatic steps is affected by its purity. Both the amount of nucleic acid isolated and its purity are affected by the type of tissue that it is isolated from, the amount of tissue used, and the isolation technique used. In this lab exercise, you will evaluate your plasmid quality and quantity using:

1. Spectrophotometric analysis of ultraviolet absorption at A_{260} will provide concentration, and a ratio of A_{260}/A_{280} (contaminating protein) and A_{260}/A_{230} (contaminating carbohydrates) will give you an idea of quality.
2. Qualitative analysis using agarose gel electrophoresis is one of the most common methods that provide valuable information about size, quality, other contaminants, and relative concentration.

The Spectrophotometer

The de facto method for quantitating nucleic acids that all other methods rely on is ultraviolet absorption. When other methods are used, a nucleic acid standard is prepared based on its absorbance at 260 nm, measured by a spectrophotometer. An advantage in the use of a spectrophotometer in nucleic acid quantitation lies in its high precision and the fact that sample is not destroyed by the assay and can be put to further analysis after quantitation. Spectral analysis is very fast, another major reason for its routine use in a molecular lab for quantitating nucleic acids.




Figure 7.1: Simple diagram of a spectrophotometer. Reprinted with permission, <https://commons.wikimedia.org/wiki/File:Spetrophotometer-en.svg>

A spectrophotometer makes use of the transmission of light through a solution to determine the concentration of a solute within the solution. This is accomplished by placing a lamp on one side of a sample and a photocell or detector on the other side. All molecules absorb radiant energy at one wavelength or another, depending on the chemical types of functional groups they are comprised. Those that absorb energy from within the visible spectrum are known as pigments. Proteins and nucleic acids absorb light in the ultraviolet range.

The design of the single beam spectrophotometer involves a light source, a prism or grating that separates light into different colors or wavelengths, a sample holder and a photocell. Connected

to each are the appropriate electrical or mechanical systems to control the illuminating intensity, the wavelength, and for conversion of energy received at the photocell into a voltage fluctuation. The voltage fluctuation is then displayed on a meter scale, is displayed digitally, or is recorded via connection to a computer for later investigation.

The concentration of colored solute in a solution is directly proportional to the intensity of its color, which in turn is proportional to the amount of absorbance of light at the wavelength that the color absorbs. The color, or absorbance, of a solution, is also proportional to the path length that the light passes. Spectrophotometers are useful for measuring concentrations of solutions because of the relation of the intensity of the color (absorbance) in a sample to the amount of solute within the sample.

This is often expressed as the **Beer-Lambert Law or Beer's Law:**

Beer-Lambert Law

$$A = \epsilon C l$$

Where, A is absorbance at a given wavelength of light,
 ϵ is the extinction coefficient (amount of color absorbance of the solute per mole),
C is the concentration of solute in the solution (doubling the concentration doubles the amount of light absorbed)
l is the path length (if you double the width of the cuvette, you double the absorbance)

Given the geometry of a spectrophotometer, what is measured at the photocell is the amount of light energy which arrives at the photocell. The voltage meter is reading the amount of light TRANSMITTED to the photocell. Light transmission is not a linear function but is rather an exponential function. That is why the solution was APPROXIMATELY half as intense when viewed in its diluted form. Most spectrophotometers have a built-in means of direct conversion of this reading to absorbance. *Absorbance is essentially the opposite of transmittance: what light is not absorbed is transmitted. The percent transmittance is related to absorbance mathematically as:*

$$A = 2 - \log (\% T)$$

Where A is absorbance at a given wavelength of light
% T is the percent transmittance or $\frac{\text{light transmitted through a sample} \times 100}{\text{light transmitted through a blank}}$

With the aid of spectroscopy, the quantitative analysis of nucleic acids and proteins has established itself as a routine method in many laboratories. Both nucleic acids and proteins absorb in the ultraviolet range, but while nucleic acids absorb strongly at 260 nm, proteins absorb more strongly at 280 nm. The absorption of 1 OD (Optical Density or Absorbance Unit, a dimensionless quantity) is equivalent to approximately 50 $\mu\text{g/ml}$ dsDNA, approximately 33 $\mu\text{g/ml}$ ssDNA, 40 $\mu\text{g/ml}$ RNA or approximately 30 $\mu\text{g/ml}$ for oligonucleotides.

Purity determination of DNA Interference by contaminants can be recognized by the calculation of “ratio.” The ratio A_{260}/A_{280} is used to estimate the purity of nucleic acid since proteins absorb at 280 nm. *Pure DNA should have a ratio of approximately 1.8, whereas pure RNA should give a value of approximately 2.1.* Absorption at 230 nm reflects biological contaminants of the sample such as carbohydrates, peptides, salts, or proteins. Also, many chemicals commonly used in nucleic acid preparations, such as phenol, EDTA, and SDS, can be detected by their absorbance at 230 nm. In the case of pure samples, the ratio A_{260}/A_{230} should be approximately 2.2. A ratio <2.2 indicates contaminants.

MONITORING BACTERIAL GROWTH USING OPTICAL DENSITY

Optical density is a common microbiology technique used to monitor the growth of microorganisms. Although a spectrophotometer is used in this measurement (600nm), it is not the amount of light absorbed, but rather, light scattered by the microorganisms. A typical microbial growth curve can be plotted by monitoring the OD600 by time. The OD600 indicates indirectly through light scattering an approximate number of live cells. See figure below.




Figure 7-2: A typical microbial growth curve (3)

During the **lag phase**, the cells are adjusting to their new environment and may not reproduce immediately. Depending on the species this lag phase can last a few hours (typically) up to a few days as it adjusts to its new growth conditions. In the **log phase** of growth, the cells reproduce logarithmically. This is an important stage for researchers as the metabolic rate of individual cells is at its maximum, the cell wall is intact, and growth is stable and healthy. As nutrients become depleted, bacterial growth slows and enters the **stationary phase**. During this phase nutrients are depleted, waste is built up, and rate of reproduction stops (the number of dying and produced cells is equal). It's important to note that during this phase the metabolic rate of live cells slows.

During the **death phase**, the number of dying cells outnumbers the live cells substantially. There may still be some live cells in the culture. For some organisms, during this phase, they may develop resting structures such as endospores.

AGAROSE GEL ELECTROPHORESIS TO ANALYZE NUCLEIC ACID

As discussed in lab unit 4, agarose gel electrophoresis is a molecular biology technique used to separate a mixture of nucleic acid fragments according to size in an agarose matrix. DNA or RNA samples are separated by applying an electric field to move the largely negatively charged molecules toward the positive electrode. In this lab, you will analyze your plasmid DNA quality using agarose gel electrophoresis.


Since DNA molecules are double-stranded (and the same shape), they migrate through the gel matrix separating according to size alone. Smaller molecules travel faster through the gel matrix than larger molecules at a rate that is inversely proportional to the log₁₀ of the number of base pairs. Therefore, by running a known molecular weight ladder at the same time as your nucleic acid sample, you can determine the size of the nucleic acid fragments using a linear regression analysis. In addition to agarose concentration and molecular weight, the mobility of nucleic acids in agarose gels is also influenced by the molecular conformation of the nucleic acid.

The three forms of circular plasmid DNA are:

Form I: Closed, circular, negatively supercoiled DNA

Form II: “Nicked” DNA, which has been partially cut through one strand of the DNA, causing it to unwind

Form III: Linear DNA, which has been cut by restriction enzymes



Form I usually has the greatest electrophoretic mobility of all DNA forms because supercoiled DNA molecules tend to be compact. Think of it as a bullet moving through the agarose gel. The size of closed circular plasmid DNA cannot be determined on an agarose gel relative to the bands of a linear marker. To determine the size of a plasmid, you must cut the plasmid with restriction enzymes. The slowest moving DNA of all is Form II because it is an open circle with a partially unwound strand, which causes it to drag in the gel. While the linear DNA can snake through the agarose particles, reducing drag, the circular, unwound plasmid cannot. Form III, or plasmid DNA that has been cut with restriction enzymes, has decreased mobility because the linear DNA is like a long string or rod that can drag. Unlike the other forms of DNA, *linear DNA migrates through a gel at a rate that is inversely proportional to the logarithm of its molecular weight*. Therefore, the molecular weight of linear DNA can be estimated from a gel if compared to DNA fragments of known molecular weight (**markers**).

Rough handling of the plasmid DNA or endogenous endonucleases can result in nicking of the plasmid and **linearization**. Long-term storage of plasmids will result in increased amounts of Forms II and III due to the presence of endogenous endonucleases (see below). Uncut plasmid should be stored in the refrigerator to avoid the risk of ice crystals shearing, or damaging, the supercoiled DNA molecule. A cut plasmid can be kept in the freezer with no risk of damage.

References:

1. Bio-Rad. “Aurum Plasmid Mini-prep Kit.” (Cat# 732-6460)
2. Robert W Bauman “Microbiology,” 4th edition. 2015.
3. Prince George Community College. < <http://academic.pgcc.edu/>>. April 2014.

Lab Unit 7-A: Plasmid Mini-preparation

SAFETY: What are safety considerations when working with *E. coli* and plasmids? Think about both handling these materials AND waste disposal.

PROTOCOL

Day 1: Overnight cultures

This may have been done for you. Ask your instructor.

1. Obtain your plate of *E. coli* pGLO transformants on LB/amp/ara and a labeled tube containing 5 mL Luria broth with ampicillin (LB/amp).
2. Use a sterile inoculating loop to pick up a single transformant colony glowing green.
3. Place the loop into the LB and swirl to completely transfer the bacteria into the broth.
4. Place in a 37°C, vigorously shaking incubator 18-24hrs. It's important that this culture be fresh – alternative procedures include 37°C /48hrs or 25°C/48hrs (highest yield). It is not recommended to store cultures at 4°C before extracting this may result in a low yield.

Day 2: Plasmid Isolation & NanoDrop Analysis of DNA

MATERIALS: For parts, 1 & 2 below, create a comprehensive materials list. Compare with your lab partner before you begin!

--	--

PART I: Determining the Optical Density of the bacterial culture

Twelve OD•mL is typically recommended to obtain a reasonable yield of plasmid DNA to perform a restriction digest. If you have too few bacteria, your plasmid yield will be low, if you have too many bacteria, your lysis will be inefficient, and that too will result in a low yield of plasmid. To calculate the OD•mL, measure the optical density of the bacterial culture in a spectrophotometer, set to 600nm (OD600), and multiply the culture volume in mL. Why measure at 600nm?!

Calculate the volume of culture needed to prepare 12 OD•mL of liquid culture:

1. Warm up spectrophotometer and collect materials.
2. It is recommended to dilute your culture 10-fold in warm LB broth. You will need 1mL final volume for the cuvettes.
3. Set the spectrophotometer to 600nm and blank with LB broth.
4. Measure the diluted sample: _____ OD600
5. Calculate:

$$12 \text{ OD}\cdot\text{mL} = \text{_____} (\text{OD}) \times \text{dilution factor } 10 \times (\text{X mL culture})$$

Solve for X: _____

PART II: Plasmid Isolation

PROTOCOL	EXPERIMENTAL NOTES
<ol style="list-style-type: none"> 1. Pellet 12 OD•mL of liquid culture in a labeled microcentrifuge tube or recommended) a total volume of 4mL of the 5mL culture. <ol style="list-style-type: none"> a. Centrifuge 4mL of culture into one 2mL microcentrifuge tube by centrifuging 2mL of culture at 12,000rpm for one minute, b. Remove the supernatant by decanting into a liquid biohazard container c. Add into the same tube another 2mL of culture and repeat centrifugation step. d. Pour off all of the supernatants. 2. Add 250 µL of the Resuspension Solution and vortex until the cell pellet is completely suspended. 3. Add 250 µL of the Lysis Solution and mix by GENTLY inverting the capped tube 6-8 times. The solution should become viscous and slightly clear if the cell lysis has occurred. 4. Add 350 µL of the Neutralization Solution and mix by GENTLY inverting the capped tube 6-8 times. A visible precipitate (consisting of cellular debris) should form. 	<p>OD₆₀₀: _____</p> <p>Calculate OD•mL:</p>

5. Pellet the cell debris for five minutes at 12,000 rpm in a microcentrifuge. A compact white pellet will form along the side or bottom of the tube. The clear supernatant in this step contains the plasmid DNA.
6. Insert a plasmid mini-column into one of the 2 ml cap-less wash tubes supplied.
7. Transfer the supernatant from step 5 to the column. Centrifuge at 12,000 rpm for one minute. The purpose of this step is to bind the plasmid DNA to the column.
8. Remove the spin column from the wash tube, discard the filtrate at the bottom of the wash tube, and replace the column in the same wash tube.
9. Add 750 μ L of Wash Buffer and centrifuge at 12,000 rpm for one minute. The wash buffer contains ethanol and washes away impurities from your sample.
10. Remove the spin column from the wash tube, discard the filtrate at the bottom of the wash tube and replace the column in the same wash tube. Centrifuge for an additional minute. It is important to spin twice to remove residual traces of ethanol.
11. Remove the spin column and discard the wash tube. ***Place the column in a clean elution tube.***
12. Add 50 μ L of elution solution directly on top of the white membrane in the middle of the column. Let sit at room temperature for one minute to saturate the membranes on the column.
13. Elute the plasmid DNA from the membrane by centrifuging at 12,000 rpm for one minute.
14. Transfer your plasmid DNA preparation to a clean, labeled 1.5mL microcentrifuge tube.
15. If you are not analyzing your plasmid preparation today, store your plasmid at 4°C in a storage box provided.

Lab Unit 7-B: Spectrophotometric Analysis of Plasmid DNA

Determine the concentration of your eluted DNA on the low-volume (NanoDrop) Spectrophotometer. Include this worksheet in your lab report!

1. Using the SOP provided in your SOP booklet set the NanoDrop to DNA (nucleic acid).
2. Blank the NanoDrop with 1.5ul of elution solution.
3. Measure 1.5ul of plasmid DNA solution. If possible, print out the spectral graph from the NanoDrop and include with your report.
4. Record plasmid DNA concentration as well as the A260/A280 ratio:

Concentration: _____ ng/ μ L


A260/A280: _____ A260/A230: _____

5. What does the A260/A280 ratio tell you? What is the optimal ratio for plasmid DNA? How does your ratio compare?

6. Obtain a storage box for the whole class to use. Label the box with the class name and instructor name. Ensure your tubes are labeled, so they are not mixed up with your classmate's tubes. Store the box at -20°C until next lab exercise.

****OPTIONAL** ANALYZE YOUR PLASMID ON AN AGAROSE GEL- Use protocol in Lab Unit 4.**

Frequently researchers will analyze their plasmid DNA isolations on an agarose gel to determine the quality of the plasmid preparation. Fully intact supercoiled DNA migrates quickly in one band on an agarose gel – but not according to size – why? When there are small nicks in the plasmid, you may get several different populations that migrate at differently.





4. In this lab, you added Cell Lysis solution to your bacterial pellet. What is the purpose of this solution and what do you think would happen to the results of your experiment if you left out this step?

5. What is the purpose of the spin filter? How does DNA bind to the filter? How do you get DNA to elute from the filter?

CONCLUSION:

LAB UNIT 8: RESTRICTION ENZYME DIGEST

Jack O’Grady, M.S., Linnea Fletcher, Ph.D., Angela Wheeler, M.S., Patricia Phelps, Ph.D.

OBJECTIVES

Your performance will be satisfactory when you can:

- ✓ understand how restriction enzymes are used in a biotechnology lab
- ✓ demonstrate safe and proper lab skills of working with recombinant DNA material
- ✓ demonstrate the safe and proper use of lab equipment including micropipettes, microcentrifuge, and electrophoresis equipment
- ✓ demonstrate knowledge of restriction enzymes to cut double stranded DNA
- ✓ prepare agarose gel, and prepare and load DNA samples onto an agarose gel
- ✓ analyze electrophoretic results


INTRODUCTION

Restriction enzyme digests are a tool used almost daily by molecular biologists. Recall that restriction digests can be used to cut DNA at specific nucleotide sequences, allowing for subsequent manipulation of the DNA, and often to generate recombinant DNA molecules. Restriction digests are also important as a diagnostic tool. Creating recombinant DNA molecules requires numerous steps involving different DNA molecules.

In this lab, you are continuing to work with the plasmid pGLO DNA you just isolated. Recall that in the previous lab you transformed bacterial cells with the pGLO plasmid. Bacterial colonies from this transformation were grown up in LB nutrient broth containing ampicillin as a selective agent. In this lab, you isolated and purified plasmid DNA from these bacterial cells. At this point, it is wise to confirm that the plasmid DNA isolated from your bacterial “overnight” is indeed the pGLO plasmid, and not chromosomal DNA or other foreign DNA. To do this, we will first presume that we have in fact isolated the correct DNA from the transformed bacteria.

Verifying Plasmid Insert


Most of the time researchers know the sequence and map of their plasmid DNA. Using this map, you can select the appropriate restriction enzyme to determine if you, in fact, have the plasmid with the insert. Restriction analysis can even be used to determine the orientation of the insert in the plasmid. In the example below, digestion with enzyme RE1 will linearize the 7200bp plasmid into one single 7200bp fragment. Alternatively, digestion with both RE1 and RE2 will result in two bands; the 1200bp insert and the 5000bp backbone.



pGLO Plasmid

The pGLO contains several genes that enable replication of the plasmid DNA and expression of the fluorescent trait (phenotype) in bacteria following transformation. Some of the essential genes and other elements include:

- **GFP:** The jellyfish gene that codes for the production of Green Fluorescent Protein
- **amp^r:** A gene that encodes the enzyme beta-lactamase, which breaks down the antibiotic ampicillin.
- **Ori:** The origin of pGLO plasmid DNA replication
- **araC:** A gene that encodes the regulatory protein that binds to the pBAD promoter. Only when arabinose binds to the araC protein is the production of GFP switched on
- **pBAD Promoter:** A specific DNA sequence upstream of the GFP gene, which binds araC-arabinose and promotes RNA polymerase binding and transcription of the GFP gene
- **Multiple Cloning Sites:** Regions of known restriction (*NdeI*, *HindIII*, *EcoRI*, etc.) sites that permit insertion or deletion of the gene of interest.



Restriction Map of pGLO plasmid DNA

How to choose percentage of Agarose Gel

When performing restriction analysis, you will obtain varying sizes of DNA. Depending on the size of DNA fragments you will choose a concentration of agarose gel that will help resolve those fragments. The following table is a recommended agarose percentage for size separation.

Other factors affecting the resolution of bands include voltage, using a wider comb, and loading less DNA. Note, you would also choose an appropriate marker (DNA ladder) to correspond to the

size fragments you will obtain. In our lab, we use a broad ladder and only consider the linear portion of that ladder to determine the size of the DNA fragment.

% Agarose	Size Fragments (bp)	Recommended Voltage
0.8	2,000 to 30,000 (Genomic DNA and large plasmids)	80-125V or lower for genomic DNA
1.0	500 to 10,000 (most general applications, one fragment in lane)	100-125V or lower
1.5	200 to 3,000 (better resolution for smaller fragments)	100V or lower
2	50 to 1,000 (better resolution for small fragments)	80V or lower
2.5	50 to 800 *better resolution of small fragments from each other in the same lane	80V or lower

Lab Unit 8-A: Bioinformatics: Restriction Enzyme Digest of Plasmid DNA

PROCEDURE

Part I: NEB Cutter

As part of your pre-lab exercise, complete the following short bioinformatics exercise.

1. Open up NEB cutter: <http://tools.neb.com/NEBcutter2/>
2. Insert the Accession number for this plasmid (Accession #: U62637).
3. Click on "All commercially available" and "circular" DNA. Press submit.
4. The default is single cutters, but you can toggle between single cutters and multiple cutters by pressing the button at the bottom of the screen "Display."
5. Copy and paste the single cutter map into an MS word document. Submit this with your report.
6. Based on this map, predict the size of fragments if you cut this plasmid with the following enzymes:
 - EcoRI: _____
 - Double digest: EcoRI & EcoRV: _____
7. From this map, can you see one or two restriction enzymes that will cut out GFP? What are they?
8. Confer with your classmates. As a class, decide the best enzyme(s) to use to determine if you have GFP in your plasmid: _____
9. Ask your instructor if you have these enzymes available in the lab.
10. Predict the size of restriction fragments you expect to see on the agarose gel after restriction digestion with these enzymes. You may use NEB cutter to help you out!

11. What percentage of agarose gel do you recommend to resolve these fragments from each other? See introduction!

1. It is ideal to digest a maximum of 1 μg of DNA in these reactions to see DNA fragments on an agarose gel. Using the concentration of the isolated plasmid DNA, calculate the volume needed to have 1 μg in each tube. The final volume in each tube after every reagent is added should be 20 μL .

NOTE: If your concentration is below 100ng/ul, you will need to use the maximum volume of your plasmid.

Calculate actual mass used in each reaction:

Concentration DNA: _____ ng/ul x _____ ul = _____ ng

2. Using the concentration of the buffer (10X), calculate how much should be added so that when the final volume is 20 μL , the buffer concentration is 1X.

Calculations:

3. Calculate the amount of water needed in each tube, after accounting for the volumes of the other reagents, to make a final volume of 20 μL . Record in the table.
4. Add the reagents to the tubes in the order listed below. **Keep DNA and enzymes on ice. Use a new tip for each solution.** Tube 1 is the control and will contain only DNA, buffer, and water, *NO* enzyme! Why not?
 - a. Water should be added first to each tube (if needed).
 - b. The buffer, which should be added second, maintains the proper pH for enzymatic activity. This must be thawed and vortexed to resuspend buffer components.
 - c. Then DNA should be added. Make sure you use a clean tip for each DNA sample.
 - d. The enzyme should always be added LAST. Pipette the enzyme directly into the solution at the bottom of the microcentrifuge tube, and pipet up and down to mix. The enzyme should be retrieved directly from the freezer and returned immediately to the freezer. It is in glycerol and does not need to be thawed. Pop-spin before opening, KEEP ON ICE, DO NOT VORTEX. Return to freezer immediately.
5. Gently flick to mix, pop spin all the tubes in a picofuge to ensure all components are in the bottom of the tube.
6. Place all tubes, with lids closed, (in a floating tube rack) in a 37° C water bath to incubate for at least 60 minutes and up to 3 hours. However, for this procedure, *it may be beneficial to incubate overnight.*

Lab Unit 8-C: Agarose Gel Electrophoresis of Restriction Digest DNA Fragments

Part I: Prepare Agarose Gel with SYBRsafe

Materials required:

Agarose	Balance, weigh boats
1X gel running buffer	SYBR Safe stock solution, 10,000X
Gel-casting apparatus	Ultraviolet (UV) light box (or imaging system)
Gel electrophoresis box	Erlenmeyer flask
Power supply	Saran Wrap
Microwave oven or hot plate	Micropipettes and tips
Molecular weight DNA ladder (1kB DNA ladder, Promega G7541)	


1. **Set up a gel casting apparatus** as directed by your instructor. Ensure the bumpers are snug in place to avoid leaking during casting. Insert a comb into the tray, choosing a comb by the size and number of wells it will create. Remember to include your molecular weight marker when deciding on the comb size. Set the apparatus on a flat surface that will be undisturbed while the gel is solidifying.
2. **Prepare 400 mL of 1X TAE buffer.**

Calculation: How will you prepare 400mL of 1XTAE from the 50XTAE stock you prepared earlier in the semester?

3. **Prepare 30mL of agarose in 1X TAE buffer.** What percentage of gel did you decide to prepare based on the size of your restriction enzyme DNA fragments? _____ %
4. Weigh out the required amount of agarose and add it to the appropriate amount of 1X gel running buffer in an Erlenmeyer flask. The flask should be at least twice the volume of the buffer and no more than 5X the volume. For example, to prepare 30 mL of a 1% agarose gel, add 0.3 g of agarose to 30 mL of the 1X buffer in a 125-mL flask.


Calculation: How much agarose will you need to make 30mL of a _____% agarose gel?

5. Place the flask on a level surface and carefully mark the glass at the fluid level with a permanent marker (DO NOT mark the white marking area with a permanent marker! It cannot be removed).
6. Heat the mixture until all agarose has dissolved. A hot plate or microwave oven can be used; using a hot plate will take more time. Interrupt the heating at regular intervals and swirl the container to mix the contents. The solution should be brought to a boil, but do not allow the solution to boil over. Microwave for a minimum amount of time to avoid buffer evaporation, which will cause a dramatic increase in the percentage of the gel. For example, a 30mL solution will require less than 1min total in a typical microwave.



CAUTION! The flask will be hot! Use hot hands to remove flask and swirl. Do not hold over the top of your face while swirling.

7. When the agarose is completely dissolved, the solution will be completely clear and homogeneous; that is, you will not observe any granules or threads of non-dissolved agarose in the solution. Observe the fluid level about the mark you made on the flask. If a significant amount of water has evaporated, carefully add water to return to the level of the mark and swirl the solution.
8. Cool the solution to 50 – 60°C. This is typically cool to the touch. Do not cool too long or it will solidify in the flask.



NOTE: If you make a mistake with this solution DO NOT dump out down the sink. The agarose will solidify in the sink and clog it. Ask your instructor where to dispose of the solution. Usually, the trash can is fine.

9. **Add SYBR Safe** stock solution (10,000X) to a final concentration of 1X. This is a 10,000-fold dilution, so for each 100 mL of agarose solution, add 10 μ L of stock SYBR Safe solution.

Calculation: How much SYBRsafe will you add to your flask to achieve a final concentration of 1X starting with a 10,000X stock?

10. Pour the gel immediately into your casting tray. Allow the gel to form completely; typically, 20 minutes at room temperature is sufficient. Remove the comb and bumpers from the gel, place the gel in the electrophoresis chamber, and add a sufficient volume 1X TAE gel running buffer to just cover the surface of the gel.

⇒ NOTE: Do NOT pour molten agar down the sink. Leftover agar can be kept in the fridge or allowed to solidify in the flask and disposed of in the trash.

11. Each sample and marker should contain loading dye to a final concentration of 1X. Gel loading solution is usually supplied at a 2-10X concentration and contains several substances. Glycerol makes loaded samples heavier than the buffer, allowing them to sink more quickly into the wells and prevents them floating away. The solution is buffered to prevent degradation of the nucleic acids it is mixed with. Finally, a visible tracking dye (usually blue, green, or orange) is present to allow the technician to monitor the progress of electrophoresis without removing the gels and observing with UV light.

Add 5ul of 5-6X load dye to each of the restriction digests. Flick to mix and pop spin in a picofuge.

12. The molecular weight marker should already be prepared for you. Collect an aliquot to run on the gel with your samples.

13. Load the total volume in each of the molecular weight marker (6ul) and samples (25ul) into the wells. By convention, the first left lane is the molecular weight marker. Do not skip lanes when loading samples. Record the order of samples in your electrophoresis documentation form.

14. Place the cover on the electrophoresis chamber in the correct orientation and connect leads to a power supply. Set the power supply to approximately 80 volts and allow to electrophorese until the tracking dye is approximately 3/4th the way to the bottom of the gels (approximately 1 hour for 2% minigels).

15. Fill out a gel electrophoresis documentation form

16. Visualize DNA bands by placing gels on a UV transilluminator.

17. Capture two sets of images; low exposure and high exposure. The low exposure is to ensure you see all the molecular weight marker bands, the cut, and uncut plasmid. The high exposure is to see the insert fragments which may be lighter.

- a. Print out both exposures and affix both to the same gel electrophoresis form. Turn this in with your lab report.
- b. Save a PDF copy to your removable drive. Print out a large image (half the size of a page). This will be used to analyze your data.

18. Save the remaining plasmid DNA in a labeled microcentrifuge tube. Plasmid DNA should be stored at 4°C.

PART II: Data Analysis

Include this worksheet with your lab report!

1. Using the large printed out copy of your images, measure the distance migrated by each band in your molecular weight marker (in mm!). Record this data in a table below. It is recommended you look up an image of the marker used online (Promega G7541). The 1000, 3000, 8000, and 10,000 bp bands are brighter intensity making them easier to identify.

Molecular Weight Marker (Kilobases)	Log MW	Migration Distance (mm)
250		
500		
750		
1000		
1500		
2000		
2500		
3000		
4000		
5000		
6000		
8000		
10,000		

2. Using MS Excel, graph your molecular weight standard curve. Plot log molecular weight on the y-axis and migration distance on the x-axis. Label you axis and give the graph and appropriate descriptive title (NOT logMW vs. migration distance!).
3. Using MS Excel, determine the equation of the line for the linear points of the curve. The highest MW marker may not be in the linear range, so you can exclude it if it falls off the linear portion of your graph.
4. Equation of the Line: _____
5. Determine the R^2 value for your graph. This is a correlation coefficient that will tell you how well your data correlates linearly to each other. An R^2 value >0.95 is acceptable, but >0.98 is preferred.
6. R^2 Value: _____
7. Measure the distance of each of the DNA bands in your lanes. Use the equation of the line to determine the molecular weight of each of your samples. Record in the table below.

Experimental Samples	Migration Distance (mm) of each band	Log MW of each band	MW (kilobases) of each band
Tube 1: Undigested control plasmid			
Tube 2: Restriction digested sample pGLO plasmid			

8. Turn a copy of your graph with your completed lab report.
9. Did your plasmid migrate at the expected size? Why or why not? See your NEB cutter predictions.



Lab Unit 8 – Plasmid Isolation Assignment

1. Include your NEB restriction digest map with your report.
 2. Based on this map, predict the size of fragments if you cut this plasmid with the following enzymes:
 - o EcoRI: _____
 - o Double digest: EcoRI & EcoRV: _____
 3. Which enzyme did the class choose to cut out GFP from this plasmid? Why?

 4. Predict the size of restriction fragments you expect to see on the agarose gel after restriction digestion with these enzymes. You may use NEB cutter to help you out!
-
5. What percentage of agarose gel did the class choose to run the digested fragments? Why?

6. Can you accurately determine the size of an uncut plasmid on an agarose gel? Why or why not? (HINT: READ INTRO)

7. Analyze your data. In your report include the gel documentation form and a labeled picture of your gel.
- o Label on your gel the molecular weight bands
 - o Label on your gel the vector and insert of the restriction digested plasmid.
 - o Determine the molecular weight of your vector and insert and label this on your gel.

MW Size of Insert: _____

MW Size of Vector: _____

8. Did you obtain the predicted size of your insert and vector? Why or why not?

CONCLUSION:

LAB UNIT 9: PROTEIN PURIFICATION USING COLUMN CHROMATOGRAPHY

Jack O'Grady, M.S., Linnea Fletcher, Ph.D., Angela Wheeler, M.S., Patricia Phelps, Ph.D.,
Adapted from Bio-Rad's GFP Chromatography kit

OBJECTIVES

Your performance will be satisfactory when you can

- ✓ demonstrate the safe handling of the biohazardous material
- ✓ extract proteins from cells
- ✓ separate a complex protein mixture using hydrophobic interaction chromatography
- ✓ discuss the role of green fluorescent protein

INTRODUCTION

Chromatography is a very powerful method for separating complex mixtures of biomolecules into separate components. There are many types of chromatography, but in each case, the separation of components of a mixture is based on differences in chemical and physical properties of the components.

In all types of chromatography, the separation takes place between two different phases: the **stationary phase** that does not move and the **mobile phase** that moves steadily past the stationary phase. Different components of a solution will separate due to their differential affinity for the stationary, compared with the mobile, phase. The stationary phase can be a flat sheet (as in paper or thin layer chromatography) or a column of material (as in liquid or gas chromatography). The separations can be based on molecule size (as in size exclusion or gel permeation chromatography), by charge and polarity (as in ion exchange chromatography), or by specific binding (as in affinity chromatography).

This lab exercise involves hydrophobic interaction chromatography, in which components bind, or **adsorb**, to the stationary phase due to hydrophobicity. The stationary phase is made of insoluble particles, also called a resin, of polysaccharide beads called Sepharose™. These beads are small (40 to 165 μm diameter), and made of agarose that has been chemically cross-linked to make the beads less likely to be crushed in a large column. While Sepharose is hydrophilic, cross-linking phenyl groups make it into the hydrophobic phenyl sepharose. Proteins with patches of hydrophobic (literally, “water-fearing”) amino acids on their surfaces will be attracted to the phenyl groups on the resin. Higher salt concentrations and higher temperatures strengthen these hydrophobic attractions. Under high salt conditions, even the least hydrophobic proteins will bind to the phenyl sepharose beads, but at low salt conditions, only the most highly hydrophobic proteins will remain bound to the phenyl sepharose beads. You may use an equilibration buffer with high salt concentration to bind the proteins in a mixture to phenyl sepharose and use elution buffers with successively lower concentrations of salt to separate them from each other according to their relative hydrophobicity.

Your goal in this lab is to use hydrophobic interaction chromatography to purify GFP from a bacterial cell lysate. Proteins are long chains of amino acids, some of which are very hydrophobic. GFP has many patches of hydrophobic amino acids, which collectively make the entire protein hydrophobic. Moreover, GFP is much more hydrophobic than most other bacterial proteins. We can take advantage of the hydrophobic properties of GFP to purify it from the other, less hydrophobic (more hydrophilic or "water-loving") proteins.

First, you will obtain a liquid bacterial culture that was prepared by isolating a single green transformant *E. coli* colony, adding it to a tube of the liquid medium, and then incubating with vigorous shaking overnight. You will process this culture to lyse the cells and release the proteins contained therein. You will load the cell lysate onto a HIC column in a high salt buffer. The salt causes the three-dimensional structure of proteins to change so that the hydrophobic regions of the protein move to the exterior of the protein and the hydrophilic ("water-loving") regions move to the interior of the protein. This will cause most proteins to adsorb to the column. As the salt concentration of the buffer is decreased, the three-dimensional structure of proteins changes again so that the hydrophobic regions of the proteins move back into the interior and the hydrophilic ("water-loving") regions move to the exterior.


These four buffers comprise the separation scheme:

Equilibration Buffer—A high salt buffer (2 M $(\text{NH}_4)_2\text{SO}_4$)

Binding Buffer—A very high salt buffer (4 M $(\text{NH}_4)_2\text{SO}_4$)

Wash Buffer—A medium salt buffer (1.3 M $(\text{NH}_4)_2\text{SO}_4$)

Elution Buffer—A very low salt buffer (10 mM Tris/EDTA)



*****ASSIGNMENT***** Before coming to class, read the entire detailed protocol and complete the missing pieces to the short protocol you will use in class to perform the experiment.

MATERIALS

Thoroughly evaluate this procedure and come up with a detailed list of materials (reagents & equipment) you will need to complete this exercise. Be specific! If you need micropipette tips, what size? If you need microcentrifuge tubes, how many? What equipment do you need?

Materials & Reagents needed per Person	Equipment & Materials to share with class

SAFETY PRECAUTIONS

UV radiation is harmful to your eyes. Be very careful when handling the hand-held UV light and always point down or away from you or your lab partner (not up towards your face and eyes).

Always use eye protection when using UV light.

What other safety precautions do you need to consider in this lab?

PROCEDURE

Day One

1. Examine your LB/AMP and LB/AMP/ARA plates from the transformation lab with and without UV light. To prevent damage to your skin or eyes, avoid exposure to the UV light. Never look directly into the UV lamp. Wear safety glasses whenever possible.
2. Obtain a culture tube containing 5 mL of sterile LB/AMP/ARA medium and one culture tube containing 5mL of LB/AMP. Using a sterile inoculation loop, lightly touch the "loop" end to a single green colony and gently scoop up the cells without gouging the agar. Immerse the loop in the LB/AMP/ARA culture tube. Spin the loop between your index finger and thumb to disperse the entire colony.
3. Using a new sterile loop, repeat and immerse it in the LB/AMP culture tube. What do you predict will be the result tomorrow when you look at both tubes with a UV light?
4. Place them in a 37° C incubator with vigorous shaking for 18-24 hours. After incubation store at 4°C until extraction.

PROCEDURE DAY TWO – ISOLATING GFP USING HIC CHROMATOGRAPHY

The detailed procedure for part II is on the next few pages. Using this procedure create a quick guide protocol for part II on the next page. Part I is already done for you as an example.

PROTOCOL	EXPERIMENTAL NOTES
<p>PART I: CELL LYSIS</p> <ol style="list-style-type: none"> 1. Collect culture tubes and invert to mix. 2. Observe tubes in normal room lighting and then with the UV light. Record 3. Pour a 2mL aliquot of the LB-amp-ara culture into a clear 2mL microcentrifuge tube. 4. Centrifuge the tube for 1 minute at 12,000-14,000 rpm. Discard supernatant and pour an additional 2mL aliquot to the same tube and centrifuge as above. Discard supernatant in liquid biohazard waste. 5. Add 250 µl of TE Solution and resuspend pellet thoroughly by vortexing. 6. Add one drop of lysozyme by pipet to the resuspended pellet. Cap and mix by inversion. 7. The lysozyme will start digesting the bacterial cell wall. Observe the tube under the UV light. 8. Place the microcentrifuge tube in a -80 °C freezer for 5 min. 9. Thaw for 1 min at 37°C. 10. Repeat freezing two more times. 11. Centrifuge for 5 minutes at 12,000-14,000 rpm. 12. Transfer 50 µl of the supernatant into a new tube labeled "Protein Extract." 13. Transfer the remaining 200 µl of the supernatant into a new tube labeled "supernatant." 	<p>Culture incubation conditions: _____</p> <p>LB-amp-ara: _____</p> <p>LB-amp: _____</p> <p>Pellet observations: _____</p> <p>Lysis Observations: _____</p> <p>Time 1: Time 2: Time 3: *prepare column while waiting</p> <p>Observations with UV light:</p>



Assignment: Before class! Read the detailed procedure on the following page and write a short protocol below

PART II: COLUMN CHROMATOGRAPHY

Collection Tube Number	Observations Under UV Light
Tube 1 Sample in Binding Buffer	
Tube 2 Sample with Wash Buffer	
Tube 3A Sample with TE Buffer	
Tube 3B Sample with TE Buffer	
Tube 3C Sample with TE Buffer	

DETAILED PROTOCOL

1. Remove your two liquid cultures from the incubator (or fridge). If they have been sitting invert to mix.
2. Observe tubes in normal room lighting and then with the UV light. Note any color differences that you observe. Is this what you predicted?
3. Transfer the entire contents of the GFP glowing liquid culture by pipet into a clear, 2 mL microcentrifuge tube labeled GFP. You will need to do this in 2 to 3 aliquots into the same tube as follows:
 - a. Mix bacterial culture tube by inversion a few times and pour a 2mL aliquot into a clear 2mL microcentrifuge tube.
 - b. Centrifuge the tube for 1 minute at 12,000-14,000 rpm. Be sure to balance the tubes in the centrifuge.
 - c. After centrifugation observe the pellet and supernatant under UV light. Note your observations in your notebook.
 - d. Discard supernatant and pour an additional 2mL aliquot to the same tube and centrifuge as above.
4. After centrifugation, open the tube and slowly pour off the liquid supernatant into a liquid biohazard waste container.
5. Pour more of the liquid bacterial culture into the SAME tube and centrifuge as above. It is recommended to use all 5mL of culture.
6. Add 250 μ L of TE Solution to the tube. Resuspend the bacterial pellet thoroughly by vortexing (don't use a micropipette, the bacteria will get stuck in the tip).
7. Add 1 drop of lysozyme by pipet to the resuspended pellet. Cap and mix the contents by flicking the tube with your index finger. The lysozyme will start digesting the bacterial cell wall. Observe the tube under the UV light. Place the microcentrifuge tube in the freezer. Freezing will cause the cells to rupture completely. Freezing *slowly* at -20°C for 30 min-1 hr works efficiently at cell rupture. If time is limiting, you can use a -80°C freezer for 15-30 minutes. When placing in freezer, do not put in a rack, it will prevent freezing.
Recommended: 2x20min freeze-thaw cycles at -20°C . If there is not enough time for this, 3x5min freeze-thaw cycles at -80°C also work very well.
8. Remove your tube from the freezer and thaw using hand warmth. Place the tube in the centrifuge and pellet the insoluble cell debris by spinning for 5 minutes at 12,000-14,000 rpm. Label a new microcentrifuge tube.
9. While your sample is freeze-thawing, prepare the chromatography column. Shake the column to resuspend the beads. If the beads have dried out, add 1mL of Equilibration Buffer and shake.

10. Shake the column down to consolidate the beads at the bottom. Remove the top cap and snap off the tab bottom of the chromatography column. Allow the liquid buffer to drain from the column into a beaker (this will take 3–5 minutes).
11. Equilibrate the column by adding 2 mL of Equilibration Buffer to the top of the column, 1 mL at a time, being careful not to disturb the bead bed at the top. Drain the buffer.
NOTE: If the column is not dripping notify your instructor. Sometimes the tip (2mm) may need to be cut off with a pair of scissors to open up the hold at the bottom of the column.
12. After the centrifugation, immediately remove the microcentrifuge tube from the centrifuge. Examine the tube with the UV light. The cell debris should be visible as a pellet at the bottom of the tube. The liquid that is present above the pellet is called the supernatant. Note the color of the pellet and the supernatant.
13. Transfer 50 µl of the supernatant into a new tube labeled “Protein Extract.” We will run this aliquot on the SDS-PAGE gel to determine how well our purification worked to separate GFP from the other proteins in the extract.
14. Transfer the remaining 200 µl of the supernatant into a new tube labeled “supernatant.”
15. Transfer 200 µl of Binding Buffer to this tube containing the supernatant. Mix by pipetting up and down a few times.
16. Obtain five clear 1.5mL collection tubes and label them as shown in the table below. Include the following data table in your lab notebook:

Collection Tube Number	Prediction under UV light	Observations Under UV Light
Tube 1 Sample in Binding Buffer		
Tube 2 Sample with Wash Buffer		
Tube 3A Sample with TE Buffer		
Tube 3B Sample with TE Buffer		
Tube 3C Sample with TE Buffer		

17. When the last of the equilibration buffer has drained from the HIC column bed, move the column to collection tube 1.

18. Carefully load all of the supernatant (in Binding Buffer) into the top of the column by resting the pipette tip against the side of the column and letting the supernatant gently drip down the side of the column wall. Examine the column using the UV light. Note your observations in the data table. Let the entire volume of supernatant flow into the column and tube 1.
19. Move the column to collection tube 2. VERY CAREFULLY, slowly add 250 μ l of Wash Buffer and let the entire volume flow into the column. As you wait, predict the results you might see with this buffer. Examine the column using the UV light and record your results.
20. Move the column to tube 3A. VERY CAREFULLY, slowly add 250 μ l of TE buffer (Elution Buffer) and let the entire volume flow into the column. Again, make a prediction and then examine the column using the UV light. List the results in the data table.
21. Move the column to tube 3B. VERY CAREFULLY, slowly add 250 μ l of TE buffer (Elution Buffer) and let the entire volume flow into the column. Again, make a prediction and then examine the column using the UV light. List the results in the data table.
22. And finally, move the column to tube 3C. VERY CAREFULLY, slowly add 250 μ l of TE buffer (Elution Buffer) and let the entire volume flow into the column. Again, make a prediction and then examine the column using the UV light. List the results in the data table.
23. Examine all of the collection tubes using the UV lamp and note any differences in color between the tubes. Record in your notebook.
24. Cap the collection tubes and label with your initials. Place tubes (and pre-column aliquot) in a labeled freezer box. Store at -20°C until next lab.

LAB UNIT 10: DENATURING POLYACRYLAMIDE GEL ELECTROPHORESIS

Jack O'Grady, M.S., Linnea Fletcher, Ph.D., Angela Wheeler, M.S., Patricia Phelps, Ph.D.
Adapted from Bio-Rad's pGLO SDS-PAGE extension

OBJECTIVES

Your performance will be satisfactory when you can:

- ✓ Explain how polyacrylamide gel electrophoresis can be used to analyze proteins
- ✓ Analyze protein samples using PAGE; load, run and stain a polyacrylamide gel
- ✓ Analyze and interpret the results of a PAGE
- ✓ Plot molecular weight marker with MS Excel and using linear regression to determine the size of GFP in chromatography fractions

INTRODUCTION

General Principles of Protein Electrophoresis and SDS-PAGE

Electrophoresis ("to carry with electricity") is the migration of charged molecules in an electric field toward the electrode with the opposite charge. This technique is widely used in molecular biology research to answer a variety of questions by examining proteins. For example:

- What proteins are in my sample?
- What are the molecular weights of the proteins?
- How pure is my protein of interest?
- How much protein do I have?

Ulrich Laemmli developed his system of polyacrylamide gel electrophoresis with two gel phases so that all of the proteins in a gel begin separating or resolving, at the same time. Since sample volumes can vary from lane to lane, forming vertically narrow or broad bands in the wells, all of the proteins in a sample do not enter the stacking gel zone simultaneously. However, the low percentage (4%) of the stacking gel allows the proteins to migrate rapidly and accumulate at the edge of the denser resolving gel, regardless of their sizes. The samples of mixed proteins are thus concentrated into uniformly thin bands in each lane before they move into the denser (5-20%) resolving gel and begin to separate according to their weights.

There is no obvious visual border between the stacking and resolving zones of a commercially prepared gel, but if you watch your samples immediately after turning on the power supply, you will see the protein samples being focused into a narrow band at the interface. Prestained protein markers first stack into a tight band, and then the individual prestained proteins become distinct as the proteins begin to separate according to their molecular weights.

Why are we using polyacrylamide gels, and not agarose gels, to analyze proteins?

The gel matrix formed by polyacrylamide is much tighter and able to resolve much smaller molecules than agarose gels. Polyacrylamide gels have pore sizes similar to the sizes of proteins. Nucleic acids are orders of magnitude larger than proteins, and agarose is usually the preferred medium for these molecules. However, when separating very small fragments of DNA, for



example during DNA sequencing, polyacrylamide is the matrix of choice.




Figure 1. Precast gels are very thin polyacrylamide gels sandwiched between clear plates. Each gel has two separate zones, the stacking gel and the separating gel, which is also known as the resolving gel. In polyacrylamide gel electrophoresis, samples are loaded into wells at the top of the stacking gel, and the proteins move downward toward the positively charged electrode.

The Chemistry and Physics behind Electrophoresis

The size of biomolecules is expressed in Daltons (D), a measure of molecular weight. One Dalton equals the mass of a hydrogen atom, which weighs 1.66×10^{-24} gram. Most proteins have masses on the order of thousands of Daltons, so the term kilo-Dalton (kDa) is used for molecular weight. Proteins range in size from several kD to thousands of kD. In contrast, the nucleic acids we study are often larger than 1000 base pairs, or 1 kilobase (kb), and each kilobase pair has a mass of approximately 660 kD. For example, when cloning DNA, a 2 kb fragment of DNA can be inserted into a plasmid vector of 3 kb, giving a total plasmid length of 5 kb. The mass of this 5-kb plasmid would be approximately 3.3 million Daltons or 3,300 kD, much larger than the average protein!

A molecule's electrical charge and its mass affect its mobility through a gel during electrophoresis. The ratio of charge to mass is called charge density. Since every protein is made of a unique combination of amino acids, each of which may have a positive, negative, or neutral charge, the net charge of each protein is naturally different. The inherent charges of proteins must be removed as a factor affecting migration for polyacrylamide electrophoresis to be effective as a method of protein molecular weight determination.

The intrinsic charges of proteins are obscured by placing a strongly anionic (negatively charged) detergent, SDS, in both the sample buffer and the gel running buffer. SDS coats the proteins with negative charges and also keeps them denatured as linear chains. In this form, proteins migrate in a polyacrylamide gel as if they have equivalent negative charge densities, and mass becomes the only variable affecting the migration rate of each protein. This technique is called SDS-PAGE.

Polyacrylamide Acts As a Molecular Sieve: The degree of sieving within a gel can be controlled by adjusting the polyacrylamide concentration. Higher concentrations of polyacrylamide resolve smaller molecular weight ranges. For example, a 5% polyacrylamide gel separates large proteins of 100 to 300 kD, while an 18% polyacrylamide gel is better for separating smaller proteins in the range of 5 to 30 kD.

For this lab, we will use a premade 12-15% polyacrylamide gel, which provides excellent separation of proteins in the range of 10 to 100 kD. Our attention will be focused on variations among the smaller proteins, in the range of 15 to 50 kD, since it is easier to discern differences among these proteins. Smaller proteins migrate further through the gel and are better resolved than proteins of high molecular weights.

Running a Polyacrylamide Gel: Polyacrylamide gels are pre-cast in a plastic cassette. The gel cassette is inserted into a vertical electrophoresis apparatus, and the running buffer is added until each well is covered with buffer. Samples, controls, and molecular weight markers are loaded into the wells. A lid is placed on the apparatus, and leads are plugged into a power supply. A current is applied at a constant voltage, bubbles rise from the electrodes, and the loading dye and proteins in the samples begin to enter the gel.

Sample Preparation – Disrupting Protein Structure




Figure 2: Secondary, tertiary and quaternary protein structure must be disrupted, or denatured to separate proteins by size.

To effectively determine the molecular weights of proteins, the secondary (2°), tertiary (3°), and quaternary (4°) structures of the protein complexes within a protein extract are disrupted before electrophoresis. This process of structural disruption is called denaturation.

1. Primary structure = order of amino acids
2. Secondary structure = domains of repeating structures, such as β -pleated sheets and α -helices as a result of H bonding between peptide backbone
3. Tertiary structure = 3-dimensional shape of a folded polypeptide, maintained by disulfide bonds, electrostatic interactions, hydrophobic effects, H bonding of the R groups on the amino acids
4. Quaternary structure = several polypeptide chains associated together to form a functional protein

Secondary, tertiary, and quaternary structures are disrupted by the combination of heat and

SDS. A reducing agent, such as β -mercaptoethanol (BME) or dithiothreitol (DTT), may be added to ensure complete breakage of disulfide bonds. These three factors – heat, ionic detergent, and reducing agent – completely disrupt the 2°, 3°, and 4° structures of proteins and protein complexes, resulting in linear chains of amino acids which allow the molecules to snake through the gel at rates proportional to their molecular masses.




Figure 3: A quaternary protein complex denatured with reducing agents, heat and SDS can be separated into individual proteins and resolved by size using SDS-PAGE.

Visualizing the Proteins: After electrophoresis is complete, the gel is stained so that blue-colored protein bands appear against a clear background.

Molecular Weight Standards: Electrophoresis protein standards, or molecular weight markers, consist of a mixture of proteins of known molecular weight. They are available in some protein size ranges. The markers to be used should correspond to the sizes of the proteins of interest. Molecular weight standards are available either prestained or unstained. Unstained markers are not visible until the gel is stained with a protein stain, such as Bio-Safe™ Coomassie stain. The prestained Kaleidoscope standards used in this lab are visible as they separate on the gel. The dyes bound to the Kaleidoscope marker proteins affect the migrations of the proteins, and the actual sizes of the dyed molecules differ slightly from batch to batch. Please refer to the size chart that comes with each vial of Kaleidoscope prestained standards for the calibrated molecular weights of each of the dyed proteins.




Figure 4: Kaleidoscope prestained standards


Identifying Proteins in Polyacrylamide Gels

It is not possible to definitively identify unknown proteins in an SDS-PAGE gel without additional analysis. In an experiment like this one, each protein extract contains a complex mixture of proteins. The different proteins appear as distinct blue-stained bands on the gel. From the positions and intensities of these bands, we can determine the size and relative abundance of the proteins, but we can only make educated guesses about the identity of each protein, based on available references.

Even when the molecular weight of a protein is known, and used as a criterion for identification, there are two possible sources of error. First, bands that migrate almost identically on a gel may be different proteins of very similar sizes. Second, proteins of very similar composition, function, and evolutionary origin may be different in molecular weight, because of variations acquired as they evolved. Definitive identification of a protein requires mass spectrometry, sequencing, or immunodetection. Immunodetection methods, such as western blotting, use antibodies that specifically recognize the proteins of interest. Such antibodies can provide identification when bands cannot be identified by molecular weight alone.

Using Molecular Weight Marker to Determine MW of Unknown

There is a linear relationship between the log of the molecular weight of the protein and migration distance (R_f) when proteins are separated on a denaturing polyacrylamide gel. This relationship can be exploited to determine the molecular weight of a target protein. Using the equation of the line and the migration distance of the target protein, you can determine its approximate molecular weight. This relationship works best with denaturing gels because the three-dimensional shape and charge of the protein do not influence its migration – it migrates by size alone.



PROTEIN ANALYSIS OF GFP HIC SAMPLES USING SDS-PAGE

MATERIALS

<ul style="list-style-type: none">• 1, polyacrylamide gel per person• Chromatography samples• Kaleidoscope molecular weight marker• Sample load buffer• Gel load tips & micropipettes• Screw cap tubes• Hand-held UV light	<ul style="list-style-type: none">• 95°C water bath and floater• Vertical Gel Apparatus (one per lab group)• 1000mL, 1X Denaturing Gel Running Buffer per vertical gel apparatus• Coomassie Safe stain• Ruler• Cracker to open gel
--	---

SAFETY

1. The wires connecting the cell to the power supply must be in good condition, not worn or cracked. Broken or worn wires not only cause rapid changes in resistance that adversely affects electrophoresis, but they also create an electrocution hazard.
2. An area of at least 6 inches around the power supply and cell should be bare of clutter and other equipment and dry.
3. Wear gloves while loading and handling the gels; the unpolymerized acrylamide is a neurotoxin! Most prepared gels are preserved using sodium azide which is very toxic!
4. The sample load dye has either DTT or beta-mercaptoethanol. Both of these are hazardous chemicals, handle with gloves, and avoid breathing in vapors and dispose of tips and tubes in hazardous waste containers. Keep tubes closed when aliquoting.
5. Coomassie blue will stain clothing and hands. Wear gloves when handling the staining and destaining solutions.

PROTOCOL

Part I: Preparing Protein Samples

1. Label 7 screw-cap tubes A-G. You will prepare your protein samples you collected from the column chromatography lab in these tubes. Use screw cap tubes.
 - A. Protein Extract undiluted
 - B. Column fraction 1 – proteins not bound to column
 - C. Column fraction 2 – proteins washed off column
 - D. Column fraction 3A – elution 1
 - E. Column fraction 3B – elution 2
 - F. Column fraction 3C – elution 3
 - G. *Column fraction with the most GFP (this tube will not be heated)*

2. Add 50 μ l of Laemmli sample load buffer to each of the tubes you just labeled. Laemmli buffer contains SDS and DTT to help denature the proteins in your sample as well as tracking dye and glycerol to weigh the samples down when you load them on the gel.

➡ Your sample buffer also has reducing agent DTT. Therefore you must *dispose of all tips and tubes in a waste beaker.*

⇒ ***WARNING! DTT is harmful*** if absorbed through the skin. It irritates skin, eyes, and respiratory tract. Wear gloves and wash hands, avoid from directly inhaling. Keep tubes CAPPED. Do not breathe in fumes.

3. Add 50 µl of each sample you collected to the corresponding tube you labeled in step 1.
4. Heat tubes A-F ONLY (***NOT tube G***) for 5 minutes at 95°C (a boiling water bath works best for this).
5. ***Molecular Weight Marker***: You most likely will be using Kaleidoscope prestained standards. This marker already has load dye. ***Do NOT heat this marker***, it will cause it to degrade! Warm to room temperature on your bench top before using. You will need 10µl per gel.
6. Prepare 1000ml, of 1X Tris/glycine Running Buffer using the 10X stock.

Calculations:

7. Your instructor will demonstrate how to set up your gel apparatus.
⇒ ***WARNING! Acrylamide is a neurotoxin! Always handle gels with gloves.***
 - a. Remove the comb and tape along the bottom of the pre-made polyacrylamide gels and place your gel in the chamber with the short plates facing inside as shown by your instructor.
 - b. Make sure your apparatus does not leak by pouring 1X TGS into the center and watch for leaking out the bottom. Insert into the gel box and fill to the mark “Two gels.”
 - c. Wash out the wells vigorously to remove unpolymerized acrylamide with a disposable transfer pipet or p1000 micropipette.

Part II: Separating Protein Samples on the SDS-PAGE gel:

GLP tips to load a sample into a well:

1. Using a thin gel loading tip, aliquot the correct amount of your sample from your microcentrifuge tube slowly.
2. Insert the micropipette tip into the top of the well to at least four mm of the bottom of the well in between the two plates. It helps to lean the tip on the tall plate and using both hands to steady tip.
3. Ensure the pipette tip is between the two plates and *very slowly and gently* expel the solution from the micropipette tip into the well while holding the micropipette steady. The blue solution should fall to the bottom of the well, gradually filling it.
4. Do not press the micropipette to the second stop – it is important to avoid blowing air bubbles into the well.

5. Do not release your thumb until you have slowly withdrawn the micropipette tip from the well so that you avoid removing the sample that you have so carefully loaded!

PROTOCOL:

1. Load 10 μ l of the Kaleidoscope prestained standard in the far left lane. The MW band sizes are 250, 150, 100, 75, 50, 37, and 25, 20, 15, 10 kDa. This marker does NOT get heated!
2. Load 30 μ l of each protein sample (A-G) into a separate well. Remember to fill in gel-documentation form!
3. Put the lid on the tank and insert leads into the power supply, matching red to red and black to black.
4. Electrophorese for 30 minutes at a constant voltage of 200 V or until the dye front reaches the bottom of the gel. Alternatively, you can use 250V for 20 min with special fast gels from Bio-Rad (ask your instructor).
5. When the loading dye has almost reached the bottom of the gels, stop the power supply and disconnect the leads. Remove the gel cassettes. Rinse with water and wipe with a paper towel. Using a hand-held UV light, look to see if you can see your GFP protein glowing. If you can, take a sharpie and mark on the plate where you see it. Measure with a ruler the migration distance of the glowing GFP and note what lane it is in. Also, Measure the blue dye front. You will need this for the graphing procedure in part IV.

- a. GFP glowing in lane(s):
 - b. The migration distance of GFP (from the bottom of the well):
 - c. The migration distance of the dye front (from the bottom of the well):
6. Lay a gel cassette flat on the bench with the short plate facing up. Carefully pry apart the gel plates using a spatula. The gel will adhere to one of the plates.
7. Transfer the plate with the gel adhering to it to a tray containing dH₂O. Rinse for 15 min with dH₂O to remove excess running buffer that may interfere with the Bio-Safe Coomassie blue stain. This is best done a three, 5-minute washes.
8. Using the UV light box, capture an image of your gel BEFORE staining. Do you see the GFP glowing in any of your lanes? Which lane(s)?


Part III: Staining SDS-PAGE Gel with Coomassie Safe Stain

1. Pour out water and add just enough Coomassie safe stain to cover the gel. Stain the gels for one hour (to overnight) with gentle agitation.
 → NOTE: when pouring out water washes, hold gel with your finger lightly, or else it will get dumped down the sink!
2. After the gels have been stained, rinse the stain down the sink with lots of water.
3. Cover the gel with large volumes of dH₂O, changing several times until most of the blue background disappears and you can visualize discrete bands on the gel. This can take several hours to overnight.
4. The GFP protein is approximately 27 kDa. Do you see a band in this range for any of your samples?

Part IV: Gel Analysis

1. Measure the migration distance (in mm!) of each band in the molecular weight marker lane. Record this in the data table below/

Distance of Dye front: _____ mm



2. Calculate the Rf value by dividing the migration distance of the band by the migration distance of the blue dye front. Record this in the table.
3. Calculate the log of the molecular weight of each band. Record in the table.
4. Plot the log of the molecular weight (y-axis) versus the Rf value (x-axis).
5. Determine the equation of the line and the R² value.
6. Measure the distance of the predicted GFP in each well. Create a table to record your sample data.
7. Using the equation of the line determine the molecular weight of the major protein bands in each lane. Record this in the same table.



LAB UNIT 11: DNA BARCODING - AMPLIFY DNA BARCODE USING PCR

Jack O'Grady, M.S. Adapted from Cold Spring Harbor DNA Learning Center

This lab has been modeled from in part or whole by the Cold Spring Harbor DNA Learning Center DNA Barcoding 101 Manual: "Using DNA Barcodes to Identify and Classify Living Things." 2012 <http://www.dnabarcoding101.org/>

OBJECTIVES

Your performance will be satisfactory when you can:

- ✓ set up a PCR reaction with all required components
- ✓ operate and program a Thermocycler
- ✓ analyze DNA on an agarose gel
- ✓ interpret the gel data and evaluate the quality of amplification as well as determine the size of the amplicon

INTRODUCTION

Polymerase chain reaction (PCR) is a technique for rapidly creating multiple copies of a segment of DNA utilizing repeated cycles of DNA synthesis. PCR has revolutionized molecular biology and forensics, allowing amplification of small quantities of DNA into amounts that can be used for experimentation or forensic testing. Kary Mullis, who later won a Nobel Prize for his work, developed PCR in 1983. The subsequent discovery of a DNA polymerase that is stable at high temperatures and the introduction of thermal cyclers, instruments that automate the PCR process, brought the procedure into widespread use in the late 1980s.

From trace amounts of the DNA used as starting material (template), PCR produces exponentially larger amounts of a specific piece of DNA. The template can be any form of DNA, and only a single molecule of DNA is needed to generate millions of copies. PCR makes use of two normal cellular activities: 1) binding of complementary strands of DNA, and 2) replication of DNA molecules by DNA polymerases.

PCR Step by Step

The strength of PCR lies in its ability to make many copies of (amplify) a single region (target) of a longer DNA molecule. For example, a researcher wanting to study a single human gene needs to amplify only that portion from the enormous human genome of approximately 3.3×10^9 base pairs! The first step is to identify and sequence areas upstream and downstream from the DNA of interest. Once this is done, short strands of DNA that are complementary to the upstream and downstream DNA are synthesized. As in cellular DNA replication, these oligonucleotide primers are used as the starting point for copying the DNA of interest, but the primers used in PCR are DNA oligonucleotides, not RNA. PCR involves a repetitive series of cycles, each of which consists of template denaturation, primer annealing (binding to the template DNA strand), and extension of the annealed primer by a heat-stable DNA polymerase. All of the components needed for PCR are mixed in a microcentrifuge tube.



PCR Reaction Components:

1. Template DNA
2. *Taq* DNA polymerase (or another thermally-stable DNA polymerase)
3. Primers — synthesized to complement a specific region on the template DNA. The two primers in a pair are designed to anneal to opposite ends of the region of interest. The primers are added in excess (that is, there are many more primer molecules than template molecules in the reaction tube)
4. Nucleotides — the four individual bases in the form of deoxynucleoside triphosphates (dNTPs), which allows them to be added to a DNA polymer. The dNTP mixture includes the same amounts of dATP, dTTP, dGTP, and dCTP.
5. Reaction buffer — prepared with the correct ionic strength of monovalent and divalent cations needed for the reaction and buffered to maintain the pH needed for enzyme activity

The microcentrifuge tubes are specialized tubes used only for PCR. PCR tubes are plastic with very thin walls, allowing rapid transfer of heat through the plastic, and the tubes usually hold only 0.2 or 0.5 ml.

The PCR reaction tubes are placed in a thermal cycler, an instrument developed in 1987 that automates the heating and cooling cycles needed during PCR. Thermal cyclers contain a metal block with holes for the PCR tubes. The metal block can be heated or cooled very rapidly. Thermal cyclers are programmable, so the PCR reaction parameters (temperatures, time at each temperature, and a number of cycles) can be stored by the instrument. This means that the user can just load the samples and push a button to run the reactions. Contrast this to the early researchers who had to sit by a series of water baths with a timer, switching the tubes from one temperature to another manually for hours!

The first step of the PCR reaction is the *denaturation step*. Since DNA polymerase can use only single-stranded DNA as a template, the first step of PCR is uncoiling and separating the two strands of the template DNA. In cells, enzymes such as helicase and topoisomerase do this work, but in PCR, heat is used to separate the strands. When double-stranded DNA is heated to 95°C, the strands separate, or denature. Since complete denaturation of the template DNA is essential for successful PCR, the first step is frequently an extended denaturation period of 2–5 minutes. The initial denaturation is longer than subsequent denaturation steps because the template DNA molecules are longer than the PCR product molecules that must be denatured in subsequent cycles. Denaturation steps in subsequent PCR cycles are normally 30–60 seconds.


The thermal cycler then rapidly cools the reactions to 40–60°C to allow the primers to *anneal* to the separated template strands. The temperature at which the primers anneal to the template DNA depends on several factors, including primer length, the G–C content of the primer, and the specificity of the primer for the template DNA. If the primer sequences match the template sequences exactly, the primers will anneal to the template DNA at a higher temperature. As the annealing temperature is lowered, primers will bind to the template DNA at sites where the two strands are not exactly complimentary. In many cases, these mismatches will cause the strands to dissociate as the temperature rises after the annealing step, but they can also result in

amplification of DNA other than the target. In the annealing step, the two original strands may reanneal to each other, but the primers are in such excess that they out-compete the original DNA strands for the binding sites.

The final step is an *extension*, in which the reaction is heated to 72°C, the optimal temperature for *Taq* DNA polymerase to extend the primers and make complete copies of each template DNA strand.

At the end of the first PCR cycle (one round of denaturation, annealing, and extension steps define one cycle), there are two new strands for each original double-stranded template, which means there is twice as much template DNA for the second cycle of PCR. As the cycle is repeated, the number of strands doubles with each reaction. For example, after 35 cycles, there will be over 30 billion times more copies of the target sequence than at the beginning. The number of cycles needed for amplification depends on the amount of template DNA and the efficiency of the reaction, but reactions are frequently run for 30-40 cycles.

PCR generates DNA of a precise length and sequence. During the first cycle, primers anneal to the original template DNA strands at opposite ends and on opposite strands. After the first cycle, two new strands are generated that are shorter than the original template strands but still longer than the target DNA, because the original template sequence continues past the location where the other primer binds. It isn't until the third PCR cycle that fragments of the precise target length are generated.



Example of thermal cycling profile. In this profile an initial denaturation step of 95°C for 5 min is followed by 40 cycles of 1-minute denaturation, 1-minute annealing, and 2-minutes extension. A final 6-minute extension time is added to ensure completion of DNA synthesis. The final hold ensures samples are kept stable until the samples are retrieved.



Watch PCR animation: <https://youtu.be/2KoLnIwoZKU>

Write the exact PCR protocol the class will program into the thermocycler. For each step, explain its purpose.

Temperature	Time	Purpose

PCR STERILE TECHNIQUE

PCR is a powerful and sensitive technique that enables researchers to produce large quantities of specific DNA from very small amounts of starting material. Because of this sensitivity, contamination of PCR reactions with unwanted DNA is always a possible problem. Therefore, utmost care must be taken to prevent cross-contamination of samples. Steps to be taken to prevent contamination and failed experiments include:

- Aerosol barrier pipet tips. The end of the barrels of micropipettes can easily become contaminated with aerosolized DNA molecules. Pipet tips that contain a filter at the end can prevent aerosol contamination from micropipettes. DNA molecules that are found within the micropipette cannot pass through the filter and contaminate PCR reactions.
- Aliquot reagents. Sharing of reagents and multiple pipetting into the same reagent tube can easily introduce contaminants into your PCR reactions. When at all possible, divide reagents into small aliquots for each team, or if possible, for each student. If only one aliquot of a reagent does become contaminated, then only a minimal number of PCR reactions will become contaminated and fail.
- Change pipet tips. Always use a new pipet tip when entering a reagent tube for the first time. If a pipet tip is repeatedly used, contaminating DNA molecules on the outside of the tip will be transferred to other solutions, resulting in contaminated PCR reactions. If you



are at all unsure if your pipet tip is clean, err on the safe side and discard the tip and get a new one. The price of a few extra tips is a lot smaller than the price of failed reactions.

- Use good sterile technique. When opening tubes or pipetting reagents, leave the tubes open for as little time as possible. Tubes that are open and exposed to the air can easily become contaminated by aerosolized DNA molecules. Go into reagent tubes efficiently, and close them as soon as you are finished pipetting. Also, try not to pick tubes up by the rim or cap, as you can easily introduce contaminants from your fingertips.
- Clean Area. Bleach at a concentration of 10% destroys DNA, so wiping down surfaces and rinsing plastic pipet barrels, mortars, and pestles with 10% bleach can get rid of any surface DNA contamination that may arise.



TEST YOUR
KNOWLEDGE!

List and explain the components in a PCR reaction tube:

Why is sterile technique important for PCR?

References:

1. Bio-Rad Cloning & Sequencing Manual
2. Cold Spring Harbor Laboratory's DNA Barcoding 101: <http://www.dnabarcoding101.org/>

Lab Unit 11-A: PCR Amplification of DNA Barcode

REAGENTS, SUPPLIES, & EQUIPMENT:

For each group of 2:

- Container with crushed ice
- Appropriate primer/loading dye mix (25 µL)* per reaction
- DNA from each lab partner specimen(s) (from plant gDNA isolation lab)*
- Micropipettes and barrier tips (1-100 µL)
- Microcentrifuge tube rack
- Permanent marker
- ready-to-go PCR Beads in 0.2- or 0.5-ml PCR tube per reaction (3 per group of 2)

To share:

- Thermal cycler

The following primers are available to use for this experiment:

Plant rbcL gene

rbclA F 5'- ATGTCACCACAAACAGAGACTAAAGC-3' (forward primer)

rbclA R 5'- GTAAAATCAAGTCCACCRCG-3' (reverse primer)

Animal coi gene

lepF1 5'- ATTCAACCAATCATAAAGATATTGG -3' (forward primer)

lepR1 5'- TAAACTTCTGGATGTCCAAAAAATCA-3' (reverse primer)

vf1f 5'- TCTCAACCAACCACAAAGACATTGG-3' (forward primer)

vf1r 5'- TAGACTTCTGGGTGGCCAAAGAATCA-3' (reverse primer)

Plants: rbcL primers (rbclA F / rbclA rev)

Fish: COI primers (VF2_t1/ FishF2_t1/ FishR2_t1/ FR1d_t1)

Insects: (LepF1_t1/ LepR1_t1)

Other animals: (LepF1_t1/ VF1_t1/ VF1d_t1/ VF1i_t1/ LepR1_t1/ VR1d_t1/ VR1_t1/ VR1i_t1)

*****ASSIGNMENT!*****

Create and type up a protocol for both 11A and 11B of this lab exercise before coming to class.

Include the following:

- a. Title page with title, name, date
- b. List of materials
- c. Safety Considerations
- d. Short quick guide protocol (with enough detail to complete the exercise)
- e. Include any data tables needed

PROTOCOL

1. Obtain 3 PCR tubes containing Ready-To-Go PCR Bead. Label one tube with your initials, one tube with your lab partners initials and one tube (-).
2. Use a micropipette with a fresh tip to add 23 μL of the primer/loading dye mix to each tube. Allow the beads to dissolve for 1 minute.
3. Use a micropipette with a fresh tip to add 2 μL of your DNA directly into the appropriate primer/loading dye mix. Ensure that no DNA remains in the tip after pipetting.
4. Use a micropipette with a fresh tip to add 2 μL of nuclease-free water to the (-) tube.
5. Store your sample on ice until your class is ready to begin thermal cycling.
6. The instructor will show you how to program the Thermocycler. Ensure the cycling is set up with the following parameters:
 - Initial Denaturing step: 94°C for 3 minutes
 - 35 cycles of the following profile:
 - Denaturing step: 94°C 30 seconds
 - Annealing step: 54°C 45 seconds
 - Elongation step: 72°C 45 seconds
 - Final extension step 72°C 4 minutes
 - Hold 4-15°C
7. After thermal cycling, store the amplified DNA on ice or at -20 °C until you are ready to continue with gel analysis.

Lab Unit 11-B: Analyze PCR by Gel Electrophoresis

REAGENTS, SUPPLIES, & EQUIPMENT

Materials required:

Agarose	Balance, weigh boats
1X TAE gel running buffer	SYBR Safe stock solution, 10,000X
Gel-casting apparatus	Ultraviolet (UV) light box (or imaging system)
Gel electrophoresis box	125mL Erlenmeyer flask
Power supply	PCR reactions
Microwave oven	
Micropipette and tips (1–100 μ l)	
DNA ladder already prepared in load dye (Promega, G7541)	
Gel electrophoresis chamber and power supply	
Gel-casting tray and comb	

Part I: Analyzing PCR on Agarose Gel

1. Set up a gel casting apparatus as directed by your instructor. Ensure the bumpers are snug in place to avoid leaking during casting. Insert a comb into the tray, choosing a comb by the size and number of wells it will create. Remember to include your molecular weight marker when deciding on the comb size. Set the apparatus on a flat surface that will be undisturbed while the gel is solidifying.
2. Prepare enough 1X running buffer for the gel and the electrophoresis chamber. This will depend on the specific apparatus you will use. An easy way to determine this is to fill the chamber with water and measure the volume. Verify with your instructor that you will need a total volume of 400mL of 1XTAE.

Calculation: How will you prepare 400mL of 1XTAE from the 50XTAE stock you prepared earlier in the semester?

3. Weigh out the required amount of agarose and add it to the appropriate amount of 1X gel running buffer in an Erlenmeyer flask. The flask should be at least twice the volume of the buffer and no more than 5X the volume. For example, to prepare 30 mL of a 2% agarose gel, use a 125-mL flask.

Calculation: How much agarose will you need to prepare 30mL of a 2% agarose gel?

4. Add SYBR Safe stock solution (10,000X) to a final concentration of 1X.

Calculation: How much SYBR safe will you add to your flask of 30mL of 1XTAE to achieve a final concentration of 1X using a 10,000X stock?


5. Place the flask on a level surface and carefully mark the glass at the fluid level with a permanent marker (DO NOT mark the white marking area with a permanent marker! It cannot be removed).
6. Heat the mixture until all agarose has dissolved. A hot plate or microwave oven can be used; using a hot plate will take more time. Interrupt the heating at regular intervals and swirl the container to mix the contents. The solution should be brought to a boil, but do not allow the solution to boil over. Microwave for a minimum amount of time to avoid buffer evaporation, which will cause a dramatic increase in the percentage of the gel. For example, a 30mL solution will require less than 1min total in a typical microwave.

⇒ CAUTION! The flask will be hot! Use hot hands to remove flask and swirl. Do not hold over the top of your face while swirling.

7. When the agarose is completely dissolved, the solution will be completely clear and homogeneous; that is, you will not observe any granules or threads of non-dissolved agarose in the solution. Observe the fluid level about the mark you made on the flask. If a significant amount of water has evaporated, carefully add water to return to the level of the mark and swirl the solution.
8. Cool the solution to 50 – 60°C. The solution must be cooled to below 60°C to prevent damage to the plastic casting trays. This is typically cool to the touch. Do not cool too long or it will solidify in the flask.

⇒ NOTE: If you make a mistake with this solution DO NOT dump out down the sink. The agarose will solidify in the sink and clog it. Ask your instructor where to dispose of the solution. Usually, the trash can is fine.

9. Pour the gel immediately into a level gel casting stand. Allow the gel to form completely; typically, 20 minutes at room temperature is sufficient. Remove the comb and bumpers from the gel, place the gel in the electrophoresis chamber, and add a sufficient volume 1X gel running buffer to cover the surface of the gel (use the remaining 1XTAE buffer).

10. Place the gel in an electrophoresis chamber and cover the gel with 1X TAE buffer. The buffer level should be about a few millimeters above the top of your gel, and it should fill all the wells.
11. DNA Ladder (molecular weight marker): This should already be prepared for you. Thaw, pop-spin in a microfuge and store in a rack on your bench. Load 6ul of the DNA ladder in the first lane of your gel.
12. PCR amplicon samples: Your PCR reactions already have sample load dye in them. They are ready to go. Load 5ul of each of the PCR reaction directly into subsequent lanes of your gel (Lane 2, 3, & 4). Do not skip lanes, and remember your negative control!
 KEEP YOUR PCR AMPLICON! This is your sample you will send to sequencing!
13. Place the cover on the chamber in the correct orientation and connect leads to a power supply. Set the power supply to approximately 80 volts and allow to electrophorese approximately 1 hour for 2% minigels.
14. Fill out a gel electrophoresis documentation form.
15. Capture an image on the UV trans illuminator and print out a copy. Affix a copy of your gel picture to the gel electrophoresis form and turn this in with your lab report. Save a copy as a PDF ("print to pdf") to print out a larger copy for your molecular weight standard curve.

PART II: ANALYZING GEL

1. Measure the distance migrated by each band you can distinguish individually (only) in your molecular weight marker lane. Note that the 1000 & 3000 bp bands are bright. Record this data in a table below.

Molecular Weight Marker (Kb)	Log MW	Migration Distance (mm)
250		
500		
750		
Bright Band - 1000		
1500		
2000		
2500		
Bright Band - 3000		
4000		
5000		
6000		
8000		
10,000		

- Using MS Excel, graph your molecular weight standard curve. ***Plot Log₁₀ molecular weight on the y-axis and migration distance on the x-axis.*** Label you axis and give the graph and appropriate descriptive title (NOT log MW vs. migration distance!).
- Using MS Excel, determine the equation of the line for the linear points of the curve. The highest MW marker may not be in the linear range, so you can exclude it if it falls off the linear portion of your graph.

Equation of the Line: _____

- Determine the R² value for your graph. This is a correlation coefficient that will tell you how well your data correlates linearly to each other. An R² value >0.95 is acceptable, but >0.98 is preferred.

R² Value: _____

- Measure the distance of your PCR amplicon(s). Use the equation of the line to determine the molecular weight of each of your samples.

Show Calculations:

Size of PCR amplicon: _____ bases

- Generate a table with the class data. In the table include the researcher's name, the plant used, a photo of the plant, the size of the PCR amplicon.
- Include your tables and graph in your lab report.**

Lab Unit 11-C: Preparing Plates for Sequencing

You can only send your amplicon for sequencing if you observed only a SINGLE amplicon on your agarose gel. You will not need to add any sequencing reagents or primers. The sequencing reagents and primers used (M13 forward and reverse) will be added by the sequencing company. Your instructor will help you fill the sequencing plate and key.

1. Thaw your amplicon at room temperature, pop spin.
2. You will need 20 ul. Using a p20 micropipette, measure and ensure there is 20ul. If you are short, add the difference with the sterile water provided.
3. The class must create a plate key (diagram) and assign a place where each student will aliquot their PCR reactions. Each student must prepare two aliquots; 10ul for the forward sequencing reaction, 10ul for the reverse reaction.
4. Aliquot 10ul of your PCR reaction into each of two plate strips as determined by the plate key the class designed.
5. Cover the strips and give the sequencing strips and the plate key to the lab technician. Keep a copy of the key for your records.

Lab Unit 11 – DNA Barcoding: PCR Amplification

1. What is the name of the plant you performed the DNA barcoding experiment?
2. Was your DNA barcoding experiment successful? How do you know?
3. What is the size of your PCR amplicon? Is this what you expected? Explain.
4. Are the PCR amplicons for the whole class the same size? Is this what you expect? Why or why not?

CONCLUSION:



LAB UNIT 12: PCR-BASED VNTR ANALYSIS OF HUMAN DNA

Jack O'Grady, M.S., Linnea Fletcher, Ph.D., Angela Wheeler, M.S., Patricia Phelps, Ph.D.
Adapted from Edvotek's PCR-based VNTR Human DNA Typing (Cat#334) <http://www.edvotek.com/334>

OBJECTIVES

Your performance will be satisfactory when you can:

- ✓ Demonstrate safe handling of biohazard material
- ✓ Explain how PCR can be used in identifying individuals or alleles
- ✓ extract DNA from your own cells
- ✓ set up a PCR reaction with all required components and operate a Thermocycler
- ✓ analyze amplicon using agarose gel electrophoresis
- ✓ interpret and evaluate the quality of amplification as well as the size of the amplicon

INTRODUCTION

Although human DNA from separate individuals is identical in more places than it is unique, many regions of the human genome exhibit a great deal of diversity. Such sequences are termed **polymorphic** (having many forms) and are used for diagnosis of genetic disease, forensic identification, and paternity testing. Many polymorphisms are located in the estimated 98% of the human genome that does not code for proteins. Since no genes that encode proteins are found in these regions, changes, or **mutations**, in these regions do not have an effect on the individual and are more likely to be passed on to offspring. Mutations in protein coding regions are far more likely to be detrimental to the health and longevity of the individual.

CODIS: In 1990 the Federal Bureau of Investigation (FBI) established the Combined DNA Index System (CODIS). It is the core of the national DNA database and developed specifically to enable forensic DNA laboratories to create searchable DNA databases of authorized DNA profiles. The CODIS software permits laboratories throughout the country to share and compare DNA data. This system permits comparison of crime scene DNA to DNA profiles in a convicted offender and a forensic (crime scene) index. The **NDIS** (National DNA Index System) is the national level of CODIS containing the DNA profiles contributed by federal, state, and local participating forensic laboratories. The NDIS contains over 13,041,408 offender profiles, 2,860,423 arrestee profiles, and 804,902 forensic profiles as of September 2017. As of September 2017, CODIS has produced over 392,684 hits assisting in more than 377,507 investigations and over 66,000 arrests just in Texas! Learn more here, at the CODIS fact sheet:

<https://www.fbi.gov/services/laboratory/biometric-analysis/codis/codis-and-ndis-fact-sheet>

ISOLATING DNA: The first step in forensic DNA fingerprinting is the collection of human tissue from the crime scene or victim. These tissues include blood, hair, skin, and body fluids. The sample is treated with a detergent to lyse cell membranes and obtain DNA for further analysis.

In forensics, the polymerase chain reaction (PCR) is now used to amplify and examine highly polymorphic DNA regions. These are regions that vary in length from individual to individual and fall into two categories: 1) Variable Number of Tandem Repeats (**VNTR**) and 2) Short Tandem Repeats (**STR**). A VNTR is a region that is variably composed of a 15-70 base pair sequence,



typically repeated 5-100 times. An STR is similar to a VNTR except that the repeated unit is only 2-4 nucleotides in length. By examining several different VNTRs or STRs from the same individual, investigators obtain a unique DNA fingerprint for that individual which is unlike that of any other person (except for an identical twin).

D1S80 is a VNTR region present on chromosome 1 and contains a 16-nucleotide sequence which is variably repeated between 16 and 40 times. An individual who is homozygous for the D1S80 genotype will have equal repeat numbers on both homologs of chromosome 1, displaying a single PCR product following gel analysis. More commonly, a person will be heterozygous, with differing D1S80 repeat numbers. Amplification of DNA from heterozygous individuals will result in two distinct PCR products. In this experiment, PCR will be used to amplify a short DNA sequence from human chromosome 1 at a point called the D1S80 locus that is a variable insertion (one of the CODIS loci). The primers used to start the amplification were designed to flank the DNA region of the D1S80 insertion site. The amplicon size (position on the gel after electrophoresis) will reveal the length of the insertion.

PCR: The primer mixture you will use contains a 25 bp forward primer that starts copying one strand and a 26 bp reverse primer that starts copying the complementary strand. These primers match only one site on human DNA so only the DNA fragment between the two primers is copied. You also need *Taq* DNA polymerase, buffer, KCl, MgCl₂, and dNTP (nucleotides with each of the four bases – A, T, C, and G) in your reaction mixture to achieve amplification. All of these are supplied in a single-use, solid bead or pellet. These components, your DNA, and the primers are all of the ingredients needed to perform the amplification. PCR beads must be stored desiccated at room temperature, or they will absorb water from the air, and the enzyme will be degraded.

The thermal cycler must be programmed so that it is preheated and ready to run when your samples are ready. The program will begin with an initial 5-minute cycle at 94°C, which makes sure that the DNA completely denatures (the complementary strands pull apart). Five minutes is enough time for the 3 billion base pairs in the human genome to denaturing. Then the machine will cycle through the steps below 35 times. Each cycle doubles the amount of DNA that was produced in the previous cycle.

94°C /30 seconds: denatures amplicon DNA

65°C /30 seconds: primer anneals to complementary sequence of DNA

72°C /30 seconds: DNA polymerase begins synthesizing new DNA

After 35 cycles, the temperature is held at 72° C for 4 minutes to allow the polymerase to back-fill any amplicons that were not amplified clear to the end. After programming, make sure the thermal cycler is preheating so that it will be warm enough to start the cycles immediately.

Lab Unit 12-A: Isolating DNA & PCR of VNTR Loci


ASSIGNMENT!

Create and type up a protocol for this lab exercise – 12A - before coming to class. Include the following:

- Title page with title, name, date
- List of material
- Safety Considerations
- Protocol

PROCEDURE: PART I: ISOLATING DNA FROM CHEEK SAMPLES

Module I: Isolation of DNA from Human Cheek Cells



- LABEL** a 1.5 ml screw top microcentrifuge tube and a cup with your lab group and/or initials.
NOTE: Saline solution MUST be used for cheek cell wash. Sports drinks will inhibit amplification of DNA by Polymerase Chain Reaction in Module II.
- RINSE** your mouth vigorously for 60 seconds using 10 ml saline solution. **EXPEL** the solution into cup.
- SWIRL** the cup gently to resuspend the cells. **TRANSFER** 1.5 ml of solution into the labeled tube.
- CENTRIFUGE** the cell suspension for 2 minutes at full speed to pellet the cells. **POUR** off the supernatant, but **DO NOT DISTURB THE CELL PELLETT!** **REPEAT** steps 3 and 4 once more.
- RESUSPEND** the cheek cells in 140 µl lysis buffer by pipetting up and down or by vortexing vigorously.
- CAP** the tube and **PLACE** in a waterbath float. **INCUBATE** the sample in a 55° C waterbath for 5 minutes.
- MIX** the sample by vortexing or by flicking the tube vigorously for 20 seconds.
- INCUBATE** the sample in a 99° C waterbath for 5 minutes. *NOTE: Students MUST use screw-cap tubes when boiling DNA isolation samples.*
- CENTRIFUGE** the cellular lysate for 2 minutes at full speed.
- TRANSFER** 80 µl of the supernatant to a clean, labeled microcentrifuge tube. **PLACE** tube in ice.
- PROCEED** to Module II: Amplification of the D1S80 Locus.

STEP 4:
if cell pellet size is not large enough, REPEAT steps 3 - 4 until students can easily see the pellet.



OPTIONAL STOPPING POINT:

The extracted DNA may be stored at -20°C for amplification at a later time.

PART II: PCR OF DNA SAMPLES

For every PCR reaction you perform, it is essential to include controls to help interpret the data. For example, for this experiment, it is ideal to include a positive control with DNA that you know the primers will produce an amplicon if your DNA extraction was successful. Also, it is ideal to include a no template control (negative control) to show that the reagents are not contaminated with DNA. Depending on reagent availability, your instructor may ask you to set up positive and negative controls for your group, or the instructor may set up a positive and negative control for the class. Ask your instructor how to proceed with sample & control set up.

1. Each student PAIR will need 4 Ready-to-Go PCR beads. Label the sides of the tube as shown in the table below (not the lids!).
2. Add 20 uL D1S80 primer mix to each PCR tube and flick to dissolve. Pop-spin and place on ice. Once the bead is dissolved, it's very important to keep on ice until the samples are transferred to the PCR machine.
3. Add to each PCR tube the DNA or Nuclease-free water as outlined in the table below. Mix gently by pipetting up and down or flick and pop-spin in the picofuge and return to the ice. Make sure the bead is completely dissolved.

PCR Tube	Sample	PCR Bead	Volume Primer mix	Volume DNA	Nuclease-free Water
S1	Student 1 gDNA	1 bead	20ul	5 ul DNA	0 ul
S2	Student 2 gDNA	1 bead	20ul	5ul DNA	0ul
+	Positive Control DNA	1 bead	20ul	5ul Control DNA	0 ul
-	Negative (No template)	1 bead	20 ul	NO DNA added!	5 ul

4. Set up the PCR machine as directed by your instructor. The PCR program for this lab is as follows:
 - 94°C for 5 min
 - 35 cycles of: 94°C for 30sec, 65°C for 30sec, 72°C for 30sec
 - 72°C for 4 min
5. Keep PCR reactions on ice and wait until the thermocycler has reached 94°C. Press pause and carefully load the samples ensuring the lids are snapped in all the way, and the tubes are in the plate holes. Press resume.
6. When the program has completed, store tubes at -20°C until ready to begin the next step.

Lab Unit 12-B: Analysis of VNTR Loci PCR Amplicons

*****ASSIGNMENT***** Before coming to class write a STANDARD OPERATING PROCEDURE (SOP) for analyzing your VNTR samples on an agarose gel. You may use the following to help guide your SOP writing. See your previous gel procedures! Have your lab partner verify your SOP BEFORE class! *Use correct SOP format.*

Part I: Analyzing the VNTR PCR amplicons on an agarose gel

Prepare a 2.5% agarose gel in 1X TAE buffer. We will use a high percentage of agarose gels to ensure proper resolution of multiple alleles (amplicons) close in size that are relatively small.

Calculation: How do you prepare 400mL of 1XTAE from 50XTAE?

Calculation: How much agarose do you need for 30mL of a 2.5% concentration?

Calculation: What volume of 10,000X SYBR safe do you need

SAFETY: Are there any safety considerations?

MATERIALS: List all the materials needed – be specific!

PROCEDURE

1. Place the gel in an electrophoresis chamber and cover the gel with 1X TAE buffer. The buffer level should be well above the top of your gel, and it should fill all the wells.
2. Prepare PCR amplicon samples: Add 5 ul of 6XGel Loading solution to the PCR amplicon. Vortex briefly, pop-spin and store on a rack on your bench.
3. Prepare DNA marker: This is already prepared for you. Thaw, pop-spin in a picofuge and store in a rack on your bench.
4. Fill in gel electrophoresis documentation form. Load gels with 6uL of DNA ladder and 30ul of PCR amplicons, taking note of where each sample is loaded in the gel.
5. Place the cover on the chamber in the correct orientation and connect leads to a power supply. Set the power supply to approximately 80 volts and allow to electrophorese until the blue tracking dye is approximately 3/4th the way to the bottom of the gels (approximately 1 hour or longer for 2.5% minigels).
6. Visualize DNA bands by placing gels on a UV transilluminator and photograph gels using the gel documentation system. Affix a copy of your gel picture to the gel electrophoresis documentation form and turn this in with your lab report.

PART II: ANALYZING GEL

1. Measure the distance migrated of **ONLY each band you can resolve** in the molecular weight marker lane. Note that the 1000 & 3000 bp bands are bright. Record this data in a table below.

Molecular Weight Marker (Kb)	Log MW	Migration Distance (mm)
250		
500		
750		
Bright Band - 1000		
1500		
2000		
2500		
Bright Band - 3000		
4000		
5000		
6000		
8000		
10,000		

- Using MS Excel, graph your molecular weight standard curve. ***Plot Log₁₀ molecular weight on the y-axis and migration distance on the x-axis.*** Label you axis and give the graph and appropriate descriptive title (NOT logMW vs. migration distance!).
- Using MS Excel, determine the equation of the line for the linear points of the curve. The highest MW marker may not be in the linear range, so you can exclude it if it falls off the linear portion of your graph.

Equation of the Line: _____

- Determine the R² value for your graph. This is a correlation coefficient that will tell you how well your data correlates linearly to each other. An R² value >0.95 is acceptable, but >0.98 is preferred.

R² Value: _____

- Measure the distance of your PCR amplicon(s). Use the equation of the line to determine the molecular weight of each of your samples.

Show Calculations here:

Size of PCR amplicon: _____ bases

- Generate a table with the class data. In the table include the researcher's name and the size of the PCR amplicon(s).
- Include your tables and graph in your lab report.

LAB UNIT 13: DNA BARCODING SEQUENCING & BIOINFORMATICS

Jack O'Grady, M.S., & Cold Springs Harbor

This lab has been adapted from the Cold Spring Harbor DNA Learning Center DNA Barcoding 101 Manual: "Using DNA Barcodes to Identify and Classify Living Things." 2012 (<http://www.dnabarcoding101.org/>).

The introduction and background for this lab provided by BioRad's Cloning & Sequencing Kit Instruction Manual (Cat#166-5000EDU) and exercise instructions were provided by Cold Spring Harbor's DNA Barcoding 101 Manual.

OBJECTIVES

Your performance will be satisfactory when you can:

- ✓ Understand the principle behind the Sanger sequencing method
- ✓ Determine the quality of sequencing data
- ✓ Locate scientific publications and biological databases of DNA and protein sequences on the Internet (GenBank, BLAST, NEBcutter, ClustalW)
- ✓ Retrieve and compare sequence information from databases
- ✓ Compare evolutionary relatedness and draw phylogenetic trees from sequence comparisons
- ✓ Study protein structure and function from database comparisons

INTRODUCTION

Sequencing means determining the exact order of nucleotide bases (guanine (G), adenine (A), thymine (T), and cytosine (C)) in a DNA molecule. DNA sequencing began in the 1970s when two research groups developed different methods for sequencing, the Maxam-Gilbert method, and the Sanger method, at almost at the same time. Today, most researchers send their samples to core laboratory facilities where, for a nominal charge, their DNA is sequenced for them using an automated sequencer user the Sanger method. The researchers receive the sequence data the next day.

Sanger Sequencing Method


In Europe, Sanger and Coulson developed the chain termination method for DNA sequencing, or as they called it, the "plus and minus" method (Sanger et al. 1977). Since the mid-1980s, chain termination has been the predominant method used for sequencing, in large part because the technique could be automated. Frederick Sanger received a Nobel Prize for his work.

The steps in Sanger Sequencing, Chain Termination Method are:

1. Prepare a template of the DNA to be sequenced.
2. The following are added to a reaction tube in a microtiter plate well along with target DNA.
 - a. Sequencing primer that will start DNA synthesis at the area to be sequenced.
 - b. DNA polymerase that is heat stable



- c. Labeled nucleotides — these are deoxynucleotide triphosphates (dNTPs: dGTP; dATP; dTTP and dCTP), and they are always present in excess in the reaction.
 - d. Modified nucleotides called a dideoxynucleotide (ddNTP) (dye-terminator nucleotides) are added at a low concentration.
3. Allow DNA synthesis to proceed – this is similar to a PCR reaction: Cycling of denaturing DNA, annealing primer to target, elongation from primer by DNA polymerase. During synthesis, almost all of the nucleotides that are incorporated into the new DNA strand are labeled dNTPs, because the dNTPs are in excess. However, when a ddNTP is incorporated, DNA synthesis will stop on that strand, as there is no 3'-hydroxyl to form the next phosphodiester bond.
 4. Use capillary gel electrophoresis and fluorometry to separate the labeled DNA fragments by size, and a computer reads the fluorescent tag, processes the data and provides a sequencing chromatogram.
 5. The result is a graph called a chromatogram or electropherogram, on which the bases are represented by a sequence of colored peaks. The peak height indicates the intensity of the fluorescent signal. The automated sequencer interprets the results, assigning G, A, T, or C to each peak. If the software cannot determine which nucleotide is in a position, it will assign the letter N to the unknown base.



Example of a sequencing chromatogram

To learn more about DNA Sequencing watch this animation: <https://youtu.be/ONGdehkB8jU>

Whole Genome Sequencing. From Wikipedia (https://en.wikipedia.org/wiki/Whole_genome_sequencing)

The DNA sequencing methods used in the 1970s and 1980s were manual, for example, [Maxam-Gilbert sequencing](#) and [Sanger sequencing](#). The shift to more rapid, automated sequencing methods in the 1990s finally allowed the sequence of whole genomes. The first organism to have its entire genome sequenced was *Haemophilus influenzae* in 1995. After it, the genomes of other bacteria and some archaea were first sequenced, largely due to their small genome size. *H. influenzae* has a genome of 1,830,140 base pairs of DNA.

In 1999, the entire DNA sequence of human chromosome 22, the shortest human autosome, was published. The first plant genome - that of the model organism [Arabidopsis thaliana](#) - was also fully sequenced by 2000. By 2001, a draft of the entire human genome sequence was published. In 2004, the [Human Genome Project](#) published an incomplete version of the human

genome. To date, researchers have sequenced the complete genomes of thousands of organisms. For a complete list: https://en.wikipedia.org/wiki/Lists_of_sequenced_genomes

Analysis of DNA Sequences Using Bioinformatics Tools

The ability to determine the exact DNA sequence of genes emerged in the late 1970s and a technique to synthesize large quantities of target regions of DNA using polymerase chain reaction (PCR) was developed in the early 1980s. An electronic repository for the many genes being discovered was created in the late 1980s. This database, called **GenBank**, is operated by the **National Center for Biotechnology Information (NCBI)** and funded by the U.S. National Institutes of Health (**NIH**).

GenBank is accessible via the Internet to scientists, teachers, and students worldwide free of charge. Major efforts to completely sequence entire genomes were initiated in the 1990s and have now been completed for humans as well as for numerous model organisms studied by scientists. The capacity for isolating and sequencing genes has grown so quickly that the number of submissions to GenBank has doubled every two years since 1993 (NCBI, 2005). The challenge of analyzing all of the DNA sequences deposited in GenBank spurred the development of numerous computer programs for interpreting DNA and protein sequence data. This computer-aided analytical approach is called **bioinformatics**. In addition to GenBank, other databases housing sequence information is available, as are a wide range of software programs and tools designed to obtain, analyze, and organize this information.

DNA Sequencing Data

Once the sequencing reaction has been performed, and the samples have been analyzed on a sequencing instrument, the result is a data file that contains a chromatogram. A chromatogram is a representation of the DNA molecules generated from the Sanger chain termination sequencing protocol, where the sequence of peaks represents the sequence of bases. A chromatogram provides information on the peak intensities, the time course in which they eluted, and the base calls that the instrument made for these peaks. The data can be analyzed manually by opening the data file in a reader-style program such as **DNA Subway**, which can be found for free on the Internet.

In this laboratory exercise, we will analyze the DNA Barcoding sequences using a software program called “DNA Subway” found here: <https://dnasubway.cyverse.org/> Using this software program, you will evaluate your forward and reverse sequences, trim off any unusable portion, and then pair them up to establish a consensus of your DNA barcode. You will continue to use DNA Subway to **BLAST** your DNA barcode to determine if there is a similar sequence in the NCBI database. This will allow you to possibly identify your organism or determine if you have a novel sequence. Finally, DNA Subway will also permit you to evaluate the relatedness of your DNA Barcode to other published sequences or classmate sequences.

REFERENCES

1. BioRad Cloning & Sequencing Kit Instruction Manual (Cat#166-5000EDU)
2. Cold Spring Harbor DNA Barcoding 101 laboratory: <http://www.dnabarcoding101.org/>



PROCEDURE: **It is important that you create a DNA subway account BEFORE you come to class.

1. Log into DNA Subway at <https://dnasubway.cyverse.org/>
2. Register for an account. Do this only once. It may take up to 24hrs before it lets you log back on. If it doesn't let you on immediately wait for 24hr. Do not try to create another account. Ask your instructor for assistance if you have any difficulty. Remember to write down your password! 😊

For Part, I and II below, create an MS-word office document to capture each of the steps below. Save your MS word document as you go through the instructions.

Part I: Use BLAST to Find DNA Sequences in Databases

➤ **Your instructor may have you perform part I as a pre-lab exercise!**

Perform a BLAST search as follows:

1. Do an Internet search for "NCBI BLAST."
2. Click on the link for the result *BLAST: Basic Local Alignment Search Tool*. This will take you to the Internet site of the National Center for Biotechnology Information (NCBI).
3. Under the heading "Basic BLAST," click on "nucleotide blast."
4. Enter the primer set you used into the search window (below). These are the query sequences. Omit any non-nucleotide characters from the window because they will not be recognized by the BLAST algorithm.

The following primers were used in this experiment:

Plant rbcl gene

rbclAf 5'- ATGTCACCACAAACAGAGACTAAAGC-3' (forward primer)

rbclA rev 5'- GTAAAATCAAGTCCACCRCG-3' (reverse primer)

Animal coi gene

lepF1 5'- ATTCAACCAATCATAAAGATATTGG -3' (forward primer)

lepR1 5'- TAAACTTCTGGATGTCCAAAAATCA-3'(reverse primer)

vf1f 5'- TCTCAACCAACCACAAAGACATTGG-3' (forward primer)

vf1r 5'- TAGACTTCTGGGTGGCCAAAGAATCA-3' (reverse primer)

5. Under "Choose Search Set," select "NCBI Genomes (chromosome)" from the pull-down menu.
6. Under "Program Selection," optimize for "Somewhat similar sequences (blastn)."
7. At the bottom of the page, click on "Algorithm parameters" to expand more options.
8. Change the "Expect Threshold" to "1000"
9. Change "Word size" to "15"
10. Change "Match/Mismatch Scores" to "1,-2"
11. Click on "BLAST."

This sends your query sequences to a server at the National Center for Biotechnology Information in Bethesda, Maryland. There, the BLAST algorithm will attempt to match the primer sequences to the DNA sequences stored in its database. A temporary page



showing the status of your search will be displayed until your results are available. This may take only a few seconds or more than 1 minute if many other searches are queued at the server.

12. The results of the BLAST search are displayed in three ways as you scroll down the page.
- First, a Graphic Summary illustrates how significant matches, or “hits,” align with the query sequence. Capture a print screen of the Graphic summary (or try right click copy on the graphic and see if it copies over). **Copy into your MSword document.**

Why are some alignments longer than others?

- This is followed by Descriptions of sequences producing significant alignments, a table with links to database reports. **Copy and paste the first 5 in the table into the MSword document.**

What is the most significant alignment?

- The accession number is a unique identifier given to a sequence when it is submitted to a database, such as GenBank. The accession link leads to a detailed report on the sequence.

What is the Assession number to the most significant alignment?

- Note the scores in the “e” column on the right.
The Expectation or E value is the number of alignments with the query sequence that would be expected to occur by chance in the database. The lower the E value, the higher the probability that the hit is related to the query. For example, an E value of 1 means that a search with your sequence would be expected to turn up one match by chance.

What is the E value of your most significant hit, and what does it mean?

What does it mean if there are multiple hits with similar E values?

What do the descriptions of significant hits have in common?

- e. Next is an Alignments section, which provides a detailed view of each primer sequence (Query), aligned to the nucleotide sequence of the search hit (Sbjct, subject). Copy the aligned area into your MSword document.

Notice that hits have matches to one or both of the primers:

	<u>Forward Primer</u>	<u>Reverse Primer</u>
rbcl	nucleotides 1-26	nucleotides 27-46
Lep or VF	nucleotide 1-25	nucleotides 26-53

13. Predict the length of the product that the primer set would amplify in a PCR reaction (in vitro).

a. In the Alignments section, select a hit that matches both primer sequences.

b. Which nucleotide positions do the primers match in the subject sequence?

c. The lowest and highest nucleotide positions in the subject sequence indicate the borders of the amplified sequence. Subtracting one from the other gives the difference between the coordinates. However, the PCR product includes both ends, so add one nucleotide to the result to determine the exact length of the fragment amplified by the two primers:

d. What value do you get if you calculate the fragment size for other species that have matches to the forward and reverse primer? Why is this so?

14. Determine the type of DNA sequence amplified by the primer set:

a. Click on the accession link (beginning with “ref”) to open the data sheet for the hit used in Question above.

b. The data sheet has three parts: • The top section contains basic information about the sequence, including its basepair (bp) length, database accession number, source, and references to papers in which the sequence is published. • The bottom section lists the nucleotide sequence. • The middle section contains annotations of gene and regulatory FEATURES, with their beginning and ending nucleotide positions (“xx..xx”). These features may include genes, coding sequences (cds), regulatory regions, ribosomal RNA (rRNA), and transfer RNA (tRNA).

c. Identify the feature(s) located between the nucleotide positions identified by the primers, as determined above.

Part II: Identify Species and Phylogenetic Relationships Using DNA Subway

The following directions explain how to use the Blue Line of DNA Subway to analyze novel DNA sequences generated by a DNA sequencing experiment. If you did not sequence your own DNA sample, follow these directions using DNA sequences produced by other students.

➤ **Remember to copy your work into the MSword document as you move through the analysis.**

1. Create a DNA Subway Project and Upload DNA Sequences

- a. Log into DNA Subway at www.dnasubway.org. If you do not have an account, you will need to register first to save and share your work.
- b. Select “Determine Sequence Relationships” (Blue Line) to begin a project.
- c. Select “rbcl” or “COI” from the “Select Project Type” section. (rbcl (plant) sequences must be analyzed separately from COI (animal) sequences.)
- d. “Select Sequence Source” provides several ways to obtain sequences for barcode analysis:
 - Upload sequence(s) in ab1 (files ending with .ab1) or FASTA format. Click “Browse” to navigate to a folder on your desktop or drive containing your sequence(s). Select a sequence by clicking on its file name. Select more than one sequence by holding down the ctrl key while clicking file names. Once you have selected the sequences you want, click “Open.”
 - Enter a sequence in FASTA format. Below is an example of this format. The next line. >sequence name atgcccccttaatatgcctt.....
 - Import a sequence/trace from the DNALC. Click on your tracking number. Select one or more files from the list. Click to “Add” selected files.
 - Select a sample sequence.
- e. • Upload sequence(s) in ab1 (files ending with .ab1) or FASTA format. Click “Browse” to navigate to a folder on your desktop or drive containing your sequence(s). Select a sequence by clicking on its file name. Select more than one sequence by holding down the ctrl key while clicking file names. Once you have selected the sequences you want, click “Open.”
 - Enter a sequence in FASTA format. Below is an example of this format. The next line. >sequence name atgcccccttaatatgcctt.....
 - Import a sequence/trace from the DNALC. Click on your tracking number. Select one or more files from the list. Click to “Add” selected files.
 - Select a sample sequence.
- f. Provide a title in the Name Your Project section.
- g. Write a short description of your project in the Description section (optional).
- h. Click “Continue.”

2. View and Build Sequences

- a. On the Assemble Sequences branch line, click “Sequence Viewer.” Click on a sequence name to view an electropherogram with a bar graph with quality scores for each nucleotide.
- b. The DNA sequencing software measures the fluorescence emitted in each of four channels – A,T,C,G – and records these as a trace, or electropherogram. In a good sequencing reaction, the nucleotide at a given position will be fluorescently labeled far more than the background. Thus, peaks in the electropherogram correlate to nucleotide positions in the DNA sequence. *Capture a print screen of a sample of great sequence in the electropherogram. Label it. Capture a print screen of a poor sequence, and copy it into your MS word document. Label it.*

Describe the quality of your sequencing results:

- c. A software program called Phred analyzes the sequence file and “calls” a nucleotide (A, T, C, G) for each peak. If two or more nucleotides have relatively strong signals at the same position, the software calls an “N” for an undetermined nucleotide.
- d. Phred also examines the peaks around each call and assigns a quality score for each nucleotide. The quality scores use a logarithmic scale to describe the probability that each nucleotide call is wrong, or, conversely, accurate.

Phred Score Error Accuracy

10	1 in 10	90%
20	1 in 100	99%
30	1 in 1,000	99.9%
40	1 in 10,000	99.99%
50	1 in 100,000	99.999%

- e. The bar below each nucleotide is the Phred score for that nucleotide. The horizontal line equals a Phred score of 20, which is generally the cut-off for high-quality sequence. Thus any bar at or above the line is considered a high-quality read.

What are the error rate and accuracy associated with a Phred score of 20? Discuss your results.

- f. Every sequence “read” begins with nucleotides (A,T,C,G) interspersed with Ns. In “clean” sequences, where experimental conditions were near optimal, the initial Ns will end within the first 25 nucleotides. The remaining sequence will have very few if any, internal Ns. Then, at the end of the read, the sequence will abruptly change over to Ns.
- g. Large numbers of Ns scattered throughout the sequence indicate poor quality sequence. Sequences with average Phred scores below 20 will be flagged with a “Low-Quality Score Alert.” You will need to be careful when concluding analyses made with the poor quality sequence. What do you notice about the electropherogram peaks and quality scores at nucleotide positions labeled “N”?
- h. Click on “Sequence Trimmer” to automatically remove Ns from the 5’ and 3’ ends of selected sequences. Once the program is finished, click again to view the trimmed sequences.

Why is it important to remove excess N’s from the ends of the sequences?

3. Pair Forward and Reverse Reads

- a. If you have a good quality forward and reverse reads (i.e. the sequence generated using the reverse primer) for any sample, click on “Pair Builder” to associate a forward read with its corresponding reverse read.
- b. Check the boxes for two sequences you wish to pair, and confirm your selection in the pop-up.
- c. Click on the “F” to the right of the reverse sequence. The entry will change to R, indicating that the sequence has been transformed into its reverse complement.
- d. Click on “Save” to save your pair assignments.
- e. Click on “Consensus Builder” to align the paired forward and reverse reads. Then, click on a forward-reverse pair to view its consensus sequence.

Why is the consensus sequence longer than the forward and reverse reads? Why does this occur?

- f. Positions highlighted in yellow mark differences in nucleotide calls between the forward and reversed reads.

Do differences tend to occur in certain areas of the sequence? Why?

- g. Large numbers of yellow mismatches – especially in long blocks – may indicate that you have incorrectly paired sequences from two different sources (organisms), or that you failed to reverse complement the reverse strand.
 - Return to Pair Builder to check your pairs and reverse complements.
 - Click on the red “x” to redo a pairing, and toggle “F” and “R” settings, as needed.
- h. A large number of mismatches in properly paired and reverse complemented sequences indicate that one or both sequences are of poor quality. Often, one of the sequencing reactions produces a high quality read that can be used on its own. To determine this:
 - Examine the distribution of Ns to see if they are mainly confined to one of the two sequences.
 - Examine the electropherograms to see if one of the two sequences is of good quality.
 - If one of the sequences seems of good quality, return to Pair Builder, and click the red x to undo the pairing.

Describe & explain your results:

- i. Few or no internal mismatches indicate good quality sequence from forward and reverse reads. If you like, you can check the consensus sequence at yellow mismatches and override the judgment made by the software:
 - Click on a highlighted mismatch to see the electropherograms and graphic summarizing Pfred scores for each read. Remember that the horizontal line equals a Phred score of 20, the cut-off for the high-quality sequence.
 - Click on the desired nucleotide in the black rectangle to change the consensus sequence at that position. You should only change the consensus if you have a strong reason to believe the consensus is wrong.
 - Click the button to “Save Change(s).”

4. **BLAST Your Sequence:** A BLAST search can quickly identify any close matches to your sequence in sequence databases. In this way, you can often quickly identify an unknown sample to the genus or species level. It also provides a means to add samples for phylogenetic analysis.

- a. On the Add Sequences branch, click on “BLASTN.” Then, click on the “BLAST” button next to the sequence you want to query against DNA databases.
- b. The returned list has information about the 20 most significant alignments (hits). For the top hit, copy the following into your MS word document:

- a. Accession number, a unique identifier given to each sequence submitted to a database. Prefixes indicate the database name – including gb (GenBank), emb (European Molecular Biology Laboratory), and dbj (DNA Databank of Japan).

What is the accession number for the top hit?

- b. Organism and sequence description or gene name of the hit. Click on the genus and species name for a link to an image of the organism, with additional links to detailed descriptions at Wikipedia and Encyclopedia of Life (EOL). Include an image and name of the plant which is the top hit in your MS word document.

What is the plant?

- c. Several statistics are shown in the window allow comparison of hits across different searches. The number of mismatches over the length of the alignment gives a rough idea of how closely two sequences match. The bit score formula takes into account gaps in the sequence; the higher the score, the better the alignment. The Expectation or E value is the number of alignments with the query sequence that would be expected to occur by chance in the database. The lower the E value, the higher the probability that the hit is related to the query. For example, an E value of 1 means that a search with your sequence would be expected to turn up one match by chance.

Why do the most significant hits typically have E values of 0? (This is not the case with BLAST searches with primers.) What does it mean when there are multiple BLAST hits with similar E values?

- d. Add BLAST sequence data to your phylogenetic analysis by checking the box (es) above any accession number(s), then clicking on “Add BLAST hits to project” at the bottom of the BLAST results window.

5. Add Sequences to Your Analysis

- a. Click on “Upload Data” to include additional data. Either upload data in ab1 or FASTA format or import data from other sources.
- b. Click on “Reference Data” to select data that will let you compare your barcode sequence in an appropriate phylogenetic context.

6. Analyze Sequences: Select and Align. Many unknown species can be rapidly identified by a BLAST search. In this case, a phylogenetic analysis adds depth to your understanding by showing how your sequence fits into a broader taxonomy of living things. If your BLAST search fails to identify your sequence, phylogenetic analysis can usually identify it to at least the family level.

- a. Click on “Select Data” on the “Analyze Sequences” branch. Then check the boxes to select any or all of the sequences you have uploaded from your own sequencing projects, from BLAST searches, and from reference data sets. Click on Save.
- b. Click on “MUSCLE” to align your sequences. When the program is finished, click again to view the alignment in Jalview.
 - a. Scroll through your alignments to see similarities between sequences. Nucleotides are color-coded, and each row of nucleotides is the sequence of a single organism or sequencing reaction. Columns are matches (or mismatches) at a single nucleotide position across all sequences. Dashes (-) are gaps in the sequence, where nucleotides in one sequence are not represented in other sequences.
 - b. Note that the 5’ (leftmost) and 3’ (rightmost) ends of the sequences are usually misaligned, due to gaps (-) or undetermined nucleotides (Ns). What causes these problems?
 - c. Note any sequence that introduces large, internal gaps (-----) in the alignment.
- c. This is either poor quality or unrelated sequence that should be excluded from the analysis. To remove it, return to “Select Data,” uncheck that sequence, and save your change. Then click on “MUSCLE” to recalculate.
- d. Trim Unaligned Ends of the Sequences

- a. Identify the leftmost point at which all or most sequences show corresponding nucleotide color bars. (There should be few or no gaps in the vertical column of nucleotides at this point.)
- b. Click in the nucleotide coordinate bar directly above this nucleotide in the first sequence. This will activate a red cursor and a pop-up menu.
- c. Click on “Remove left” to trim the leftmost sequences to this nucleotide position.
- d. Repeat first two steps of 6.c. above and click “Remove right” to trim the rightmost sequences.
- e. You can return to “Select Data” (in step b. above) to remove any sequence that has large sequence gaps. Why is it important to remove sequence gaps and unaligned ends?
- f. Click “Submit trimmed alignment.”

7. Analyze Sequences: Create a Phylogenetic Tree

NOTE: The Phylogenetic Tree software has been updated and may no longer be interactive. These instructions may not be accurate. Ask your instructor if you will complete this section.

- a. Click on “PHYML ML” to generate a phylogenetic tree using the maximum likelihood method. A tree will open in a new window, and the MUSCLE alignment used to produce it will open in another window.
- b. A phylogenetic tree is a graphical representation of relationships between taxonomic groups. In this experiment, a gene tree is determined by analyzing the similarities and differences in DNA sequence.
- c. Look at your tree.
 - a. The branch tips are the DNA sequences of individual species or samples you analyzed. Any two branches are connected to each other by a node (¶), which represents the common ancestor of the two sequences.
 - b. The length of each branch is a measure of the evolutionary distance from the ancestral sequence at the node. Species or sequences with short branches from a node are closely related; those with longer branches are more distantly related.
 - c. A group formed by a common ancestor and its descendants is called a clade. Related clades, in turn, are connected by nodes to make larger, clades.
 - d. Click on a node (¶) to highlight sequences in that clade. Click the node again to deselect the clade. What assumptions are made when one infers evolutionary relationships from sequence differences?
 - e. Generally, the clades will follow established phylogenetic relationships ascending from genus > family > order > class > phylum. However, gene and phylogenetic trees do disagree on some placements, and much research is focused on “reconciling” these differences.

Why do gene and phylogenetic trees sometimes disagree?

- d. Find and evaluate your sequence's position in the tree.
 - a. If your sequence is closely related to any of the references or uploaded sequences, it will share a single node with those species.
 - b. If your sequence is identical to another sequence, the two will diverge directly from the node without branches.
 - c. If your sequence is distantly related to all of the species in your tree, your sequence will sit on a branch by itself – with the other sequences grouping together as a clade.
 - d. To identify the smallest clade that includes your sequence, click on the node that is directly connected to your sequence. The sequences that are highlighted are the closest relatives of your sequence in the tree.
 - e. Look at the scientific names of sequences within the most closely associated clade. If all members share the same genus name, you have identified your sequence as belonging to that genus. If different genus names are represented, check and see if they belong to the same family or order.
- e. Return to the menu, and click on "PHYLIP NJ" to generate a phylogenetic tree using the neighbor-joining method. Include this in your report!

How does it compare to the maximum likelihood tree? What does this tell you?

- f. If neither tree places your sequence within an identifiable clade -- or if that clade is only at order level – you will need to add more sequences that may increase the resolution of your analysis. Return to Step 5, and add more reference sequences or obtain sequences within the order or family clade that contained your sequence. Then repeat Steps 6-7 to select, align, and generate trees from your refined data set.
- g. Print out your MSword document.
- h. **Staple your report in the following order:**
 - 1. Title page
 - 2. Manual up to and including Part I
 - 3. MS office work for Part I
 - 4. Manual of Part II
 - 5. MS office work for Part II
 - 6. Q&A pages below.



Lab Unit 13 – DNA Barcoding: Bioinformatics

- 1 What is DNA sequencing?
- 2 List all of the components of a Sanger Sequencing reaction and what they are for in the reaction.
- 3 What does a Phred score of 30 tell you?
- 4 Copy and paste your consensus sequence in an MSword document. Print it out and include it with your report.
- 5 What is the first (number 1) BLAST result? Is this result what you were expecting? Explain.

LAB UNIT 14: PREPARING FOR A CAREER IN BIOTECHNOLOGY

Jack O'Grady, M.S., & Linnea Fletcher, Ph.D.

OBJECTIVES

Your performance will be satisfactory when you can:

- ✓ Explore potential careers in biotechnology
- ✓ Create a resume

Exciting Biotechnology Careers

Biotechnology is a diverse career requiring a full range of academic and workforce skills.

Biotechnology offers career opportunities in:

- ✓ Bioscience: Medical, Agricultural, Environmental
- ✓ Applied chemistry: Testing
- ✓ Physics & Engineering: Biomedical devices
- ✓ Computer Science: Bioinformatics

Biotechnology Industry in Austin, Texas

The biotechnology industry has also been steadily growing in the Austin area. Today, Austin's bioscience community encompasses over 100 companies that employ more than 7000 people in the areas of research, diagnostics, pharmaceuticals and medical devices. Some of these businesses include Xbiotech, Agilent, Asuragen, Perkin Elmer, Fischer Scientific, Pfizer, CPL, and PPD to name a few!

Austin is also a major contributor to academic research in the biological sciences, at the University of Texas, Texas State University, and the M.D. Anderson Cancer Research Center in nearby Bastrop. Additionally, many students obtain jobs at the Texas Department of State Health Services where they test seasonal flu samples! To find out more biotech opportunities in Texas, explore here: <https://texaswideopenforbusiness.com/industries/biotechnology-life-sciences>
Explore Biotech jobs in Austin: <https://www.indeed.com/q-Biotech-l-Austin,-TX-jobs.html>

Biotechnology is a global endeavor! Many companies have research and development, manufacturing, or administrative facilities in several different countries. Bio-Link is a website that has an excellent biotechnology careers exploration tool. Go here to learn about potential careers in Biotechnology across the nation: <http://biotech-careers.org/>



PART I: Writing a Resume for a Job in the Biotechnology Industry

The purpose of a resume is to get an interview! Your resume will open doors for you in your career, so it is important you take time to create a memorable resume. It is a marketing brochure, and you are the one that you are marketing. While writing your resume, remember that your intention is to summarize the major highlights of your career in a way that will emphasize your professional accomplishments that are appropriate for the specific job you are applying.

It is required that you prepare a resume in advance of the internship course you will take at the end of the program. That is why you are being provided an opportunity to work on your resume in this course. You will continue to update your resume during every course with the new techniques you have learned. When you are ready to register for the internship program, you will have an accurate and professional resume ready to go! After instructor approval, this resume will be given to prospective mentors and employers.

CAREER WORKSHOPS: There are many free and excellent workshops available to you. Workforce Solutions has an office in the Eastview Campus, building 1000 and they have weekly workshops ranging from job hunting, resume building to mock interviews. It is recommended you visit their offices and see what it is they do there. You will be visiting them when you graduate to help you find a job. If you have never prepared a resume, the 'Basic Resume Lab' workshop is recommended. If you already have a working resume, the 'Advanced Resume Workshop,') is recommended. You can find a calendar of their workshops at their website:

<http://www.wfscapitalarea.com/>

RESUME BUILDING HIGHLIGHTS

1. **Targeted Resumes.** You must target each resume to each job you are applying. Use the job description to your advantage!
2. **Keep it brief.** Only include information that will promote you and give you an edge in getting an interview. Keeping it brief will keep it neat, organized, and easy to read. If you have no work history target one full page. If you have extensive work history, use two pages maximum.
3. **Make it easy to rapidly visually scan** for salient information by organizing your thoughts in outline form and bulleting the most important points in lists, so they are easy to find. Use large standard font type and large margins. Helvetica and Arial 11-12 are most recommended for resumes. Do not use times 10.
4. **List your strongest points first,** since many readers will only read the beginning of your resume. This is also true of your bulleted lists: list the most important things first.
5. **If you lack experience and accomplishments,** make more of the specific techniques that you have been trained. Focus attention on soft skills listed in the job advertisement (teamwork, organization skills, etc.) that you may have in your non-science job history.
6. **Use action-oriented verbs.** A strong action-oriented verb in the past tense is more interesting and makes an impression. For instance, use verbs such as "established," "directed," "created," "launched," "designed." Avoid weak verbs like "assisted," "served," "completed," "contributed." Also, avoid the vague phrases such as "was responsible for."

7. **Emphasize accomplishments.** An accomplishment is any way that you have brought value to your organization. You can do this in one of three ways: by helping the organization make money, by helping the organization save money, or by developing something new, innovative, or unique. Example: *“Coordinated the installation of a new inventory management system which accurately tracked inventory resulted in \$2000 saving in supplies each month.”*
8. **Don’t use excessive jargon.** Acronyms and jargon should be avoided. Use scientific jargon only if you can be quite certain that everyone who will read it will be familiar with the terms. NOTE: *Jargon used in the job description is encouraged*; the first people to scan resumes are typically not scientists but HR personnel. If you use it, spelling and using scientific jargon correctly is essential on a resume!
9. **Don’t be shy.** It is important for you to tactfully boast about how you have helped organizations, advanced scientific knowledge, and contributed to an important cause. Write your resume as if your parent wrote it!
10. **Proofread.** Nothing makes a worse impression than typos, misspellings, incorrect spacing, or inconsistencies in the use of bolding, italics, capitalization, and underlining. Proof your resume and have at least one other person review it.

RESUME CHECKLIST

- ✓ Make sure that your **name, one email, and one phone number** is prominently displayed at the top. No need to put your home address. Call yourself and ensure the voice mail is accessible and has a personal and professional message.
- ✓ **Summary.** It is preferable to provide a summary section (rather than an objective). A summary tells the company what you offer them. Tailor the summary to the job posting. The summary should include information that sets you apart from other applicants.
- ✓ **Your educational background.** List where you have attended schools and what degrees you have earned, and year awarded. If you are currently working on a degree, post-date your graduation. Bold the degree, not the university. No need to list your GPA.
- ✓ **List employers chronologically.** Unless there is a work history that relates to the job you’re applying for, list most recent first. Each entry should include the organization’s name, the department, years you were employed, and positions held. Emphasize any job experience that relates in any way to the job that you are now interested in, but include any job you have worked at for long periods of time even if they don’t have much to do with the job that you are now applying for. These are great opportunities to highlight applicable soft skills.
- ✓ **Use other categories if appropriate.** Any information that shows your strengths or gives a compelling reason for hiring you should be included
 - **Specific job skills and training:** As a student, your job experience in biotechnology may not be as impressive as your training. If this is the case, emphasize your training.

- **Awards, Scholarships, Publications & patents:** If you have more than a few awards and publications, you should emphasize them in their own section.
- **Presentations & posters:** Include workshops, presentations, seminars, and training programs you have delivered. Even if they are not Biotech-related, include them.
- **Professional Associations:** Do you belong to the ACC Biotech Club or the ASQ Quality association? This shows dedication to the field.
- **Community Service (volunteer work):** Have you spent time helping feed the homeless, or tutoring underprivileged children? This is a good place to show your employer not just a bit of your personality, but your character.
- **Avoid:** Avoid listing your hobbies unless the hobby relates to the job, then feel free to list it!

LAB UNIT 14: ASSIGNMENT

1. PRE-LAB EXERCISE – Before you come to class, perform the activities in the introduction and the following resume
 - Find your current resume and bring it to class. If you only have a paper copy, bring this to class. If you have an electronic copy, bring that as well, as you will use it to update your resume. If you don't have one, not to worry, you will create one in class!
 - Fill in the "Resume Worksheet" with new experiences to add to your resume, or if you do not have a working resume, use it to collect all the information you'll need to put on your resume. Bring this to class.
 - Using the "Resume Critique Sheet" provided, critique your current resume (if you have one). Bring this to class.
2. Based on the information provided in this lab, and from your instructor today, create a biotechnology-specific resume.
3. Print out your draft resume, and write DRAFT on the top.
4. Swap resumes with your classmate and 'edit' them using a red pen. Fill in the "Resume Critique Sheet" for your classmate. Return the edited resume and critique sheet to your classmate. Include this with your report.
5. Make corrections based on your classmate's comments if you see fit. Print out a final resume to include with this assignment.
6. Staple your final resume to the top of your assignment. Include all activities, rubrics, versions of your resume, and worksheets in order.

Sample Resume: ACC Student

[Email address](#) • 512.555.5555

Summary

Highly motivated and experienced biotechnician with a Bachelor's of Science in Biochemistry and intensive training in a biotechnology certification program seeks employment in the biotechnology industry. Experience in experimental protocol planning, execution, and documentation with attention to detail. Confident leadership and team work skills developed in the military.

Technical Skills Summary

- Protein: Protein purification using column chromatography (FPLC, HPLC); protein analysis using SDS-PAGE, western blot, ELISA, gas chromatography, and protein quantification with Bradford Assay and spectrophotometry.
- Nucleic acid technology: RNA isolation, genomic and plasmid DNA purification, PCR (end-point and real-time), agarose gel electrophoresis, mammalian cell culture, recombinant DNA technology, Southern and northern blotting, and VNTR analysis.
- Extensive media and solution preparation including multi-component solutions and buffers.

Education

Advanced Technical Certificate, Biotechnology, Austin Community College (2018)
BS, Biochemistry, Texas A&M, Corpus Christi (2009)

Work Experience

X Biotech, Research Intern

May 2017 – Sept 2017

- Performed ELISA, FPLC, SDS-PAGE, and western blotting to determine purity, identity and activity of recombinant antibodies.
- Performed PCR, plasmid preparation, and agarose gel analysis.
- Experience in ELISA and FPLC protocol development.

Austin Organics, Greenhouse Manager

2010 – 2016

- Managed the construction and day-to-day operations of the first organic certified hydroponic lettuce greenhouse in Texas.
- Handled direct customer sales, ordering materials, planting and harvesting schedule.
- Worked with Microbial Earth to culture aerobic microorganisms for composting.

United States Navy, Aviation Electronics Technician

1999 – 2003

- Worked as an electronics technician repairing and maintaining mission-critical avionics and automated test equipment.
- Instrumental in instituting new repair protocols.
- Responsible for maintaining documentation in accordance with a quality system.



RESUME WORKSHEET

To create your resume, it is helpful to first reflect upon your past and current experiences. When you begin you will quickly realize there is a large amount of information you will need to research to create a resume. Do not leave this to the last minute before class! You will need access to personal records; email, pay stubs, job descriptions, transcripts etc...

- 1 In the worksheet below list all work experiences, training and achievements. Try to identify what you did in that experience (see example).
- 2 After making your list, identify those experiences that you've enjoyed and/or learned something that you would like to use in the future. Circle or highlight those.
- 3 Highlight in a different color (or circle in red pen) those items that will require further research.
- 4 Include this worksheet in your report.

Work	Volunteer	Education	Other
<p><u>Example:</u> Shift supervisor at McDonald's -helped customers -handled \$\$\$ -supervised staff</p>	<p><u>Example:</u> Habitat for Humanity - 30 hours of service in 2012 - Recruited 3 new volunteers</p>	<p><u>Example:</u> Biotechnology certificate - Complete in 2014 - Favorite classes: BITC1240, BITC2441</p>	<p><u>Example:</u> -SSTEM Scholarship -Biotechnology club officer</p>



RESUME CRITIQUE

A great way to really learn about good resume preparation is reading other resumes. Swap your final resume with a classmate and critique according to the table below. Include your resume and the critique sheet from your lab partner in your final lab report.

RESUME CRITIQUE SHEET

Name: _____

Score: _____

CHARACTERISTICS	SCORE		FEEDBACK
	1.....5.....10 Weak Strong		
Overall appearance <ul style="list-style-type: none"> ✓ Easy to read & inviting ✓ Professional appearance 			
Format <ul style="list-style-type: none"> ✓ Contains contact information ✓ Complete ✓ Appropriate length ✓ Consistent layout/format 			
Organization <ul style="list-style-type: none"> ✓ Strongest qualities first ✓ Emphasizes important skills, credentials & accomplishments 			
Content <ul style="list-style-type: none"> ✓ Stresses skills, education, & accomplishments ✓ Targeted for job – supports summary/objective ✓ Appropriate content for resume 			
Language <ul style="list-style-type: none"> ✓ Appropriate tense ✓ Action verb- accomplishments ✓ Spelling ✓ Grammar 			

ADDITIONAL FEEDBACK:





APPENDIX





EMPLOYABILITY SKILLS EVALUATION FORM

AUSTIN COMMUNITY COLLEGE BIOTECHNOLOGY PROGRAM

Student Name: _____

Score: _____

Course: _____

Semester: _____

Skills	Needs Improvement		Average	Excellent	
	1	2		4	5
DEPENDABLE ATTENDANCE AND PUNCTUALITY, COMPLETES WORK IN A TIMELY MANNER	1	2	3	4	5
ORGANIZATIONAL SKILLS	1	2	3	4	5
WORKING WITH OTHERS, GOOD TEAM WORKER	1	2	3	4	5
SAFE & SECURE WORK HABITS IN A REGULATED ENVIRONMENT	1	2	3	4	5
COMMUNICATION SKILLS	1	2	3	4	5
RESOURCEFULNESS, ABLE TO WORK INDEPENDENTLY	1	2	3	4	5
DOCUMENTATION IN A REGULATED ENVIRONMENT (LAB NOTEBOOK, SOP, FORMS, BATCH BINDERS, LOG BOOKS)	1	2	3	4	5
KEEPS WORK AREA CLEAN AND ORDERLY	1	2	3	4	5
PRODUCES QUALITY WORK	1	2	3	4	5
PROBLEM SOLVING SKILLS	1	2	3	4	5

EVALUATOR: _____

DATE: _____

Additional Comments:





SOLUTION PREPARATION FORM

Control # _____

Name of Solution/Media: _____

Amount prepared: _____

Preparation Date: _____

Technician(s): _____

Component	Vendor/ lot #/ Control #	Date Received	Storage conditions	MW or initial concentration	Mass used	Final Concentration

Balance used	Calibration status	
pH meter used	Calibration status	
Initial pH	Final pH	Adjusted pH with
Preparation temperature	Sterilization procedure/ sterility testing	Solution storage conditions & location

Calculations/Comments:



For every solution prepared in the laboratory, you must fill out a solution preparation form. You will find the forms in the file cabinet. A sample formed filled out is below. Never leave a blank space in a form – every field must be filled in, even if it was ‘not performed’ state so.

SOLUTION PREPARATION FORM

Control # Use label SOP

Name of Solution/Media: Complete name, include concentration, pH

Amount prepared: _____

Date: _____

Preparers(s): If two students prepared the solution, both names are recorded here. Both students must submit a copy of this (same) Solution Prep form in their lab report

Component	Vendor/ lot #/ Control #	Date Received	Storage conditions	MW or initial concentration	Mass used	Final Concentration
<i>Water is not a listed component</i>						

Balance used Balance number	Calibration status Did you calibrate it? Or did you verify it?	
pH meter used pH meter number	Calibration status Did you calibrate it? Or did you verify it?	
Initial pH Always record the pH the solution started at	Final pH Always record your final pH after you BTV	Adjusted pH with Chemical and concentration you used to adjusted the pH
Prep temperature Record the actual temperature the solution was prepared – NOT “RT”	Sterilization procedure Was this solution autoclaved or filter sterilized?	Storage conditions Where is this solution stored now?

Calculations/Comments:

ALL calculations used to prepare this solution are recorded here.



ELECTROPHORESIS DOCUMENTATION FORM

Date _____ buffer _____ gel__ % agarose/acrylamide (**circle one**)

Voltage _____ start time _____ stop time _____ stain _____

Analyst(s): _____

Experimental Summary/Description:

Lane	Sample description	Concentration	Sample Volume Analyzed	Mass loaded
1				
2				
3				
4				
5				
6				
7				
8				
9				
10				
11				
12				

INSERT CLEARLY LABELED IMAGE (S) HERE

TAPE ALL 4 EDGES, AND INITIAL ACROSS THE TAPE

Photographic settings (exposure time, aperture setting, etc):





DEVIATION REPORT FORM

Date: _____ DR#: _____

Detailed name & Model # of Equipment: _____

Serial Number or Equipment Number: _____

Where is the Equipment now? _____

Detailed Description of Deviation:

Recommended or Actual Correction of Deviation:

Date: _____ Name: _____ Signature: _____

QA Approval (Faculty Member or Laboratory Technician)

Date: _____ Signature: _____

Instructions: Keep a copy for your lab notebook, keep original with the equipment and notify lab technician where the equipment is.





BIOL1414: LABORATORY REPORT RUBRIC



Name _____ Date _____

Lab _____ Score _____

Formatting:

Title Page	5 pts _____
Overall neatness/legibility	5 pts _____
Grammar and Mechanics	5 pts _____

Work Book Exercises: Complete all exercises in the Work Book.

Complete all exercises in the workbook	20 pts _____
--	--------------

Lab Work:

Team Work: Active participation in lab exercise	10 pts _____
Thorough observations made in (written in ink)	10 pts _____
Safety: Wear PPE, lab clean up, operate equipment	10 pts _____

Results & Analysis:

Tables, Graphs, Calculations, Forms	15 pts _____
Lab Analysis Questions	15 pts _____
Conclusion	5 pts _____

Total:	100 pts _____
---------------	----------------------

Feedback:



SOP TEMPLATE & INSTRUCTIONS

Title: How to Write a Standard Operating Procedure (SOP)	
Institution: Austin Community College	Department: Biotechnology Program
Approved By:	Date of Approval: 10/07/2017
Prepared By: Jack O'Grady	Revision Number: 002

1.0 SCOPE AND APPLICATION

The purpose of this Standard Operating Procedure (SOP) is to train users to write an SOP for use in ACC Biotechnology labs and teaching curriculum. SOPs provide consistency each time a procedure or process is performed, serve as reminders to ensure that work is completed correctly, can be used to train employees the correct way to complete a new task, and also reduce the possibility of failure by enabling the employee to complete any function described in the SOP.

2.0 SUMMARY OF METHOD

Briefly, summarize the procedure.

3.0 WORKFLOW

An overview of the workflow. Create a flowchart of the procedure.

4.0 RESPONSIBILITIES

Include qualification user should have to complete task satisfactorily. A technician or student is responsible for adhering to all applicable duties outlined in this SOP. The instructor is responsible for overseeing all activities, ensuring all work satisfies the specific tasks described in this SOP.

5.0 DEFINITIONS

Identify specialized terms, abbreviations, or acronyms used in this SOP.

6.0 SAFETY GUIDELINES

6.1 Health and Safety Warnings

Indicate operations that could result in injury. Explain what will happen if the procedure is not followed correctly. List warnings here, and at critical steps, in the procedure.

6.2 Precautions

Indicate activities that may result in equipment damage, sample degradation, or possible invalidation of results. List here and at critical steps in the procedure.

7.0 WASTE MANAGEMENT

Hazardous materials disposal guidelines.



8.0 INTERFERENCES

Describe process components that may interfere with accuracy of final product

9.0 QUALITY CONTROL

This section dedicated to allowing self-verification of the quality and consistency of the work. Describes preparation of appropriate QC procedures and QC material for calibration, performance evaluations, standard preparation, and the frequency at which they should occur. Also, details limits/criteria for QC results and actions required when data exceed limitations.

10.0 REAGENTS, MATERIALS, AND EQUIPMENT

10.1 List reagents: name, manufacturer, CAS, catalog number (if known).

10.2 List specific instruments/software used in this SOP: vendor and model number. Include SOPs to use and calibrate this equipment.

11.0 CALCULATIONS

Show mathematical steps to be followed for preparation of procedure.

12.0 PROCEDURE

Use this space for detailing the experiment procedure of each step listed in the workflow. Must be numbered steps.

1. Assign a title for the SOP that briefly and clearly states what it will describe.
 - a. The title reflects a full description of the purpose of the procedure.
 - b. The title is written directly into Title Box at the top of an SOP.
2. Assign an ID number and revision number to the SOP.
 - a. The ID number reflects the type of procedure. Example: SOP-HPLC-xxx refers to SOPs for HPLC instruments. SOP-HPLC-001 refers to a specific SOP written for an HPLC procedure. This section will also list the date when the document was approved.
 - b. Revision numbers indicate how many times the document has been revised. The writer begins with the letter A and proceeds alphabetically each time a new version is approved. This section will also list the date when the document was revised and approved again.

Use this space to take notes during your experiments or other items to note during the use of this SOP

SOP number is inserted into the header of this template.

<p>c. SOP ID number is placed at the top right-hand corner of the SOP labeled “SOP#.”</p> <p>3. In Section 1.0 Scope and Application, write a statement of the purpose and scope of the procedure. Scope should be brief but descriptive and specific, include situations when this SOP is used.</p> <p>4. In Section 2.0 Summary of Method, briefly, summarize the procedure.</p> <p>5. In Section 3.0 Workflow, briefly summarize the workflow, using a list of steps or a graphical representation such as a flowchart.</p> <p>6. In Section 4.0 Responsibilities, define the responsibility of the procedure. Include qualifications and credentials user should have to complete task satisfactorily.</p> <p>7. In Section 5.0 Definitions, include any specialized terms that will be needed to follow the SOP.</p> <p>8. In Section 6.0 Safety Guidelines, define as follows:</p> <ol style="list-style-type: none"> Identify any hazardous chemicals and how they might enter the body or otherwise cause danger. Example: <i>“Concentrated sulfuric acid released strong fumes that cause lung, eye, and skin irritation and can result in serious burns if spilled on the body. Only work with concentrated sulfuric in a fume hood while wearing nitrile gloves, lab coat, and safety glasses.”</i> Indicate what personal protective equipment (PPE) is required for performing the procedure. Examples: <i>“Sodium dodecyl sulfate causes lung, eye, and skin irritation. Wear a respirator while weighing the dry chemical.”</i> Identify other physical hazards (flammability, electrical, mechanical) and how to work safely to reduce the risk of injury. Example: <i>“Hydrogen gas is extremely flammable if concentrated. Avoid any open flame around hydrogen tank.”</i> Indicate proper disposal of biohazards, chemicals, and contaminated items. Example: 	
--	--

“Place used agar plates in a red biohazard bag, tie and deposit in large cardboard biohazard collection area. Place pipet tips contaminated with 2-mercaptoethanol in the designated collection bag inside the fume hood.”

9. In Section **7.0 Waste Management**, include instructions for how to safely dispose of any hazardous materials used in, or generated by, the procedure.
10. In Section **8.0 Interferences**, describe any known processes that may interfere with the accuracy of the final product.
11. In section **9.0 Quality Control**, list instructions related to quality control of the procedure, instrument, or materials used for the SOP.
12. In Section **10.0 Reagents, Materials, and Equipment**, list materials are needed to perform the procedure.
 - a. These can include equipment, supplies, chemicals, and facilities.
 - b. Be specific about the manufacturer and model of products if a specific model is required.
 - c. Specify temperatures for water baths and incubators.
 - d. Include SOP of associated equipment use and calibration.
13. In Section **11.0 Calculation**, include an example of the calculation needed to perform this task and provide space for the user to write any calculations needed.
14. In section **12.0 Procedure**, write a procedure in numbered steps that provides instructions to the user exactly how to perform the operation.
 - a. Steps should be written as commands in the present tense.
 - b. Steps should be placed in chronological order, particularly when an advanced preparation is required. Example: *“One hour before starting procedure, place solution A at room temperature to equilibrate.”*

<p>c. Each step should provide only the details required for that step, minimizing background information or explanation.</p> <p>d. If a common problem arises during the procedure, include a step for how to prevent or correct the problem. Example: <i>“The pellet may become dislodged from the tube if it was not dense enough. Centrifuge again for two minutes, increasing speed to 10,000rpm.”</i></p> <p>15. In Section 13.0 Data and Records Management, list instructions for how to manage data from, or records related to, the procedure being followed. Provide guidelines for how to document the procedure was performed; this may be a form that was filled out, log book, or recording the action in a lab notebook.</p> <p>16. In Section 14.0 References, cite sources you used to write the SOP; this is often the user manual for a piece of equipment, a product insert, or another SOP or procedure.</p> <p>17. Submit the written SOP to the supervisor for approval.</p> <p>18. Destroy or limit access to previous versions of the SOP. Always keep a copy of the old versions for historical reference.</p>	
---	--

13.0 DATA AND RECORDS MANAGEMENT

The records section should include post-data calculations or presentation to be performed, forms to be used, and include data and record storage information.

14.0 REFERENCES

List any references, associated documents, and forms.

Seidman LA, Moore CJ. 2009. *Basic laboratory methods for biotechnology: textbook and laboratory reference*. Upper Saddle River, NJ: Prentice Hall.



